

# Passion in Health. Since 1988.

Product range

## Meditrade UK

Meditrade UK is a member of the Meditrade Group located in central London. Meditrade Group has always had a large number of UK customers, however Brexit meant that the quality of our oversea relationships were going to decline, naturally. Because of this Meditrade made the decision to open an office on the ground in the UK, which would allow Meditrade UK to offer favourable business solutions to the existing UK customer base. In addition, we can offer Irish customers services from Europe (Germany – HQ).

At Meditrade UK we have developed a logistic concept that offers a flexible and consistent supply of our products into the UK and Irish market. Part of our logistics concept is the multi-stage logistics hub, which is being built in Gainsborough, Lincolnshire. The strategically central geographical placement establishes an efficient distribution network which allows us to promptly assist our customer's demands at all times.





## Meditrade – 30 years of innovation

In 1988 Meditrade was founded under the maxim 'Passion in Health '- Meditrade Group combines more than 30 years of expertise in manufacturing, distributing and supplying medical consumables to healthcare facilities, manufacturers, distributors, retailers and wholesalers. Starting as a manufacturer of medical disposable gloves, Meditrade has now become one of the leading manufacturers of medical equipment in Europe with a wide range of products from wound care, protective clothing, operating theatre sets, incontinence products and disinfection products for the skin, hands, and surfaces.

As manufacturers and system suppliers of medical consumables, we see ourselves as partners and service providers for the basic care of all medical facilities. Quality, innovation and beneficial care options are our daily motivation to make a small but essential contribution to modern and future-oriented healthcare. The employees of Meditrade stand by this core essence and fulfil these demanding tasks for a higher standard of health with competence and professionalism. With factories all over the world and logistics centres in Germany and the United Kingdom, Meditrade meets the increasing demands and requirements of its customers and always assists its customer's requests wherever they are.



## Our product range

- Gloves
- Protective clothing
- Wound care
- Sets and operating theatre supplies
- Incontinence
- Disinfection and care
- Hygiene, cleaning and cleaning products
- Ward and patient accessories


## **Logistics**

## Our approach to logistics ensures the flexible and reliable supply of our products to the entire European market.

- Main site with storage space for 20,000 pallets
- Our expanded logistics centre in Neuenkirchen offers space for 28,000 pallets
- 12.5-metre high racking with a semi-automatic picking system
- Completion of our logistics centre in Gainsborough in 2022



## Locations




## About us

- Meditrade products are manufactured in our own facilities all over the world in order to ensure the high quality of our products
- We offer a comprehensive product portfolio to supply our customers with a full range of medical consumables
- Leading provider of PPE and medical consumables in Europe, with over thirty years of experience in the industry
- A global company that stays true to the values of a family-run business
- >150 employees



## **Our core competencies**



#### Gloves

010763	
Latex examination gloves	
Reference™	9
Gentle Skin <sup>®</sup> sensitive	9
Gentle Skin <sup>®</sup> classic (x-long)	10
Gentle Skin <sup>®</sup> grip	11
Gentle Skin <sup>®</sup> compact +	11
Gentle Tec®	12
Gentle Skin <sup>®</sup> Aloecare <sup>™</sup>	13
Gentle Skin <sup>®</sup> black	14
Nitrile examination gloves	
Nitril <sup>®</sup> 3000	15
Nitril® NextGen®	15
Nitril <sup>®</sup> Sensory <sup>®</sup>	16
Nitril® BestGen®	17
Nitril <sup>®</sup> 3000 x-long blue	18
Nitril <sup>®</sup> black	19
Nitril <sup>®</sup> Viola	20
Nitril <sup>®</sup> Magenta	21
Nitril <sup>®</sup> LIOX	22 / 23
Vinyl examination gloves	
Vinyl 2000	24
Vinyl Blue Gen PF	25
Vinyl 2000 Stretch PF	25
Cotton and hygienic gloves	
Yarn glove	26
Froggy™	27
Finger cots	
Latex finger cots	28
Household gloves	
MyClean <sup>®</sup> easy touch	29
Sterile examination gloves	
Copolymed <sup>®</sup>	30
Gentle Skin <sup>®</sup> steril	31
Nitril <sup>®</sup> steril	32
Surgical gloves	
Reference <sup>™</sup> OP	35
Gentle Skin <sup>®</sup> Premium OP™	36
Gentle Skin <sup>®</sup> Superior OP	37
Gentle Skin <sup>®</sup> Smooth OP	38
Gentle Skin <sup>®</sup> Micro OP <sup>®</sup>	39
Maxitex <sup>®</sup> UG	40
Gentle Skin <sup>®</sup> Securitex <sup>®</sup> OP	41
Gentle Skin <sup>®</sup> Isopretex <sup>®</sup> OP	42
Neopretex <sup>®</sup> OP	43
Finessis <sup>®</sup> Aegis <sup>®</sup>	45
Finessis <sup>®</sup> Zero <sup>®</sup>	46
Finessis <sup>®</sup> Corium <sup>®</sup>	47

#### Protective clothing

Masks	
Suavel <sup>®</sup> Protec Plus	50
Suavel <sup>®</sup> Protec	51
Suavel <sup>®</sup> Protec Black	51
Suavel <sup>®</sup> Comfort Plus	52
Suavel <sup>®</sup> Comfort	52
Suavel <sup>®</sup> Antifog	53
Suavel <sup>®</sup> Antifluid	53
Suavel <sup>®</sup> Sensima	54
Caps	
Suavel <sup>®</sup> Apollo	55
Suavel <sup>®</sup> Astrid	55
Scrubs / gowns / aprons	
BeeSana <sup>®</sup> PP gown	56
BeeSana <sup>®</sup> PP/PE gown	57
BeeSana <sup>®</sup> SMS gown	57 / 58
Suavel <sup>®</sup> Yeti Coverall	58
PE disposable apron	59
Suavel <sup>®</sup> OP PRO Standard / Expert	60 / 61

#### Wound care

Gauze64BeeSana® Gauze cutting64Gauze Swabs & Dressing pads65BeeSana® Sit gauze swab68BeeSana® Sit gauze swab69BeeSana® Sit non-woven swab70ABE® drainal71BeeSana® Sit non-woven swab72BeeSana® Sit non-woven swab73ABE® tracheal71BeeSana® Sit non-woven swab72BeeSana® Dep pad73ABE® nose bandage74Swabs and tamponades75BeeSana® Cauze balls76BeeSana® Cauze balls76BeeSana® Pointed swabs76ABE® nose bandage77Alcohol swabs77Alcohol swabs77Vooden tongue depressor78Wooden tongue depressor78Wooden tongue depressor78BeeSana® Cauze balis82Fixing plasters/non-woven materials82ABE® derm plus82Fixing plasters/non-woven materials83ABE® andula retention dressing83ABE® andula ges self-adhering, crepe88BeeSana® SafeCut85ABE® last drepe86ABE® last drepe88ABE® dideil 10091ABE® added bandage92ABE® lastic93ABE® dideil tondage94ABE® extra adhesive95ABE® last drepice96ABE® last chere pressure bandage96ABE® last drepice96ABE® last drepice96		4
BeeSana® Medical gaze64Gauze Swabs & Dressing pads65 - 67BeeSana® Gauze swab68BeeSana® Silt gauze swab69BeeSana® Silt non-woven swab70ABE® drainal71ABE® drainal71BeeSana® Silt non-woven swab70BeeSana® Silt non-woven swab71BeeSana® Silt non-woven swab72BeeSana® Absorbent dressing pad73BeeSana® Cauze balls75BeeSana® Cauze balls75BeeSana® Cellulose swabs76ABE® draine care and accessories77Cotton buds78Wound plasters and bandages88ABE® injection plaster80ABE® dram plaster short adhesive dressing81ABE® dram plasters/nor-woven materials82Fixing plasters/nor-woven materials83ABE® accold88ABE® derm83ABE® accold88ABE® dict fixing bandages, self-adhering, crepeBeeSana® SafeCut88ABE® fix color88ABE® fix color88ABE® fix color89ABE® fix color89ABE® fix color89ABE® fix color88ABE® fix color88ABE® fix color88ABE® fix color88ABE® fix color88ABE® fix color88ABE® fix color89ABE® fix color89ABE® fix color88ABE® fix color88ABE® fix color<		
Gauze Swabs & Dressing pads65 - 67BeeSana® Gauze swab669BeeSana® Silt gauze swab669BeeSana® Non-woven swab70ABE® 'tracheal71BeeSana® Silt non-woven swab70ABE® 'tracheal71BeeSana® Lye pad73ABE® 'tracheal71BeeSana® Lye pad73ABE® 'tracheal71BeeSana® Cauze balls75BeeSana® Cellulose swabs76ABE® 'trachol swabs77Alchol swabs77Alchol swabs77Alchol swabs78Wooden tongue depressor79Wound plasters and bandages88ABE® 'libit and basive dressing81ABE® 'alchou' woven materials82ABE® 'alchou' woven fabric83ABE® 'alchou'85ABE® 'alchou'85ABE® 'alchou'85ABE® 'alchou'85ABE® 'alchou'85ABE® 'alchou'85ABE® 'alchou'85ABE® 'alchou'85ABE® 'alchou'89Oropression & Support bandages89Oropression & Support bandage90Compression & Support bandage94ABE® 'alchir's bandage94ABE® 'alchir's bandage96ABE® 'alchir's bandage96ABE® 'alchir's bandage97ABE® 'alchir's bandage97ABE® 'alchir's bandage97ABE® 'alchir's bandage96ABE® 'alchir's bandage		-
BeeSana® Gauze swab65 - 67BeeSana® Silt qauze swab69BeeSana® Silt non-woven swab70ABE® drainal71ABE® tracheal71BeeSana® Lipe pad73ABE® drainal71BeeSana® Lipe pad73ABE® drainal71Swabs and tamponades76BeeSana® Cauze balls76BeeSana® Cauze balls76BeeSana® Cauze balls76BeeSana® Cauze balls77BeeSana® Cauze balls77BeeSana® Cauze balls78Wooden tongue depresor79Wound platers and bandages81ABE® injection plaster80ABE® derm plus81ABE® derm plus82Fixing plasters/non-woven materials83ABE® non-woven fabric83ABE® ancolucages, self-adhering, crepe88BeeSana® SafeCut87ABE® quize bandage86ABE® fix color89ABE® fix color89ABE® lastic90Compression & Support bandages90Compression & Support bandage94ABE® lastic93ABE® lastic93<		64
BeeSana® Slit gauze swab68BeeSana® Slit non-woven swab69BeeSana® Slit non-woven swab70ABE® drainal71ABE® drainal71BeeSana® Absorbent dressing pad72BeeSana® Cauze balls73BeeSana® Gauze balls75BeeSana® Context State76BeeSana® Context State76BeeSana® Context State77ABE® ribbon gauze76BeeSana® Cellulose swabs77Alcohol swabs77Alcohol swabs77Cotton buds78Woond platers and bandages78ABE® injection plaster80ABE® derm81ABE® derm plus82Fixing plasters/non-woven materials83ABE® derm plus83ABE® derm plus83ABE® derm plus84BetSana® SafeCut85ABE® dict fixing bandages, self-adhering, crepeBeeSana® SafeCut88ABE® fix color88ABE® fix color88ABE® fix color89ABE® fix color89ABE® fix color89ABE® lastic93ABE® lastic93ABE® lastic93ABE® lastic94ABE® plastic93ABE® plastic93ABE® plastic93ABE® lastic93ABE® lastic93ABE® lastic93ABE® plastic93ABE® plastic93ABE® plastic		65 67
BeeSana® Non-woven swab69BeeSana® Silt non-woven swab70ABE® tracheal71BeeSana® Absorbent dressing pad72BeeSana® Eye pad73ABE® nose bandage74Bwabs and tamponades76BeeSana® Cauze balls75BeeSana® Cellulose swabs76ABE® ribon gauze76BeeSana® Cellulose swabs77Alchol swabs77Alchol swabs77Alchol swabs77Alchol swabs77Mound plasters and bandages80ABE® injection plaster80ABE® injection plaster80ABE® derm82ABE® derm plus81ABE® derm plus81ABE® derm plus81ABE® derm plus83ABE® silk84Elastic fixing bandage, self-adhering, crepe83Besana® SafeCut85ABE® last87ABE® last87ABE® last87ABE® lasti87ABE® lasti dahesive90ABE® lasti dahesive92ABE® lasti dahesive93ABE® lasti dahesive94ABE® lasti dahesive95ABE® lasti dahe		-
BeeSana® Slit non-woven swab70ABE® drainal71ABE® tracheal71BeeSana® Absorbent dressing pad72BeeSana® Eye pad73ABE® nose bandage74Swabs and tamponades75BeeSana® Gauze balls76BeeSana® Cellulose swabs76ARE® ribbon gauze76BeeSana® Cellulose swabs77Alcohol swabs77Treatment care and accessories77Cotton buds78Wooden tongue depressor78Wooden tongue depressor80ABE® light80ABE® light80ABE® derm plus82Fixing plasters'non-woven materials82ABE® derm plus83ABE® derm novoven fabric83ABE® silk84Elastic fixing bandages, self-adhering, crepe86Bestana® SafeCut85ABE® last87ABE® last crepe88ABE® last crepe88ABE® last crepe88ABE® lastic90Compression & Support bandages92ABE® ven95ABE® ven95ABE® heart catheter pressure bandage96ABE® heart catheter pressure bandage98ABE® heart catheter pressure bandage99ABE® heart catheter pressure bandage99ABE® heart catheter pressure bandage99ABE® heart catheter pressure bandage99ABE® heart catheter pressure bandage90ABE® hea		-
ABE® drainal71ABE® tracheal71BeeSana® Absorbent dressing pad72BeeSana® Cyep ad73ABE® nose bandage74Swabs and tamponades75BeeSana® Cauze balls75BeeSana® Collulose swabs76ABE® ribbon gauze76BeeSana® Cellulose swabs77Alcohol swabs77Alcohol swabs77Alcohol swabs78Wooden tongue depressor79Wound plasters and bandages81ABE® injection plaster80ABE® injection plaster81ABE® derm plus82Fixing plasters/non-woven materials83ABE® derm plus82Fixing plasters/non-woven materials84ABE® derm plus83ABE® date subadge86ABE® fix color89ABE® date subadge86ABE® last crepe88ABE® last crepe88ABE® last crepe88ABE® last crepe88ABE® last crepe90Compression & Support bandage92ABE® last dhesive93ABE® last dhesive93ABE® vertra94ABE® heart catheter pressure bandage96ABE® heart catheter pressure bandage98ABE® heart catheter pressure bandage98ABE® heart catheter pressure bandage98ABE® last chelage99ABE® last chelage99ABE® last chelage99ABE® last		-
ABE® tracheal71BeeSana® Absorbent dressing pad72BeeSana® Eye pad73ABE® nose bandage74Swabs and tamponades75BeeSana® Cauze balls75BeeSana® Cauze balls76BeeSana® Cellulose swabs77Alcohol swabs77Alcohol swabs77Alcohol swabs77Alcohol swabs77Alcohol swabs77Mound plasters and bandages78Wooden tongue depressor79Wound plasters and bandages81ABE® hipetchin plaster80ABE® Absorbent adhesive dressing81ABE® cannula retention dressing83ABE® cannula retention dressing83ABE® cannula retention dressing83ABE® last87ABE® last87ABE® last color89ABE® last color89ABE® lastic99Compression & Support bandages94ABE® padded bandage94ABE® non-woven fabric89ABE® last color89ABE® last color89ABE® last color99ABE® last adhesive90Compression & Support bandages94ABE® padded bandage94ABE® padded bandage96ABE® lastic93ABE® bolty bandage96ABE® bolty bandage96ABE® libh older100ABE® libh older100ABE® libh older100ABE® li		
BeeSana® Absorbent dressing pad72BeeSana® Eye pad73ABE® nose bandage74BeeSana® Gauze balls75BeeSana® Caluze balls76BeeSana® Cellulose swabs76Alcohol swabs77Alcohol swabs77Alcohol swabs77Alcohol swabs77Treatment care and accessories77Cotton buds78Wooden tongue depressor78Wound plasters and bandages81ABE® injection plaster80ABE® derm plus82Fixing plasters/non-woven materials82ABE® derm plus82Fixing plasters/non-woven materials84ABE® derm plus82BAE® gauze bandage86ABE® gauze bandage86ABE® fix color89ABE® fix color89ABE® fix color89ABE® fix color89ABE® lasti90Compression & Support bandages90ABE® lasti92ABE® lasti93ABE® lasti93ABE® hadded bandage94ABE® ven95ABE® lubular bandage96ABE® restra adhesive90ABE® lubular bandage96ABE® libholder100ABE® holded bandage97ABE® holded bandage98ABE® holded bandage97ABE® holder100ABE® holder100ABE® holder100 <trr>ABE® holder<td< td=""><td></td><td></td></td<></trr>		
BeeSana® Eye pad73ABE® nose bandage74Swabs and tamponades75BeeSana® Gauze balls75BeeSana® Cellulose swabs76ABE® ribbon gauze76BeeSana® Cellulose swabs77Alcohol swabs77Treatment care and accessories77Cotton buds78Wooden tongue depressor79Wound plasters and bandages81ABE® injection plaster80ABE® injection plaster81ABE® derm plus82Fixing plasters/non-woven materials83ABE® cannula retention dressing83ABE® duru85ABE® gauze bandage86ABE® fix color88ABE® fix color89ABE® lubal tack repe88ABE® lubal tack repe88ABE® lubal tack repe88ABE® lastic90Compression & Support bandages90ABE® lastic93ABE® lastic93ABE® lastic93ABE® lastic93ABE® lastic93ABE® padded bandage94ABE® restra adhesive99ABE® lubular net bandage96ABE® pressure bandage roll100ABE® pressure bandage roll101Stortum roll101Stortum roll101Restana® Cauze balls X-ray102 / 103BeeSana® Lap sponge X-ray102 / 103BeeSana® Lap sponge X-ray102 / 103BeeSana® Cauze balls X-r		-
ABE® nose bandage74Swabs and tamponades75BeeSana® Gauze balls75BeeSana® Onited swabs76ABE® ribbon gauze76BeeSana® Cellulose swabs77Alcohol swabs77Alcohol swabs77Alcohol swabs77Mound plasters and bandages80ABE® injection plaster80ABE® light80ABE® derm812ABE® derm82ABE® derm plus82Fixing plasters/non-woven materials83ABE® data retention dressing83ABE® cannula retention dressing83ABE® silk84Elastic fixing bandage, self-adhering, crepe88BeeSana® SafeCut85ABE® last crepe88ABE® last crepe88ABE® last crepe88ABE® last crepe90Compression & Support bandages92ABE® lastic93ABE® lastic93ABE® lastic93ABE® tubular bandage94ABE® tubular bandage94ABE® retra adhesive95ABE® tubular bandage97ABE® heart catheter pressure bandage90ABE® heart catheter pressure bandage90ABE® heart catheter pressure bandage91ABE® heart catheter pressure bandage94ABE® heart catheter pressure bandage96ABE® heart catheter pressure bandage90ABE® heart catheter pressure bandage90ABE	01	-
Swabs and tamponades75BeeSana® Gauze balls76BeeSana® Cellulose swabs76ABE® ribbon gauze76BeeSana® Cellulose swabs77Alcohol swabs77Alcohol swabs77Treatment care and accessories78Cotton buds78Wooden tongue depressor79Wound plasters and bandages80ABE® injection plaster80ABE® light81ABE® derm82ABE® derm plus82Fixing plasters/non-woven materials82ABE® cannula retention dressing83ABE® silk84Elastic fixing bandages, self-adhering, crepe86BeeSana® SafeCut886ABE® last87ABE® last crepe88ABE® last crepe88ABE® last color99ABE® last color99ABE® lastic91ABE® lastic93ABE® lastic93ABE® lastic93ABE® lastic94ABE® ven95ABE® tubular bandage96ABE® tubular bandage96ABE® tubular bandage96ABE® tubular bandage96ABE® peaded bandage96ABE® tubular bandage96ABE® tubular bandage96ABE® tubular bandage97ABE® body bandage98ABE® tubular bandage99ABE® tubular bandage96ABE® tubular bandage96A		-
BeeSana® Gauze balls75BeeSana® Pointed swabs76BeeSana® Cellulose swabs77Alcohol swabs77Alcohol swabs77Treatment care and accessories78Cotton buds78Wound plasters and bandages80ABE® injection plaster80ABE® derm82ABE® derm plus82Fixing plasters/non-woven materials83ABE® cannula retention dressing83ABE® alge color83ABE® alge color88ABE® alge color88ABE® last86ABE® last crepe88ABE® last color88ABE® last color89ABE® universal90Compression & Support bandages91ABE® extra94ABE® tubular net bandage94ABE® tubular net bandage94ABE® tubular net bandage95ABE® tubular net bandage96ABE® heart catheter pressure bandage96ABE® heart catheter pressure bandage96ABE® heart catheter pressure bandage90ABE® heart catheter pressure bandage90ABE® heart catheter pressure bandage90ABE® heart catheter pressure bandage91ABE® heart catheter pressure bandage96ABE® heart catheter pressure bandage90ABE® heart catheter pressure bandage90ABE® heart catheter pressure bandage90ABE® heart catheter pressure bandage90ABE® heart c		
ABE® ribbon gauze76BeeSana® Cellulose swabs77Alcohol swabs77Treatment care and accessories78Wooden tongue depressor79Wound plasters and bandages80ABE® injection plaster80ABE® light81ABE® derm82ABE® derm plus82Fixing plasters/non-woven materials82ABE® cannula retention dressing83ABE® silk84Elastic fixing bandages, self-adhering, crepe86BeeSana® SafeCut85ABE® lust color89ABE® lust color89ABE® last color89ABE® last color90Compression & Support bandages91ABE® padded bandage94ABE® vena95ABE® vertra95ABE® tubular bandage96ABE® tubular net bandage96ABE® tubular net bandage97ABE® body bandage98ABE® bidy bandage98ABE® ibidichrist's bandage98ABE® heart catheter pressure bandage90ABE® heart catheter pressure bandage90ABE® heart catheter pressure bandage90ABE® heart catheter pressure bandage91ABE® heart catheter pressure bandage96ABE® heart catheter pressure bandage90ABE® heart catheter pressure bandage90ABE® heart catheter pressure bandage100ABE® heart catheter pressure bandage90ABE® heart catheter p		75
BeeSana® Cellulose swabs77Alcohol swabs77Alcohol swabs77Treatment care and accessories78Cotton buds78Wound plasters and bandages80ABE® injection plaster80ABE® light81ABE® derm plus82Fixing plasters/non-woven materials83ABE® cannula retention dressing83ABE® non-woven fabric83ABE® gauze bandage86ABE® last87ABE® last crepe88ABE® last crepe88ABE® last crepe89ABE® last crepe90Compression & Support bandages90Compression & Support bandage94ABE® extra94ABE® tubular net bandage94ABE® tubular net bandage95ABE® tubular net bandage96ABE® tubular net bandage97ABE® body bandage98ABE® body bandage99ABE® heart catheter pressure bandage90Wound care X-ray100BeeSana® Gauze swab X-ray104BeeSana® Care swab X-ray104BeeSana® Preparation balls X-ray106BeeSana® Preparation balls X-ray107BeeSana® Gauze ball X-ray108BeeSana® Gauze ball X-ray108	BeeSana <sup>®</sup> Pointed swabs	76
Alcohol swabs77Treatment care and accessories78Cotton buds78Wooden tongue depressor79Wound plasters and bandages80ABE® injection plaster80ABE® light80ABE® derm82ABE® derm plus81Fixing plasters/non-woven materials82ABE® non-woven fabric83ABE® cannula retention dressing83ABE® cannula retention dressing84BeeSana® SafeCut85ABE® last87ABE® last crepe88ABE® fix color89ABE® fix color89ABE® last crepe88ABE® lastic90Compression & Support bandages90ABE® lastic91ABE® lastic93ABE® tabular bandage94ABE® extra94ABE® extra94ABE® tubular bandage95ABE® tubular net bandage96ABE® body bandage97ABE® holder100ABE® heart catheter pressure bandage91ABE® holder100ABE® heart catheter pressure bandage91ABE® heart catheter pressure bandage92ABE® heart catheter pressure bandage90ABE® heart catheter pressure bandage91ABE® heart catheter pressure bandage91ABE® heart catheter pressure bandage92ABE® heart catheter pressure bandage90ABE® heart catheter pressure bandage100<	ABE <sup>®</sup> ribbon gauze	76
Treatment care and accessories78Cotton buds78Wooden tongue depressor79Wound plasters and bandages80ABE® injection plaster80ABE® light80ABE® derm82ABE® derm plus82Fixing plasters/non-woven materials83ABE® cannula retention dressing83ABE® injection plaster83ABE® cannula retention dressing83ABE® non-woven fabric83ABE® jake schedulage86ABE® gauze bandages, self-adhering, crepe86BeeSana® SafeCut87ABE® last crepe88ABE® last crepe88ABE® last crepe88ABE® lastic90Compression & Support bandages94ABE® lastic93ABE® lastic93ABE® lastic93ABE® lastic94ABE® ven95ABE® tubular bandage96ABE® tubular net bandage96ABE® body bandage97ABE® body bandage98ABE® body bandage99ABE® beheart catheter pressure bandage91ABE® bilb holder100ABE® bilb holder100ABE® bandage bandage98ABE® bandage bandage99ABE® bandage bandage91ABE® bandage bandage92ABE® bandage bandage96ABE® bandage call100ABE® bandage call91 <trr>ABE® bandage bandage96&lt;</trr>	BeeSana <sup>®</sup> Cellulose swabs	77
Cotton buds78Wooden tongue depressor79Wound plasters and bandages80ABE® injection plaster80ABE® injection plaster80ABE® derm81ABE® derm plus81Fixing plasters/non-woven materials83ABE® cannula retention dressing83ABE® non-woven fabric83ABE® silk84Elastic fixing bandages, self-adhering, crepe86BeeSana® SafeCut85ABE® gauze bandage86ABE® last crepe88ABE® last crepe88ABE® last adhesive90Compression & Support bandages91ABE® lastic93ABE® lastic93ABE® lastic94ABE® varia94ABE® varia94ABE® tubular bandage95ABE® tubular bandage96ABE® tubular net bandage97ABE® body bandage98ABE® body bandage98ABE® inib holder100ABE® inib holder100ABE® inib holder100ABE® pressure bandage roll101Scrotum roll101Wound care X-ray102 / 103BeeSana® Gauze swab X-ray104 / 105BeeSana® Gauze balls X-ray108BeeSana® Gauze balls X-ray108BeeSana® Gauze balls X-ray108BeeSana® Gauze ball X-ray108BeeSana® Gauze ball X-ray108BeeSana® Gauze ball X-ray108 <td>Alcohol swabs</td> <td>77</td>	Alcohol swabs	77
Wooden tongue depressor79Wound plasters and bandages80ABE® injection plaster80ABE® light80ABE® derm81ABE® derm plus82Fixing plasters/non-woven materials83ABE® cannula retention dressing83ABE® non-woven fabric83ABE® silk84Elastic fixing bandages, self-adhering, crepe85BeeSana® SafeCut85ABE® Last87ABE® Last crepe88ABE® likcic loor89ABE® listic color89ABE® last crepe88ABE® last ic90Compression & Support bandages91ABE® extra92ABE® listic93ABE® istic93ABE® tubular bandage94ABE® tubular bandage95ABE® tubular net bandage96ABE® tubular net bandage97ABE® body bandage98ABE® body bandage99ABE® heart catheter pressure bandage90ABE® inib holder100ABE® pressure bandage roll101Scrotum roll101Wound care X-ray102 / 103BeeSana® Gauze swab X-ray104BeeSana® Gauze balls X-ray108BeeSana® Gauze balls X-ray10	Treatment care and accessories	
Wound plasters and bandagesABE® injection plaster80ABE® light80ABE® light81ABE® derm82ABE® derm plus82Fixing plasters/non-woven materials83ABE® cannula retention dressing83ABE® non-woven fabric83ABE® silk84Elastic fixing bandages, self-adhering, crepe86BeeSana® SafeCut86ABE® gauze bandage86ABE® last87ABE® last crepe88ABE® last crepe88ABE® last crepe89ABE® lastic90Compression & Support bandages91ABE® lastic93ABE® lastic93ABE® varia adhesive95ABE® padded bandage94ABE® extra94ABE® extra94ABE® extra95ABE® tubular bandage96ABE® tubular bandage97ABE® loadiage coll100ABE® body bandage99ABE® heart catheter pressure bandage90ABE® inb holder100ABE® inb holder100ABE® inb holder101Scrotum roll101Wound care X-ray102 / 103BeeSana® Gauze swab X-ray104BeeSana® Gauze balls X-ray108BeeSana® Gauze		-
ABE® injection plaster80ABE® light80ABE® light80ABE® Absorbent adhesive dressing81ABE® derm82Fixing plasters/non-woven materials82Fixing plasters/non-woven materials83ABE® cannula retention dressing83ABE® non-woven fabric83ABE® silk84Elastic fixing bandages, self-adhering, crepe86BeeSana® SafeCut85ABE® gauze bandage86ABE® Last crepe88ABE® lix color89ABE® lix color89ABE® last adhesive90Compression & Support bandages91ABE® lastic93ABE® lastic93ABE® lastic93ABE® extra94ABE® extra adhesive95ABE® cilchrist's bandage96ABE® tubular bandage97ABE® body bandage98ABE® body bandage99ABE® heart catheter pressure bandage90ABE® ilim holder100ABE® pressure bandage roll101Scrotum roll101Wound care X-ray102 / 103BeeSana® Lap sponge X-ray102 / 103BeeSana® Cauze swab X-ray106BeeSana® Cauze ball X-ray108BeeSana® Cauze ball X-ray108BeeSana® Cauze ball X-ray108BeeSana® Cauze ball X-ray108BeeSana® Gauze ball X-ray108BeeSana® Cauze ball X-ray108BeeSana®		79
ABE® light80ABE® Absorbent adhesive dressing81ABE® derm82ABE® derm plus82Fixing plasters/non-woven materials83ABE® cannula retention dressing83ABE® non-woven fabric83ABE® non-woven fabric84Batsers/non-woven fabric84BeeSana® SafeCut85ABE® gauze bandage86ABE® last87ABE® last crepe88ABE® fix color89ABE® fix color89Compression & Support bandages90Compression & Support bandages91ABE® lastic93ABE® lastic93ABE® ilastic93ABE® extra94ABE® extra94ABE® extra adhesive95ABE® tubular bandage96ABE® tubular net bandage97ABE® body bandage98ABE® body bandage98ABE® heart catheter pressure bandage100ABE® lastic catheter pressure bandage100ABE® ilmb holder100ABE® ilmb holder101Scrotum roll101Wound care X-ray102 / 103BeeSana® Gauze swab X-ray102 / 103BeeSana® Cauze balls X-ray106BeeSana® Preparation balls X-ray108BeeSana® Gauze balls X-ray108BeeSana® Cauze balls X-ray108BeeSana® Cauze balls X-ray108BeeSana® Cauze balls X-ray108BeeSana® Cauze balls X-ray		
ABE® Absorbent adhesive dressing81ABE® derm82ABE® derm plus82Fixing plasters/non-woven materials83ABE® cannula retention dressing83ABE® non-woven fabric83ABE® silk84Elastic fixing bandages, self-adhering, crepe86BeeSana® SafeCut85ABE® gauze bandage86ABE® Last87ABE® Last crepe88ABE® Last crepe88ABE® Last adhesive90Compression & Support bandages91ABE® lastic93ABE® lastic93ABE® lastic93ABE® lastic94ABE® extra94ABE® extra94ABE® tubular bandage96ABE® tubular net bandage97ABE® body bandage98ABE® body bandage98ABE® heart catheter pressure bandage90ABE® pressure bandage roll101Scrotur roll101Wound care X-ray102 / 103BeeSana® Lap sponge X-ray102 / 103BeeSana® Lap sponge X-ray104 / 105BeeSana® Cauze balls X-ray108BeeSana® Cauze balls X-ray108BeeSa		-
ABE® derm82ABE® derm plus82Fixing plasters/non-woven materials83ABE® cannula retention dressing83ABE® non-woven fabric83ABE® silk84Elastic fixing bandages, self-adhering, crepe86BeeSana® SafeCut85ABE® gauze bandage86ABE® Last87ABE® Last crepe88ABE® Last crepe88ABE® Last adhesive90Compression & Support bandages91ABE® Ideal 10091ABE® lastic93ABE® lastic93ABE® added bandage94ABE® extra94ABE® extra adhesive95ABE® tubular bandage96ABE® tubular bandage97ABE® body bandage98ABE® body bandage99ABE® heart catheter pressure bandage90ABE® pressure bandage roll101Scrotum roll101Wound care X-ray102 / 103BeeSana® Gauze balls X-ray104 / 105BeeSana® Cauze balls X-ray108BeeSana® Cauze balls X-ray108BeeSana® Gauze balls X-ray108BeeSana® Cauze balls X-ray108 <td< td=""><td></td><td>-</td></td<>		-
ABE® derm plus82Fixing plasters/non-woven materials83ABE® cannula retention dressing83ABE® non-woven fabric83ABE® silk84Elastic fixing bandages, self-adhering, crepe86BeeSana® SafeCut85ABE® gauze bandage86ABE® Last87ABE® Last crepe88ABE® fix color89ABE® fix color90Compression & Support bandages90Compression & Support bandages91ABE® lastic93ABE® lastic93ABE® niversal92ABE® varia adhesive95ABE® trubular bandage94ABE® tubular bandage96ABE® body bandage97ABE® body bandage99ABE® biody bandage99ABE® heart catheter pressure bandage90ABE® limb holder100ABE® pressure bandage roll101Scrotum roll101Wound care X-ray102 / 103BeeSana® Gauze swab X-ray104 / 105BeeSana® Gauze balls X-ray106BeeSana® Cauze swab X-ray107BeeSana® Freparation balls X-ray108BeeSana® Gauze ball X-ray108BeeSana® Cauze ball X-ray1	-	-
Fixing plasters/non-woven materials83ABE® cannula retention dressing83ABE® inon-woven fabric83ABE® silk84Elastic fixing bandages, self-adhering, crepe85BeeSana® SafeCut85ABE® gauze bandage86ABE® Last87ABE® Last crepe88ABE® fix color89ABE® fix color90Compression & Support bandages90Compression & Support bandages91ABE® lastic93ABE® lastic93ABE® padded bandage94ABE® extra94ABE® extra94ABE® tubular bandage95ABE® tubular bandage96ABE® body bandage97ABE® body bandage99ABE® biody bandage99ABE® heart catheter pressure bandage900ABE® heart catheter pressure bandage1000ABE® pressure bandage roll1011Vound care X-ray102 / 103BeeSana® Gauze swab X-ray104 / 105BeeSana® Gauze swab X-ray106BeeSana® Gauze balls X-ray107BeeSana® Freparation balls X-ray108BeeSana® Gauze ball X-ray108BeeSana® Gauze ball X-ray108BeeSana® Gauze ball X-ray108BeeSana® Gauze ball X-ray108		-
ABE® cannula retention dressing83ABE® non-woven fabric83ABE® silk84Elastic fixing bandages, self-adhering, crepe85BeeSana® SafeCut85ABE® gauze bandage86ABE® Last87ABE® Last87ABE® Last crepe88ABE® fix color89ABE® Ideal 10091ABE® lastic93ABE® lastic93ABE® lastic93ABE® padded bandage94ABE® extra94ABE® extra adhesive95ABE® trubular bandage96ABE® tubular bandage96ABE® tubular bandage97ABE® fix color98ABE® body bandage98ABE® tubular net bandage99ABE® heart catheter pressure bandage90ABE® pressure bandage100ABE® pressure bandage100ABE® pressure bandage coll101Scrotum roll101Wound care X-ray102 / 103BeeSana® Gauze swab X-ray104 / 105BeeSana® Cauze swab X-ray106BeeSana® Preparation balls X-ray107BeeSana® Fonited swabs X-ray108BeeSana® Gauze ball X-ray108BeeSana® Gauze ball X-ray108BeeSana® Gauze ball X-ray108BeeSana® Gauze ball X-ray108BeeSana® Cauze ball X-ray108BeeSana® Cauze ball X-ray108	•	82
ABE® non-woven fabric83ABE® silk84Elastic fixing bandages, self-adhering, crepeBeeSana® SafeCut85ABE® gauze bandage86ABE® Last87ABE® Last87ABE® Last crepe88ABE® Itx color89ABE® Ideal 10091ABE® Ideal 10091ABE® ideal 10091ABE® ideal 10091ABE® tast93ABE® tast93ABE® tast93ABE® tast94ABE® extra94ABE® extra adhesive95ABE® tubular bandage94ABE® tubular net bandage96ABE® tubular net bandage97ABE® filmb holder100ABE® pressure bandage98ABE® pressure bandage roll101Vound care X-ray102 / 103BeeSana® Gauze swab X-ray104 / 105BeeSana® Freparation balls X-ray106BeeSana® Gauze ball X-ray108BeeSana® Gauze ball X-ray108BeeSana® Gauze ball X-ray108		00
ABE® silk84Elastic fixing bandages, self-adhering, crepeBeeSana® SafeCut85ABE® gauze bandage86ABE® Last87ABE® Last crepe88ABE® Iast color89ABE® Iast adhesive90Compression & Support bandages91ABE® Ideal 100911ABE® lastic93ABE® lastic93ABE® padded bandage94ABE® extra94ABE® extra adhesive95ABE® tubular bandage96ABE® tubular net bandage96ABE® lody bandage97ABE® lody bandage98ABE® body bandage99ABE® heart catheter pressure bandage100ABE® pressure bandage roll101Scrotum roll101Wound care X-ray102 / 103BeeSana® Gauze swab X-ray104 / 105BeeSana® Freparation balls X-ray104BeeSana® Cauze balls X-ray108BeeSana® Gauze ball X-ray108BeeSana® Gauze ball X-ray108BeeSana® Gauze ball X-ray108	3	-
Elastic fixing bandages, self-adhering, crepeBeeSana® SafeCut85ABE® gauze bandage86ABE® Last87ABE® Last crepe88ABE® tast adhesive90Compression & Support bandages91ABE® Ideal 10091ABE® lastic93ABE® lastic93ABE® added bandage94ABE® extra94ABE® extra adhesive95ABE® trubular bandage96ABE® tubular bandage96ABE® tubular net bandage97ABE® finc catheter pressure bandage98ABE® body bandage99ABE® body bandage91ABE® body bandage91ABE® filth holder100ABE® pressure bandage roll101Wound care X-ray102 / 103BeeSana® Gauze swab X-ray104 / 105BeeSana® Cauze swab X-ray104BeeSana® Cauze swab X-ray108BeeSana® Gauze balls X-ray108BeeSana® Gauze ball X-ray108BeeSana® Gauze ball X-ray108BeeSana® Gauze ball X-ray108		-
BeeSana® SafeCut85ABE® gauze bandage86ABE® Last87ABE® Last crepe88ABE® Last color89ABE® Last adhesive90Compression & Support bandages91ABE® Ideal 10091ABE® lastic93ABE® padded bandage94ABE® extra94ABE® extra94ABE® extra95ABE® tubular bandage95ABE® tubular net bandage96ABE® body bandage97ABE® body bandage98ABE® heart catheter pressure bandage90ABE® pressure bandage roll101Scrotum roll101Wound care X-ray102 / 103BeeSana® Gauze balls X-ray104 / 105BeeSana® Cauze balls X-ray108BeeSana® Gauze balls X-ray108BeeSana® Gauze ball X-ray108BeeSana® Gauze ball X-ray108		04
ABE® gauze bandage86ABE® Last87ABE® Last crepe88ABE® tast color89ABE® Last adhesive90Compression & Support bandages91ABE® ldeal 10091ABE® lastic93ABE® padded bandage94ABE® extra94ABE® extra94ABE® extra95ABE® tubular bandage95ABE® tubular net bandage96ABE® libler is bandage97ABE® body bandage98ABE® body bandage99ABE® limb holder100ABE® pressure bandage roll101Scrotum roll101Wound care X-ray102 / 103BeeSana® Gauze swab X-ray104 / 105BeeSana® Cauze swab X-ray104BeeSana® Preparation balls X-ray108BeeSana® Gauze ball X-ray108BeeSana® Gauze ball X-ray108BeeSana® Gauze ball X-ray108		85
ABE® Last87ABE® Last crepe88ABE® Iast color89ABE® Last adhesive90Compression & Support bandages91ABE® Ideal 10091ABE® lastic93ABE® padded bandage94ABE® extra94ABE® extra adhesive95ABE® tubular bandage94ABE® tubular bandage95ABE® tubular net bandage96ABE® tubular net bandage97ABE® body bandage98ABE® body bandage99ABE® heart catheter pressure bandage100ABE® pressure bandage roll101Scrotum roll101Wound care X-ray102 / 103BeeSana® Lap sponge X-ray104 / 105BeeSana® Gauze swab X-ray104BeeSana® Preparation balls X-ray108BeeSana® Gauze ball X-ray108BeeSana® Gauze ball X-ray108BeeSana® Gauze ball X-ray108		-
ABE® Last crepe88ABE® fix color89ABE® Last adhesive90Compression & Support bandages91ABE® Ideal 10091ABE® universal92ABE® padded bandage94ABE® extra94ABE® extra adhesive95ABE® tubular bandage95ABE® tubular bandage96ABE® tubular net bandage97ABE® fix catheter pressure bandage98ABE® heart catheter pressure bandage90ABE® pressure bandage roll101Scrotum roll101Wound care X-ray102 / 103BeeSana® Lap sponge X-ray102 / 103BeeSana® Gauze swab X-ray106BeeSana® Preparation balls X-ray107BeeSana® Gauze ball X-ray108BeeSana® Gauze ball X-ray108BeeSana® Gauze ball X-ray108BeeSana® Gauze ball X-ray108	5	_
ABE® fix color89ABE® Last adhesive90Compression & Support bandages91ABE® Ideal 10091ABE® universal92ABE® lastic93ABE® padded bandage94ABE® extra94ABE® extra adhesive95ABE® tubular bandage95ABE® tubular bandage96ABE® tubular net bandage97ABE® fix body bandage98ABE® body bandage99ABE® heart catheter pressure bandage100ABE® pressure bandage roll101Scrotum roll101Wound care X-ray102 / 103BeeSana® Lap sponge X-ray102 / 103BeeSana® Gauze swab X-ray106BeeSana® Preparation balls X-ray107BeeSana® Gauze ball X-ray108BeeSana® Gauze ball X-ray108BeeSana® Gauze ball X-ray108		-
ABE® Last adhesive90Compression & Support bandages91ABE® Ideal 10091ABE® universal92ABE® lastic93ABE® padded bandage94ABE® extra94ABE® extra adhesive95ABE® tubular bandage95ABE® tubular bandage96ABE® tubular net bandage97ABE® filthrist's bandage98ABE® heart catheter pressure bandage90ABE® pressure bandage roll101Scrotum roll101Wound care X-ray102 / 103BeeSana® Gauze swab X-ray104 / 105BeeSana® Preparation balls X-ray107BeeSana® Gauze ball X-ray108BeeSana® Gauze ball X-ray108BeeSana® Gauze ball X-ray108BeeSana® Gauze ball X-ray108	I	89
ABE® Ideal 10091ABE® universal92ABE® lastic93ABE® padded bandage94ABE® extra94ABE® extra adhesive95ABE® tubular bandage95ABE® tubular bandage96ABE® tubular net bandage97ABE® Gilchrist's bandage98ABE® body bandage99ABE® heart catheter pressure bandage100ABE® pressure bandage roll101Scrotum roll101Wound care X-ray102 / 103BeeSana® Gauze swab X-ray104 / 105BeeSana® Freparation balls X-ray107BeeSana® Pointed swabs X-ray108BeeSana® Gauze ball X-ray108BeeSana® Gauze ball X-ray108BeeSana® Gauze ball X-ray108	ABE <sup>®</sup> Last adhesive	-
ABE® Ideal 10091ABE® universal92ABE® lastic93ABE® padded bandage94ABE® extra94ABE® extra adhesive95ABE® tubular bandage95ABE® tubular bandage96ABE® tubular net bandage97ABE® Gilchrist's bandage98ABE® body bandage99ABE® heart catheter pressure bandage100ABE® pressure bandage roll101Scrotum roll101Wound care X-ray102 / 103BeeSana® Gauze swab X-ray104 / 105BeeSana® Freparation balls X-ray107BeeSana® Pointed swabs X-ray108BeeSana® Gauze ball X-ray108BeeSana® Gauze ball X-ray108BeeSana® Gauze ball X-ray108	Compression & Support bandages	
ABE® lastic93ABE® padded bandage94ABE® extra94ABE® extra adhesive95ABE® tra adhesive95ABE® tubular bandage96ABE® tubular net bandage96ABE® Gilchrist's bandage97ABE® Gilchrist's bandage98ABE® heart catheter pressure bandage90ABE® pressure bandage roll100ABE® pressure bandage roll101Scrotum roll101Wound care X-ray102 / 103BeeSana® Gauze swab X-ray104 / 105BeeSana® Gauze balls X-ray106BeeSana® Preparation balls X-ray108BeeSana® Gauze ball X-ray108BeeSana® Gauze ball X-ray108		91
ABE® padded bandage94ABE® extra94ABE® extra adhesive95ABE® Ven95ABE® tubular bandage96ABE® tubular net bandage97ABE® Gilchrist's bandage98ABE® body bandage99ABE® heart catheter pressure bandage100ABE® pressure bandage roll101Scrotum roll101Wound care X-ray102 / 103BeeSana® Gauze swab X-ray104 / 105BeeSana® Gauze balls X-ray107BeeSana® Preparation balls X-ray108BeeSana® Gauze ball X-ray108BeeSana® Gauze ball X-ray108	ABE <sup>®</sup> universal	92
ABE® extra94ABE® extra adhesive95ABE® tva adhesive95ABE® tva adhesive96ABE® tubular bandage96ABE® tubular net bandage97ABE® Gilchrist's bandage98ABE® body bandage99ABE® heart catheter pressure bandage100ABE® pressure bandage roll101Scrotum roll101Wound care X-ray102 / 103BeeSana® Lap sponge X-ray104 / 105BeeSana® Gauze swab X-ray106BeeSana® Preparation balls X-ray107BeeSana® Pointed swabs X-ray108BeeSana® Gauze ball X-ray108BeeSana® Gauze ball X-ray108	ABE <sup>®</sup> lastic	93
ABE® extra adhesive95ABE® Ven95ABE® tubular bandage96ABE® tubular net bandage97ABE® Gilchrist's bandage98ABE® body bandage99ABE® heart catheter pressure bandage100ABE® pressure bandage roll101Scrotum roll101Wound care X-ray102 / 103BeeSana® Gauze swab X-ray104 / 105BeeSana® Gauze balls X-ray106BeeSana® Preparation balls X-ray108BeeSana® Gauze ball X-ray108BeeSana® Gauze ball X-ray109	ABE <sup>®</sup> padded bandage	94
ABE® Ven95ABE® tubular bandage96ABE® tubular net bandage97ABE® Gilchrist's bandage98ABE® body bandage99ABE® heart catheter pressure bandage100ABE® pressure bandage roll101Scrotum roll101Wound care X-ray102 / 103BeeSana® Lap sponge X-ray102 / 103BeeSana® Gauze swab X-ray104 / 105BeeSana® Preparation balls X-ray107BeeSana® Pointed swabs X-ray108BeeSana® Gauze ball X-ray108	ABE <sup>®</sup> extra	94
ABE® tubular bandage96ABE® tubular net bandage97ABE® Gilchrist's bandage98ABE® body bandage99ABE® heart catheter pressure bandage100ABE® limb holder100ABE® pressure bandage roll101Scrotum roll101Wound care X-ray102 / 103BeeSana® Lap sponge X-ray104 / 105BeeSana® Gauze swab X-ray106BeeSana® Preparation balls X-ray107BeeSana® Pointed swabs X-ray108BeeSana® Gauze ball X-ray108		95
ABE® tubular net bandage97ABE® Gilchrist's bandage98ABE® body bandage99ABE® heart catheter pressure bandage100ABE® limb holder100ABE® pressure bandage roll101Scrotum roll101Wound care X-ray102 / 103BeeSana® Gauze swab X-ray104 / 105BeeSana® Gauze balls X-ray106BeeSana® Preparation balls X-ray107BeeSana® Gauze ball X-ray108BeeSana® Gauze ball X-ray108	ABE® Ven	95
ABE® Gilchrist's bandage98ABE® body bandage99ABE® heart catheter pressure bandage100ABE® limb holder100ABE® pressure bandage roll101Scrotum roll101Wound care X-ray102 / 103BeeSana® Gauze swab X-ray104 / 105BeeSana® Gauze balls X-ray106BeeSana® Preparation balls X-ray107BeeSana® Gauze ball X-ray108BeeSana® Gauze ball X-ray108BeeSana® Gauze ball X-ray109	ABE <sup>®</sup> tubular bandage	96
ABE® body bandage99ABE® heart catheter pressure bandage100ABE® limb holder100ABE® pressure bandage roll101Scrotum roll101Wound care X-ray102 / 103BeeSana® Gauze swab X-ray104 / 105BeeSana® Gauze balls X-ray106BeeSana® Preparation balls X-ray107BeeSana® Gauze balls X-ray108BeeSana® Gauze ball X-ray108		97
ABE® heart catheter pressure bandage100ABE® limb holder100ABE® pressure bandage roll101Scrotum roll101Wound care X-ray102 / 103BeeSana® Gauze swab X-ray104 / 105BeeSana® Gauze balls X-ray106BeeSana® Preparation balls X-ray107BeeSana® Gauze balls X-ray108BeeSana® Gauze ball X-ray108BeeSana® Gauze ball X-ray109	ABE® Gilchrist's bandage	98
ABE® limb holder100ABE® pressure bandage roll101Scrotum roll101Wound care X-ray102 / 103BeeSana® Lap sponge X-ray102 / 103BeeSana® Gauze swab X-ray104 / 105BeeSana® Gauze balls X-ray106BeeSana® Preparation balls X-ray107BeeSana® Pointed swabs X-ray108BeeSana® Gauze ball X-ray109		99
ABE® pressure bandage roll101Scrotum roll101Wound care X-ray102 / 103BeeSana® Lap sponge X-ray102 / 103BeeSana® Gauze swab X-ray104 / 105BeeSana® Gauze balls X-ray106BeeSana® Preparation balls X-ray107BeeSana® Pointed swabs X-ray108BeeSana® Gauze ball X-ray109	ABE® heart catheter pressure bandage	100
Scrotum roll101Wound care X-ray102 / 103BeeSana® Lap sponge X-ray102 / 103BeeSana® Gauze swab X-ray104 / 105BeeSana® Gauze balls X-ray106BeeSana® Preparation balls X-ray107BeeSana® Pointed swabs X-ray108BeeSana® Gauze ball X-ray109	ABE <sup>®</sup> limb holder	100
Wound care X-ray102 / 103BeeSana® Lap sponge X-ray102 / 103BeeSana® Gauze swab X-ray104 / 105BeeSana® Gauze balls X-ray106BeeSana® Preparation balls X-ray107BeeSana® Pointed swabs X-ray108BeeSana® Gauze ball X-ray109	ABE <sup>®</sup> pressure bandage roll	101
BeeSana® Lap sponge X-ray102 / 103BeeSana® Gauze swab X-ray104 / 105BeeSana® Gauze balls X-ray106BeeSana® Preparation balls X-ray107BeeSana® Pointed swabs X-ray108BeeSana® Gauze ball X-ray109		101
BeeSana® Gauze swab X-ray104 / 105BeeSana® Gauze balls X-ray106BeeSana® Preparation balls X-ray107BeeSana® Pointed swabs X-ray108BeeSana® Gauze ball X-ray109	Wound care X-ray	
BeeSana® Gauze balls X-ray106BeeSana® Preparation balls X-ray107BeeSana® Pointed swabs X-ray108BeeSana® Gauze ball X-ray109	BeeSana <sup>®</sup> Lap sponge X-ray	102 / 103
BeeSana® Preparation balls X-ray107BeeSana® Pointed swabs X-ray108BeeSana® Gauze ball X-ray109		104 / 105
BeeSana® Pointed swabs X-ray108BeeSana® Gauze ball X-ray109		-
BeeSana® Gauze ball X-ray 109		-
	,	-
ABE <sup>™</sup> ribbon gauze X-ray 109		-
	ABE™ ribbon gauze X-ray	109

#### Set solutions

BeeSana <sup>®</sup> Ward catheter set	112
BeeSana <sup>®</sup> PEG bandage set	113
BeeSana <sup>®</sup> Bandage change set	114
BeeSana <sup>®</sup> Suture removal set	114
BeeSana <sup>®</sup> Ward protection set	115

#### Incontinence products

•	
Incontinence briefs	
Superslip®	118
Superslip <sup>®</sup> Air Comfort	119
BeeSana <sup>®</sup> Medi-Pants	119
Incontinence liners/pads	
Liberty Light <sup>®</sup>	120
Incontinence pads	
Superform <sup>®</sup>	120
Superform <sup>®</sup> Air	121
BeeSana <sup>®</sup> Absorbent cushion	121
BeeSana <sup>®</sup> Absorbent pad	122
BeeSana <sup>®</sup> Rectangular pad	122
BeeSana <sup>®</sup> Male incontinence pad	123
BeeSana <sup>®</sup> Female incontinence pad	123
Fixation briefs	
Light fixation briefs	124
Panty fixation briefs	125
Incontinence undersheets	
Layered undersheet	126
BeeSana <sup>®</sup> Comflock Super	127
Nursing pads	128

#### Disinfection / care

Ethasept® Ethasept® Gel	132 133 134
Ethasept <sup>®</sup> Gel	
	134
Medizid <sup>®</sup> Rapid+	
Medizid <sup>®</sup> Rapid+ / Rapid+ XXL disinfectant wipes	135
Medizid <sup>®</sup> Rapid QF Flowpack	136
Medizid® AF+	137
Medizid® Alpha+	138
Medizid <sup>®</sup> Alpha+ XXL disinfectant wipes	139
	140
Instrusol® AF+ 142 /	143
Care	
Gentle Med <sup>®</sup> wash lotion	145
Gentle Med <sup>®</sup> skin care lotion (O/W)	146
Gentle Med <sup>®</sup> cream lotion (W/O)	146
Gentle Med <sup>®</sup> active gel	147
	147
Gentle Med <sup>®</sup> skin cleansing and care foam	148
Gentle Med <sup>®</sup> moist care wipes	148
Application aids	
Metering pumps/attachments	149
	149
Canister tap/canister key	149
Dosing bottle/graduated flask	149
Disinfectant dispenser 150 /	151

#### Hygiene / cleaning / cleaning products

beopur WC fresh	154
MEDIzid <sup>®</sup> sanitary cleaning	155
Hygiene in the sanitary sector	156 / 157
Hygiene in patient wards	158 - 161
Cleaning in commercial kitchens	162 - 167

#### Ward and patient supplies

Patient and personal hygiene	
BagBath <sup>®</sup> Haircare	172
BagBath®	172
BeeSana <sup>®</sup> Wash glove	174
Protective towel	174
Protective bib	175
BeeSana® Disposable briefs	175
Anti-slip sock universal	176
Surgical sock	176
PE shoe covers	177
Cellulose sponge sheets	177
Layered cellulose	177
Rollicel <sup>®</sup> / medical crepe	178
Protective sheets	178
Protective covering	179
Multi-purpose cloths	179
Protective bed coverings	179
Universal cloth	180
Wipe cloths	180
Cosmetic wipes	181
V-fold paper hand towels	181
Sick bag (BaSick Bag)	182
Kidney dish	182
Disposable comb	182
Disposable razor	183
Urine bottle	183
Urine bottle holder	183
Diagnosis and accessories	
Nurse's stethoscope	184
Tourniquet	184
Dosicair®	185
Urine beaker	185
Ultrasound scan gel	185
Ward supplies and accessories	
Cold-Warm-Swab	186
Non-woven cover for pads	186
Rescue blanket	187
Medical beaker	187
Medication dispenser	188
Medi-Box	188
Wall holder	189
Rubbish bags	190
Kitchen towel	191
Toilet paper	191

## We know gloves like the back of our hand

When making our examination and surgical gloves, we at Meditrade understand that it is vital to provide gloves that combine mechanical strength with good touch sensitivity. Above all, our first priority is to protect the user and the patient, in accordance with regulatory requirements. We also invest a lot in product development know-how and in our production sites. Our huge range of latex glove models are complimented with extremely cost-effective gloves in vinyl and other alternative materials, including gloves from pure natural rubber and nitrile gloves made from nitrile butadiene rubber, ensuring we can always offer the right product for the user.

# **Reference**<sup>™</sup>

### Our powdered latex examination glove

#### Features

- Latex examination glove, natural (DIN EN 455 standard)
- lightly powdered
- micro-roughened
- elastic and tear-proof
- ambidextrous
- suitable for use as a food contact material



## 

### Non-sterile

Size	REF	Pack contents	Shipping unit	PZN
XS	1121XS	100 pcs	10 packs	06177515
S	11215	100 pcs	10 packs	06177521
М	1121M	100 pcs	10 packs	06177538
L	1121L	100 pcs	10 packs	06177544
XL	1121XL	100 pcs	10 packs	02243899

#### Usage

In all areas where the user is exposed to a risk of contamination. Suitable for highly infectious situations in hospitals and in general practice. Suited to use in the manufacturing sector.

# Gentle Skin<sup>®</sup> sensitive

## Our most popular latex examination glove

#### Features

- Latex examination glove, natural (DIN EN 455 standard)
- PPE Cat III
- powder-free, low protein
- micro-roughened
- elastic and tear-proof
- ambidextrous
- suitable for use as a food contact material





#### Non-sterile

Size	REF	Pack contents	Shipping unit	PZN
XS	1221RXST	100 pcs	10 packs	02243563
S	1221RST	100 pcs	10 packs	02243586
М	1221RMT	100 pcs	10 packs	02243592
L	1221RLT	100 pcs	10 packs	02243600
XL	1221RXLT	100 pcs	10 packs	02243617

#### Usage

In all areas where the user is exposed to a risk of contamination. Suitable for highly infectious situations in hospitals and in general practice. Suitable for use in nursing, by lab personnel, in the pharmaceutical industry, in cleaning work and in manufacturing.

# Gentle Skin<sup>®</sup> classic / classic x-long

## Our tough sterile latex examination glove

#### Features

- Latex examination glove, natural (DIN EN 455 standard)
- PPE Cat III
- enhanced wall thickness with 20% more latex compared with conventional latex gloves
- powder-free, low protein
- micro-roughened
- elastic and tear-proof
- ambidextrous
- suitable for use as a food contact material
- classic x-long: with long cuff (min. 29 cm)





#### Usage

In all areas where the user is exposed to a risk of contamination. Also for emergency situations, rescue services and ICU staff. Suitable for highly infectious situations in hospitals and in general practice. Usable in nursing, by lab personnel, in the pharmaceutical industry, in cleaning work and in manufacturing.

Gentle Skin <sup>®</sup> classic, non-sterile				
Size	REF	Pack contents	Shipping unit	PZN
XS	1221RXS	100 pcs	10 packs	02243698
S	1221RS	100 pcs	10 packs	02243706
М	1221RM	100 pcs	10 packs	02243729
L	1221RL	100 pcs	10 packs	02243735
XL	1221RXL	100 pcs	10 packs	02243801

Gentle Skin <sup>®</sup> classic x-long, non-sterile				
S	12235	100 pcs	10 packs	11349906
М	1223M	100 pcs	10 packs	11349935
L	1223L	100 pcs	10 packs	11349958

# Gentle Skin<sup>®</sup> grip

## Our easy-grip latex glove

#### Features

- Latex examination glove, natural (DIN EN 455 standard)
- PPE Cat III
- good, secure grip
- powder-free, low-protein
- micro-roughened
- elastic and tear-proof
- ambidextrous
- suitable for use as a food contact material



## Non-sterile

Size	REF	Pack contents	Shipping unit	PZN
XS	1221GRIP-XS	100 pcs	10 packs	02470253
S	1221GRIP-S	100 pcs	10 packs	02470282
М	1221GRIP-M	100 pcs	10 packs	02470299
L	1221GRIP-L	100 pcs	10 packs	02470307
XL	1221GRIP-XL	100 pcs	10 packs	02470313

Gentle Skin grip

#### Usage

Suitable for use in dental surgery: its easy-grip qualities enable the safe, non-slip handling of instruments and machinery.

# Gentle Skin<sup>®</sup> compact +

## Our latex glove for improved hygiene

#### Features

- Latex examination glove, natural (DIN EN 455 standard)
- our popular disposable glove for manufacturing
- good, secure grip
- powder-free
- micro-roughened
- elastic and tear-proof
- ambidextrous
- suitable for use as a food contact material



#### Usage

Usable in manufacturing and for food processing applications. Ideal for household and cleaning work.





### Non-sterile

Size	REF	Pack contents	Shipping unit	PZN
XS	1221XSI	100 pcs	10 packs	05522453
S	1221SI	100 pcs	10 packs	05522447
М	1221MI	100 pcs	10 packs	05522430
L	1221LI	100 pcs	10 packs	05522424
XL	1221XLI	100 pcs	10 packs	05521198

11

# **Gentle Tec**<sup>®</sup>

## Our long latex glove in a polybag

#### Features

- Latex examination glove, natural (DIN EN 455 standard)
- PPE Cat III
- with long cuff (min. 29 cm)
- enhanced wall thickness with 20% more latex compared with conventional latex gloves
- packed in a practical dispensing pouch for additional protection against soiling
- powder-free, low protein
- micro-roughened
- elastic and tear-proof
- ambidextrous



## Usage

In all areas where the user is exposed to a risk of contamination. Also for emergency situations, rescue services and ICU staff. Suitable for highly infectious situations in hospitals and in general practice. Usable in nursing, by lab personnel, in the pharmaceutical industry, in cleaning work and in manufacturing.

### Non-sterile

Size	REF	Pack contents	Shipping unit	PZN
S	12295	100 pcs	10 packs	06966643
М	1229M	100 pcs	10 packs	06966666
L	1229L	100 pcs	10 packs	06966672



# Gentle Skin<sup>®</sup> Aloecare<sup>™</sup>

## Our latex glove with an Aloe vera coating

#### Features

- Latex examination glove, green (DIN EN 455 standard)
- inner coating with Aloe vera
- powder-free, low protein
- micro-roughened
- elastic and tear-proof
- ambidextrous



#### Usage

LATEX

Latex glove for applications such as cosmetics. The special inner coating with Aloe vera offers benefits for users with stressed skin.

Non-stei				
Size	REF	Pack contents	Shipping unit	PZN
XS	1231XS	100 pcs	10 packs	05030230
S	12315	100 pcs	10 packs	02244516
М	1231M	100 pcs	10 packs	02244580
L	1231L	100 pcs	10 packs	02244597
XL	1231XL	100 pcs	10 packs	05030247

# Gentle Skin<sup>®</sup> black

## Our black latex glove

#### Features

- Latex examination glove, black (DIN EN 455 standard)
- powder-free, low protein
- micro-roughened
- elastic and tear-proof
- ambidextrous
- suitable for use as a food contact material



#### Usage

The latex glove for hair stylists, beauticians and tattoo/piercing artists. Usable by caterers and chefs, and for general applications in the food industry and cleaning.

Non-ste				
Size	REF	Pack contents	Shipping unit	PZN
XS	1224XS	100 pcs	10 packs	11352676
S	1224S	100 pcs	10 packs	11352682
М	1224M	100 pcs	10 packs	11352736
L	1224L	100 pcs	10 packs	11352771
XL	1224XL	100 pcs	10 packs	11352802

# Nitril<sup>®</sup> 3000

## The white nitrile examination and protective glove

#### Features

- Nitrile examination glove, white (DIN EN 455 standard)
- PPE Cat III
- powder-, protein- and latex-free
- micro-roughened
- elastic and tear-proof
- ambidextrous, with rolled edge
- suitable for use as a food contact material



### Non-sterile

Size	REF	Pack contents	Shipping unit	PZN
XS	1280XS	100 pcs	10 packs	06966703
S	12805	100 pcs	10 packs	06967967
М	1280M	100 pcs	10 packs	06968056
L	1280L	100 pcs	10 packs	06968300
XL	1280XL	100 pcs	10 packs	06968381

#### Usage

Suitable for handling cleaning products and disinfectants, and for other applications in laboratories, the pharmaceutical industry, for professional cleaning, in the catering business and for handling foodstuffs in general. Suitable for users with a diagnosed (type I) latex allergy.

# Nitril<sup>®</sup> NextGen<sup>®</sup>

## The stretchy nitrile examination and protective glove

#### Features

- Nitrile examination glove, light blue (DIN EN 455 standard)
- PPE Cat III
- powder-, protein- and latex-free
- micro-roughened
- elastic and tear-proof
- ambidextrous, with rolled edge
- suitable for use as a food contact material





#### Non-sterile

Size	REF	Pack contents	Shipping unit	PZN		
XS	1283XS	100 pcs	10 packs	09102239		
S	12835	100 pcs	10 packs	09102245		
М	1283M	100 pcs	10 packs	09102251		
L	1283L	100 pcs	10 packs	09102268		
XL	1283XL	100 pcs	10 packs	09102274		

#### Usage

Suitable for handling cleaning products and disinfectants, and for other applications in laboratories, the pharmaceutical industry, for professional cleaning, in the catering business and for handling foodstuffs in general. Suitable for users with a diagnosed (type I) latex allergy.

# Nitril<sup>®</sup> Sensory<sup>®</sup>

## Our touch-sensitive nitrile examination glove

#### Features

- Nitrile examination glove, blue, white and violet (DIN EN 455 standard)
- PPE Cat III
- powder-, protein- and latex-free
- micro-roughened
- elastic and tear-proof
- ambidextrous, with rolled edge
- suitable for use as a food contact material

#### Usage

Suitable for handling cleaning products and disinfectants, and for other applications in laboratories, the pharmaceutical industry, for professional cleaning, in the catering business and for handling foodstuffs in general. Suitable for users with a diagnosed (type I) latex allergy.



### Non-sterile

	Size	REF	Pack contents	Shipping unit	PZN
Blau	XS	2283XS	200 pcs	10 packs	14310075
	S	22835	200 pcs	10 packs	14272305
	М	2283M	200 pcs	10 packs	14272311
	L	2283L	200 pcs	10 packs	14272328
	XL	2283XL	200 pcs	10 packs	14272334
Weiß	XS	2280XS	200 pcs	10 packs	14306783
	S	2280S	200 pcs	10 packs	14272268
	М	2280M	200 pcs	10 packs	14272274
	L	2280L	200 pcs	10 packs	14272280
	XL	2280XL	200 pcs	10 packs	14272297
Violet	XS	2285XS	200 pcs	10 packs	-
	S	22855	200 pcs	10 packs	-
	М	2285M	200 pcs	10 packs	-
	L	2285L	200 pcs	10 packs	-
	XL	2285XL	200 pcs	10 packs	-

# Nitril<sup>®</sup> BestGen<sup>®</sup>

## The elastic nitrile examination and protective glove

#### Features

- Nitrile examination glove, blue (DIN EN 455 standard)
- PPE Cat III
- powder-, protein- and latex-free
- micro-roughened
- elastic and tear-proof
- ambidextrous, with rolled edge
- suitable for use as a food contact material



#### Usage

Suitable for handling cleaning products and disinfectants, and for other applications in laboratories, the pharmaceutical industry, for professional cleaning, in the catering business and for handling foodstuffs in general. Suitable for users with a diagnosed (type I) latex allergy.

Non-ste				
Size	REF	Pack contents	Shipping unit	PZN
XS	1286XS	100 pcs	10 packs	10134612
S	12865	100 pcs	10 packs	10134641
М	1286M	100 pcs	10 packs	10129545
L	1286L	100 pcs	10 packs	10134658
XL	1286XL	100 pcs	10 packs	10134670
XXL	1286XXL	100 pcs	10 packs	11349987

# Nitril<sup>®</sup> 3000 x-long blue

## Our long nitrile glove

#### Features

- Nitrile examination glove, blue (DIN EN 455 standard)
- PPE Cat III
- with long cuff (approx. 29 cm) for better protection
- powder-, protein- and latex-free
- micro-roughened
- elastic and tear-proof
- ambidextrous
- suitable for use as a food contact material



#### Usage

Thanks to its long cuff, this glove offers effective (splash) protection both to the hands and to parts of the forearms when handling liquids. This makes it suited to the handling of chemicals, cytostatic agents, detergents and disinfectants. The glove can also be used in the lab, pharmaceuticals and catering. Also for latex allergy sufferers with a diagnosed type I allergy.

Non-ste				
Size	REF	Pack contents	Shipping unit	PZN
S	12825	100 pcs	10 packs	00420825
М	1282M	100 pcs	10 packs	00420802
L	1282L	100 pcs	10 packs	00400857
XL	1282XL	100 pcs	10 packs	00400834

# Nitril<sup>®</sup> black

## Our black nitrile examination and protective glove

#### Features

- Nitrile examination glove, black (DIN EN 455 standard)
- PPE Cat III
- powder-, protein- and latex-free
- micro-roughened
- elastic and tear-proof
- ambidextrous, with rolled edge
- suitable for use as a food contact material





#### Usage

Suitable for handling cleaning products and disinfectants, and for other applications in laboratories, the pharmaceutical industry, for professional cleaning, in the catering business and for handling foodstuffs in general. Suitable for users with a diagnosed (type I) latex allergy.

071

## Non-sterile

	Size	KEF	Pack contents	Snipping unit	PZN
	S	12845	100 pcs	10 packs	12614730
_	Μ	1284M	100 pcs	10 packs	12614747
	L	1284L	100 pcs	10 packs	12614753
	XL	1284XL	100 pcs	10 packs	12614782

# Nitril<sup>®</sup> Viola

## The stretchy nitrile examination and protective glove

#### Features

- Nitrile examination glove, violet (DIN EN 455 standard)
- PPE Cat III
- powder-, protein- and latex-free
- with silk-sericin coating
- micro-roughened
- elastic and tear-proof
- ambidextrous, with rolled edge
- suitable for use as a food contact material



#### Usage

Suitable for handling cleaning products and disinfectants, and for other applications in laboratories, the pharmaceutical industry, for professional cleaning, in the catering business and for handling foodstuffs in general. Suitable for users with a diagnosed (type I) latex allergy.

## Non-sterile

Size	REF	Pack contents	Shipping unit	PZN
XS	1285XS	100 pcs	10 packs	-
S	12855	100 pcs	10 packs	-
М	1285M	100 pcs	10 packs	-
L	1285L	100 pcs	10 packs	-
XL	12855	100 pcs	10 packs	-

# Nitril<sup>®</sup> Magenta

## Our pink nitrile glove

#### Features

- Nitrile examination glove, pink (DIN EN 455 standard)
- PPE Cat III
- powder-, protein- and latex-free
- with silk-sericin coating
- micro-roughened
- elastic and tear-proof
- ambidextrous, with rolled edge
- suitable for use as a food contact material

#### Usage

A favourite in the cosmetics industry. Suitable for handling cleaning products and disinfectants, and for applications in laboratories, in the catering business and for handling foodstuffs in general. The silk-sericin inner coating moisturises dry and sensitive skin.

Non-ste	rile			
Size	REF	Pack contents	Shipping unit	PZN
S	12875	100 pcs	10 packs	11584731
М	1287M	100 pcs	10 packs	11584694
L	1287L	100 pcs	10 packs	11584783





# Nitril LIOX<sup>®</sup>

## Our first antimicrobial examination glove

#### Features

- Nitrile examination glove, indigo (DIN EN 455 standard)
- PPE Cat III
- powder-, protein- and latex-free
- micro-roughened
- elastic and tear-proof
- ambidextrous, with rolled edge
- suitable for use as a food contact material



#### Usage

An examination and protective glove for universal use, for applications in any areas where the user is exposed to a risk of contamination. Suitable for emergency situations, rescue services and ICU staff. Suitable for highly infectious situations in hospitals and in general practice.

#### **Benefits:**

- Protects against a wide variety of hospital-acquired infections
- Kills up to 99.999% of selected microorganisms
- Skin compatibility dermatologically tested

#### Mechanism of action

An active dye incorporated into the outer surface of the glove material converts conventional oxygen into a compound known as singlet oxygen <sup>1</sup>O<sub>2</sub> as soon as the glove material comes into contact with both light and oxygen. The highly reactive singlet oxygen attacks the cell membranes of gram-positive bacteria (including MRSA, VRE, etc.), which ultimately leads to a dying-off of the affected microbes. Singlet oxygen is already in widespread use within medicine.



Size	REF	Pack contents	Shipping unit	PZN
XS	2284XS	200 pcs	10 packs	15377756
S	22845	200 pcs	10 packs	15377762
М	2284M	200 pcs	10 packs	15377779
L	2284L	200 pcs	10 packs	15377785
XL	2284XL	180 pcs	10 packs	15377816

# Vinyl 2000

## Our vinyl examination glove

#### Features

- latex-free examination glove (DIN EN 455 standard)
- with rolled edge
- ambidextrous
- vinyl is suitable for use as a latex substitute for users and patients with a diagnosed (type I) latex allergy.



#### Usage

Suitable for daily patient care in retirement homes and hospitals, in general practice, for general healthcare work, and for handling foodstuffs that do not contain fat. Also suitable for use in the manufacturing sector. Thanks to its features and cost-effectiveness, it is ideal for applications with a very high consumption of disposable gloves.

### Powdered, non-sterile

Size	REF	Pack contents	Shipping unit	PZN
XS	1151XS	100 pcs	10 packs	02244002
S	11515	100 pcs	10 packs	06177691
М	1151M	100 pcs	10 packs	06177716
L	1151L	100 pcs	10 packs	06177722
XL	1151XL	100 pcs	10 packs	06177739

## Powder-free, non-sterile

XS	1251XS	100 pcs	10 packs	02243907
S	12515	100 pcs	10 packs	02243913
М	1251M	100 pcs	10 packs	02243936
L	1251L	100 pcs	10 packs	02243942
XL	1251XL	100 pcs	10 packs	02243994

# Vinyl Blue Gen PF

## Our blue vinyl examination glove

#### Features

- latex-free examination glove (DIN EN 455 standard)
- beaded Cuff
- ambidextrous
- vinyl is suitable for use as a latex substitute for users and patients with a diagnosed (type I) latex allergy

#### Usage

Suitable for daily patient care in retirement homes and hospitals, in general practice, for general healthcare work, and for handling foodstuffs that do not contain fat. Also suitable for use in the manufacturing sector. Thanks to its features and cost-effectiveness, it is ideal for applications with a very high consumption of disposable gloves.

### Powder-free, non-sterile

Size	REF	Pack contents	Shipping unit	PZN
XS	1255XS	100 pcs	10 packs	-
S	12555	100 pcs	10 packs	-
М	1255M	100 pcs	10 packs	-
L	1255L	100 pcs	10 packs	-
XL	12535L	100 pcs	10 packs	-

# Vinyl 2000 Stretch PF

## Our highly elastic vinyl examination glove

#### Features

- latex-free examination glove (DIN EN 455 standard)
- powder-free
- more stretchy and also denser when compared with conventional vinyl gloves
- with rolled edge

- ambidextrous
- vinyl is suitable for use as a latex substitute for users and patients with a diagnosed (type I) latex allergy.

# Viny 2009 Strategy and the second se

#### Usage

The inexpensive, latex-free and comfortable alternative for care work. Suitable for use in the manufacturing sector. Thanks to its features and cost-effectiveness, it is ideal for applications with a very high consumption of disposable gloves.

Non-ste	rile			
Size	REF	Pack contents	Shipping unit	PZN
S	12535	100 pcs	10 packs	09544664
М	1253M	100 pcs	10 packs	09544670
L	1253L	100 pcs	10 packs	09544687
XL	1253XL	100 pcs	10 packs	09544724





# Yarn glove

Our soft glove from pure cotton

#### Features

- pure cotton glove
- air permeable
- anatomical fit



#### Usage

For use by skincare practitioners. Especially when treating a patient with a very greasy or slow-absorbing ointment, for example for atopic dermatitis, for very dry skin or eczema and to cover gauze dressings (avoid direct contact with open wounds!). Also usable for the reprocessing of sterilised surgical instruments. Not a medical device.

Non-sterile				
Size	REF	Pack contents	Shipping unit	PZN
Size 6/XS	3611	10 pairs	30 packs	04354686
Size 7/XS	3612	10 pairs	30 packs	04354692
Size 8/S	3613	10 pairs	30 packs	04354700
Size 9/S	3614	10 pairs	30 packs	04354717
Size 10/M	3619	10 pairs	30 packs	04354568
Size 11/M	3620	10 pairs	30 packs	04354574
Size 12/L	3621	10 pairs	30 packs	04354580
Size 13/L	3622	10 pairs	30 packs	04354597
Size 14/L	3623	10 pairs	30 packs	04354657
Size 15/L	3624	10 pairs	30 packs	04354663

## **Froggy**<sup>™</sup> Our hygienic glove made from soft polyethylene

#### Features

- transparent PE hygienic glove
- powder-, protein- and latex-free
- in a handy glove dispenser for easy removal
- welded seams
- ambidextrous
- in women's and men's sizes

		M	editrade'	
Meditrade		Froggy	, made	
		für den Einmalgebrauch	E	
Froggy Hydionelianodischuhe aus PE tür den Einemalgebrauch		Damen gehämmern, sehr reißfest		
für den Exerces Herren gehämmert, sehr reil/kest	REJ 11025		2	
gehammer -	2			

Non-ster	ile			
Size	REF	Pack contents	Shipping unit	PZN
Women	11925	100 pcs	100 packs	06177745
Men	1192L	100 pcs	100 packs	06177751

#### Usage

Economical to use for simple, day-to-day activities. Suitable for protecting the hands from general environmental dirt and soiling, and applications needing to meet hygienic standards. Usable in chemists, supermarkets, bakeries, the food industry and service stations. <u>IMPORTANT</u>: Use only once when filling motor vehicles with diesel fuel. Short-term protection against slight, superficial diesel soiling: this PE disposable glove is not classed as personal protective equipment (PPE) in accordance with Regulation (EU) 2016/425. Not a medical device.

# Latex finger cots

Our non-sterile disposable product with rolled edge and a smooth surface

#### Features

- approved as safe for use as a short-term food contact material
- white natural latex finger cots, powder-, silicone- and thiuram-free
- only suitable for a minimal level of risk
- not suitable for activities involving hazardous, sharp-edged or aggressive materials



#### Usage

For covering individual fingers to insert pessaries/suppositories, apply creams or conduct simple examinations. Typically used in the healthcare, hygiene and manufacturing sectors.

## Non-sterile

Size	REF	Pack contents	Shipping unit	PZN
S	9776	100 pcs	1 pack	01333784
М	9777	100 pcs	1 pack	01333790
L	9778	100 pcs	1 pack	01333809
XL	9779	100 pcs	1 pack	01333815



# MyClean<sup>®</sup> easy touch

The reusable glove made from latex

#### Features

- reusable glove made from latex, powder-free
- fully anatomical fit
- comfortable to wear thanks to cotton flocked lining
- special, easy-grip surface with honeycomb texture on glove palms
- suitable for food contact use







#### Usage

Especially suitable for cleaning work, and for rinsing, washing, garden work and animal care. As a food contact material, the glove can also be used in the hospitality, catering and food processing industries.

Non-steril	e		
Size	REF	Pack contents	Shipping unit
S	MAM74152	1 pair	100 packs
М	MAM74153	1 pair	100 packs
L	MAM74154	1 pair	100 packs
XL	MAM76286	1 pair	100 packs

## **Copolymed**<sup>®</sup> Our sterile copolymer examination glove

#### Features

- sterile copolymer examination glove (DIN EN 455 standard)
- powder-, protein- and latex-free
- touch-sensitive
- tear-proof
- ambidextrous
- pressed on paper



#### Usage

Suitable for short-term medical interventions, for sterile activities in the operating theatre and in wards, as well as in manufacturing. Can be used either as an underglove or an overglove if additional protection is required against soiling and fluids.

Sterile, i	ndividually pa	cked		
Size	REF	Pack contents	Shipping unit	PZN
S	80925	100 pcs	12 packs	03896416
М	8092M	100 pcs	12 packs	03896221
L	8092L	100 pcs	12 packs	03896215
Sterile, p	backed in pairs	5		

S	80915	50 pairs	12 packs	11362864
М	8091M	50 pairs	12 packs	11362870
L	8091L	50 pairs	12 packs	11362887

# Gentle Skin<sup>®</sup> steril

## Our sterile latex examination glove

#### Features

- sterile latex examination glove (DIN EN 455 standard)
- powder-free, low protein
- micro-roughened
- elastic and tear-proof
- ambidextrous
- packed singly or in pairs





#### Usage

Suitable for use in any surgical application. Anywhere where wearers are working in a sterile environment. For smaller surgical interventions, changing dressings, taps, etc.

Sterile, individually packed				
Size	REF	Pack contents	Shipping unit	PZN
S	8022S	50 pcs	8 packs	06402811
М	8022M	50 pcs	8 packs	06402828
L	8022L	50 pcs	8 packs	06402834
Storilo				

Julie,	backed in pairs			
S	8021RSW	50 pairs	4 packs	02244479
М	8021RMW	50 pairs	4 packs	02244485
L	8021RLW	50 pairs	4 packs	02244491

# Nitril<sup>®</sup> sterile

## Our sterile nitrile examination glove

#### Features

- sterile nitrile examination glove (DIN EN 455 standard)
- PPE Cat III
- micro-roughened
- elastic and tear-proof
- ambidextrous



#### Usage

Usable for surgical applications. Anywhere where wearers are working in a sterile environment. Usable for smaller surgical interventions, changing dressings, taps, etc. Suitable for latex allergy sufferers with a diagnosed type I allergy.

Sterile, p	acked in pairs	5		
Size	REF	Pack contents	Shipping unit	PZN
S	80235	50 pairs	4 packs	09544658
М	8023M	50 pairs	4 packs	09544641
L	8023L	50 pairs	4 packs	09544635
XL	8023XL	50 pairs	4 packs	11362976

# Our surgical gloves set the standards for protection and safety

During surgery, protecting medical staff and patients is of prime importance. The Meditrade surgical glove range includes a wide selection of glove models made from a variety of materials. Common to all gloves is their outstanding mechanical durability, thickness, grip, comfort and biological compatibility.

Surgical gloves



# **Reference<sup>™</sup> OP**

### The powdered surgical glove made from latex

### Features

- powdered latex surgical glove, natural (DIN EN 455 standard)
- PPE Cat III
- micro-roughened on the fingers
- lightly powdered
- anatomical fit modelled on natural hand posture
- left- and right-handed, with rolled edge
- packed in pairs
- sterile (ISO 11137-2 standard)

### Usage

Suitable for use in all surgical areas that require sterility and anywhere where the user is working in a sterile environment with an identified risk of contamination. For highly infectious situations in hospitals, general practice, in laboratories and in the pharmaceutical industry. Also suitable for double-gloving.



Size	REF	Pack contents	Shipping unit	PZN
6	93216	50 pairs	6 packs	02243824
6.5	932165	50 pairs	6 packs	02243830
7	93217	50 pairs	6 packs	02243847
7.5	932175	50 pairs	6 packs	02243853
8	93218	50 pairs	6 packs	02243876
8.5	932185	50 pairs	6 packs	02243882
9	93219	50 pairs	6 packs	05030313



# **Gentle Skin**<sup>®</sup> Premium OP<sup>™</sup>

### The powder-free surgical glove made from latex

### Features

- powder-free latex surgical glove, natural (DIN EN 455 standard)
- PPE Cat III
- with polymer coating
- micro-roughened on the fingers
- anatomical fit modelled on natural hand posture
- left- and right-handed, with rolled edge
- packed in pairs
- sterile (ISO 11137-2 standard)

### Usage

Especially suited for vascular surgery, orthopaedics, abdominal surgery and in all surgical fields with sterile requirements. Anywhere there is an identified risk of contamination for the user in a sterile environment. For highly infectious situations in hospitals, general practice, in laboratories and in the pharmaceutical industry. Also suitable for double-gloving.



Size	REF	Pack contents	Shipping unit	PZN
5.5	902155W	50 pairs	4 packs	05030299
6	90216W	50 pairs	4 packs	02244114
6.5	902165W	50 pairs	4 packs	02244278
7	90217W	50 pairs	4 packs	02244290
7.5	902175W	50 pairs	4 packs	02244315
8	90218W	50 pairs	4 packs	02244338
8.5	902185W	50 pairs	4 packs	02244427
9	90219W	50 pairs	4 packs	02244433

Surgical gloves



# Gentle Skin<sup>®</sup> Superior OP

### The latex surgical glove

### Features

- powder-free latex surgical glove, natural (DIN EN 455 standard)
- PPE Cat III
- anatomical fit modelled on natural hand posture
- left- and right-handed, with rolled edge
- packed in pairs
- sterile (ISO 11137-2 standard)

### Usage

Suitable for use in all surgical areas that require sterility and anywhere where the user is working in a terile environment with an identified risk of contamination. For highly infectious situations in hospitals, general practice, in laboratories and in the pharmaceutical industry. Also suitable for double-gloving.



Size	REF	Pack contents	Shipping unit	PZN
5.5	904155	50 pairs	4 packs	-
6	904160	50 pairs	4 packs	-
6.5	904165	50 pairs	4 packs	-
7	904170	50 pairs	4 packs	-
7.5	904175	50 pairs	4 packs	-
8	904180	50 pairs	4 packs	-
8.5	904185	50 pairs	4 packs	-
9	904190	50 pairs	4 packs	-



# Gentle Skin<sup>®</sup> Smooth OP

### The latex surgical glove with a smooth texture

### Features

- powder-free latex surgical glove, natural (DIN EN 455 standard)
- PPE Cat III
- with polymer coating
- anatomical fit modelled on natural hand posture
- left- and right-handed, with rolled edge
- packed in pairs
- sterile (ISO 11137-2)

### Usage

Especially suited for vascular surgery, orthopaedics, abdominal surgery and in all surgical fields with sterile requirements. Anywhere there is an identified risk of contamination for the user in a sterile environment. For highly infectious situations in hospitals, general practice, in laboratories and in the pharmaceutical industry. Also suitable for double-gloving.



Size	REF	Pack contents	Shipping unit	PZN
6	90266	50 pairs	4 packs	11154775
6.5	902665	50 pairs	4 packs	11154781
7	90267	50 pairs	4 packs	11154798
7.5	902675	50 pairs	4 packs	11154806
8	90268	50 pairs	4 packs	11154812
8.5	902685	50 pairs	4 packs	11154829
9	90269	50 pairs	4 packs	11154835



# Gentle Skin<sup>®</sup> Micro OP<sup>®</sup>

### The touch-sensitive surgical glove made from latex

### Features

- powder-free latex surgical glove, brown (DIN EN 455 standard)
- PPE Cat III
- micro-roughened on the fingers
- anatomical fit modelled on natural hand posture
- left- and right-handed, with rolled edge
- packed in pairs
- sterile (ISO 11137-2)

### Usage

The surgical glove for handling demanding applications where particularly high touch sensitivity is needed. Generally suitable for use in all surgical areas that require sterility and anywhere where the user is working in a sterile environment with an identified risk of contamination. For highly infectious situations in hospitals, general practice, in laboratories and in the pharmaceutical industry. Also suitable for double-gloving.



Size	REF	Pack contents	Shipping unit	PZN
6	90316	50 pairs	4 packs	06972460
6.5	903165	50 pairs	4 packs	06972589
7	90317	50 pairs	4 packs	06972595
7.5	903175	50 pairs	4 packs	06972626
8	90318	50 pairs	4 packs	06972655
8.5	903185	50 pairs	4 packs	06972661



### Maxitex<sup>®</sup> UG

• powder-free latex underglove, green

• can be combined with any primary

(DIN EN 455 standard)

• anatomical fit modelled on

natural hand posture

• left- and right-handed,

with rolled edge

• packed in pairs

surgical glove

### The surgical underglove made from latex

### Features

### Usage

As a surgical underglove for use in double-gloving for greater safety during operations. Double-gloving increases the level of protection from contamination and infection by reducing the risk of perforation. Suitable for invasive interventions, especially for high-risk patients, as well as for emergency service use and interventions with a much higher risk of

pathogen transfer.

9

1202007



Sterile, pac	ked in pairs			
Size	REF	Pack contents	Shipping unit	PZN
6	1202001	50 pairs	4 packs	-
6.5	1202002	50 pairs	4 packs	-
7	1202003	50 pairs	4 packs	-
7.5	1202004	50 pairs	4 packs	-
8	1202005	50 pairs	4 packs	-
8.5	1202006	50 pairs	4 packs	_

50 pairs

4 packs

• sterile (ISO 11137-2 standard)

#### 40



# Gentle Skin<sup>®</sup> Securitex<sup>®</sup>OP

### The surgical glove set made from latex with indicator function

### Features

### Usage

- powder-free latex surgical glove, green and brown (DIN EN 455 standard)
- micro-roughened on the fingers
- anatomical fit modelled on natural hand posture
- left- and right-handed, with rolled edge
- 2x packed in pairs
- sterile (ISO 11137-2)

Suitable for surgery and interventions with a high risk of infection that require a high level of protection, in hospitals and general practice, in labs and the pharmaceutical industry. For surgical applications that require sterility and anywhere where the user is working in a sterile environment with an identified risk of contamination.



### Sterile set, 2x packed in pairs

Size	REF	Pack contents	Shipping unit	PZN
6	90516	25 sets	4 packs	11349970
6.5	905165	25 sets	4 packs	03851224
7	90517	25 sets	4 packs	03851230
7.5	905175	25 sets	4 packs	03851247
8	90518	25 sets	4 packs	03851253
8.5	905185	25 sets	4 packs	03851276



# Gentle Skin<sup>®</sup> Isopretex<sup>®</sup>OP

### The latex-free surgical glove made from polyisoprene

### Features


- powder-free polyisoprene surgical glove, beige (DIN EN 455 standard)
- PPE Cat III
- micro-roughened on the fingers
- anatomical fit modelled on natural hand posture
- Left- and right-handed, with rolled edge
- packed in pairs
- sterile (ISO 11137-2 standard)

### Usage

Latex-free glove for all surgical areas that require sterility, and anywhere where the user is working in a sterile environment with an identified risk of contamination. For highly infectious situations in hospitals, general practice, in laboratories and in the pharmaceutical industry. Also suitable for double-gloving.



Size	REF	Pack contents	Shipping unit	PZN
5.5	952255	50 pairs	4 packs	11190854
6	95226	50 pairs	4 packs	11190860
6.5	952265	50 pairs	4 packs	11190877
7	95227	50 pairs	4 packs	11190883
7.5	952275	50 pairs	4 packs	11190908
8	95228	50 pairs	4 packs	11190914
8.5	952285	50 pairs	4 packs	11190920
9	95229	50 pairs	4 packs	11190937



# Neopretex<sup>®</sup> OP

### The latex-free surgical glove made from polychloroprene

### Features

- powder-free polychloroprene surgical glove, brown (DIN EN 455 standard)
- PPE Cat III
- micro-roughened on the fingers
- anatomical fit modelled on natural hand posture
- left- and right-handed, with rolled edge
- packed in pairs
- sterile (ISO 11137-2)

### Usage

Latex-free glove for all surgical areas that require sterility, and anywhere where the user is working in a sterile environment with an identified risk of contamination. For highly infectious situations in hospitals, general practice, in laboratories and in the pharmaceutical industry. Also suitable for double-gloving.



Size	REF	Pack contents	Shipping unit	PZN
6	95216	50 pairs	4 packs	02426890
6.5	952165	50 pairs	4 packs	02426909
7	95217	50 pairs	4 packs	02426915
7.5	952175	50 pairs	4 packs	02426921
8	95218	50 pairs	4 packs	02426938
8.5	952185	50 pairs	4 packs	02426944
9	95219	50 pairs	4 packs	02426950


Finessis. An innovative surgical glove programme, that reinvents

### **PURITY, COMFORT** and **REACTIVE SAFETY**

### Made from FLEXYLON<sup>™</sup> – a high-performance elastomer

Flexylon<sup>™</sup> is a compact poured film that contains no accelerators or cross-linkers.

This approach eliminates virtually all stitch holes and other micro-defects, thereby achieving the very highest quality standards (AQL = 0.1)!



### **Reactive** protection

Outer layer made from FLEXYLON™

#### Number of HIV1 viruses transferred by a 22G hollow bore needle Reduction factor: 98 98 871 900 800 700 600 500 400 300 200 100 17 Double-glove system No glove Finessis

### Inner layer made from FLEXYLON™



Figure 1: A sharp object presses against the outer layer, but without breaking through.



Figure 2:As pressure is increased, the thin membranes tear that separate the gel drops form one another. A larger, contiguous quantity of the antimicrobial, gel-like liquid is created.



Figure 3: As soon as the sharp object breaks through, it is enveloped by the solution and the viruses or bacteria are effectively neutralised.

# **Finessis**

## Finessis. Aegis®

### **Reactive protection**

### Trauma and high-risk Orthopaedics

### **FLEXYLON<sup>™</sup>** | Without accelerators | Synthetic surgical glove

Aegis<sup>®</sup> is a revolutionary protective system. It is a surgical glove that is equipped with a reaction-triggered disinfection system. Finessis<sup>®</sup> Aegis<sup>®</sup> was designed as a glove interior featuring an antimicrobial layer between the sheath, which consists of high-performance elastomers. The disinfection system is activated and effective only IF the glove is PERFORATED. In this case, the activated gel is released in the RIGHT QUANTITY and at the EXACT LOCATION where the perforation has occurred.

### **Special characteristics:**

- powder-free surgical glove (DIN EN 455 standard)
- PPE glove (DIN EN 374 standard)
- reaction-initiated disinfection system
- WITHOUT accelerators
- antimicrobial intermediate layer, which releases activated gel in the event of a perforation
- best AQL in the industry (AQL = 0.10)
- manufactured from FLEXYLON<sup>™</sup> (layer technology)
- high electrical resistance
- anatomical fit

Size	REF	Pack contents	Shipping unit
6	90286	40 pairs	4 packs
6.5	902865	40 pairs	4 packs
7	90287	40 pairs	4 packs
7.5	902875	40 pairs	4 packs
8	90288	40 pairs	4 packs
8.5	902885	40 pairs	4 packs
9	90289	40 pairs	4 packs



# Finessis, Zero<sup>®</sup>

### Performance and safety General surgery

### **FLEXYLON™** | Without accelerators | Synthetic surgical glove

Finessis<sup>®</sup> Zero<sup>®</sup> is the next generation of surgical gloves that are manufactured without any accelerators to achieve a new level of PURITY. By using FLEXYLON<sup>™</sup> layer technology, the risk of a type I skin allergy is completely eliminated and the risk of a type IV skin allergy is reduced to a previously unattainable minimum.

### **Special characteristics:**

- powder-free surgical glove (DIN EN 455)
- PPE glove (DIN EN 374)
- NO accelerators
- best AQL in the industry (AQL = 0.10)
- manufactured from FLEXYLON<sup>™</sup> (layer technology)
- high electrical resistance
- anatomical fit

Size	REF	Pack contents	Shipping unit
6	90276	50 pairs	4 packs
6.5	902765	50 pairs	4 packs
7	90277	50 pairs	4 packs
7.5	902775	50 pairs	4 packs
8	90278	50 pairs	4 packs
8.5	902785	50 pairs	4 packs
9	90279	50 pairs	4 packs
		1	L were







# Finessis Corium®

### Microsurgery and underglove

### **FLEXYLON<sup>™</sup>** | Without accelerators | Synthetic surgical glove

Finessis<sup>®</sup> Corium<sup>®</sup> is the glove with a high tactility for double-gloving and microsurgery. The glove was developed to ensure improved touch sensitivity. Corium<sup>®</sup> works as an underglove, to achieve optimal comfort without losing grip or traction. When combined with the Finessis<sup>®</sup> Zero<sup>®</sup> for the perfect fit, the Finessis<sup>®</sup> double-gloving system presents an impermeable barrier for any surgical application.

### **Special characteristics:**

- powder-free surgical glove (DIN EN 455)
- PPE glove (DIN EN 374)
- NO accelerators
- best AQL in the industry (AQL = 0.10)
- manufactured from FLEXYLON<sup>™</sup> (layer technology)
- smooth surface
- high wearer comfort with improved tactility
- high electrical resistance
- anatomical fit

Size	REF	Pack contents	Shipping unit
6	90246	50 pairs	4 packs
6.5	902465	50 pairs	4 packs
7	90247	50 pairs	4 packs
7.5	902475	50 pairs	4 packs
8	90248	50 pairs	4 packs
8.5	902485	50 pairs	4 packs
9	90249	50 pairs	4 packs





**Finessis** 


# Improving safety in operating theatres and wards

Our products make an important contribution to stopping infection pathways and the transmission of germs. Personnel, patients and visitors should feel safe and comfortable as they interact with each other. Wearer comfort, breathability and an adjustable fit are important aspects for suitable clothing.

For protective clothing in particular, it is not enough to simply rely on general or universal products: instead, clothing needs to meet the high standards in the medical field. When making Suavel face masks, for example, our techniques include an ultrasonic sealing method that prevents stitches coming loose in the final garment.

### Improving safety in operating theatres and wards

### Features of **Suavel**<sup>®</sup> surgical masks



# Suavel<sup>®</sup> Protec Plus

# Suavel<sup>®</sup> Protec Plus

### Our high quality surgical mask with elastic ear loops, made in Germany

### Features

- Type II/II R surgical mask (EN 14683 standard)
- bacterial filter efficiency  $\geq$  98%
- high-performance filter with low respiratory resistance
- 3-ply, green exterior
- elastic ear loops
- genuine surgical non-woven
- cool and hypoallergenic
- disposable
- non-sterile



### Our surgical mask with superior filtering and elastic ear loops

### Features

- Type II surgical mask (EN 14683 standard)
- bacterial filter efficiency  $\geq$  98%
- 3-ply, green exterior
- elastic ear loops
- genuine surgical non-woven
- hypoallergenic
- disposable
- non-sterile



### Usage

For all surgical applications and for ward use, for the protection of both staff and patients.

### Non-sterile

REF	Туре	Colour	Pack contents	Shipping unit	PZN
80-956	Ш	Green	50 pcs	20 packs	16760256
80-957	II R	White	50 pcs	20 packs	_



### Usage

For all surgical applications and for ward use, for the protection of both staff and patients.

Non-sterile						
REF	Colour	Pack contents	Shipping unit	PZN		
80-955MP	Green	50 pcs	18 packs	02244605		

# Suavel<sup>®</sup> Protec

# Suavel<sup>®</sup> Protec Black

### Our quality surgical mask with elastic ear loops

### Features

- Type II/II R surgical mask (EN 14683 standard)
- bacterial filter efficiency  $\geq$  98%
- 3-ply, green or blue exterior
- elastic ear loops
- disposable
- non-sterile



### Features

- Type II surgical mask (EN 14683 standard))
- bacterial filter efficiency  $\geq 98\%$
- 3-ply, black exterior
- elastic ear loops
- disposable
- non-sterile



### Usage

For all surgical applications and for ward use, for the protection of both staff and patients.

### Non-sterile

REF	Туре	Colour	Pack contents	Shipping unit	PZN
80-900	II	Green	50 pcs	18 packs	00068771
80-901	II	Blue	50 pcs	20 packs	16731071
80-902	II R	Blue	50 pcs	40 packs	16731088



### Usage

For all surgical applications and forward use, for the protection of both staff and patients.

Non-s	sterile			
REF	Colour	Pack contents	Shipping unit	PZN
80-904	Black	50 pcs	20 packs	-

# Suavel<sup>®</sup> Comfort Plus Suavel<sup>®</sup> Comfort

### Our tie-on surgical mask with superior filtering

### Features

- Type II surgical mask (EN 14683 standard)
- bacterial filter efficiency  $\geq$  98%
- 3-ply, blue exterior
- tie-on mask
- long ties
- genuine surgical non-woven
- hypoallergenic
- disposable
- non-sterile

### Our quality tie-on surgical mask

### Features

- Type II/II R surgical mask (EN 14683 standard)
- bacterial filter efficiency  $\geq$  98%
- 3-ply, blue exterior
- tie-on mask
- long ties
- disposable
- non-sterile



### Usage

For all surgical applications and for ward use, for the protection of both staff and patients.

# Non-sterileREFColourPack contentsShipping unitPZN80-410MBlue50 pcs12 packs02244611



### Usage

For all surgical applications and for ward use, for the protection of both staff and patients.

Non-sterile						
REF	Туре	Colour	Pack contents	Shipping unit	PZN	
80-400	Ш	Blue	50 pcs	12 packs	00068788	
80-402	II R	Blue	50 pcs	20 packs	_	

# Suavel<sup>®</sup> Antifluid

### **Our surgical mask** for spectacle wearers

### **Features**

- with effective anti-fogging insulation
- Type II surgical mask (EN 14683 standard)
- bacterial filter efficiency  $\geq$  98%
- 3-ply, blue exterior
- tie-on mask
- long ties
- disposable
- non-sterile



### **Our surgical mask** with face shield

### **Features**

- High degree of protection against infection with integrated face shield for mouth, eves and nose
- liquid stop layer
- Type II R surgical mask • (EN 14683 standard)
- Bacterial filter efficiency ≥98%
- 4-ply, blue exterior
- tie-on mask
- long ties
- disposable
- non-sterile





High degree of protection against infection with integrated face shield for mouth, eyes and nose, plus liquid stop layer

### Usage

Specially designed for spectacle wearers. For all surgical applications and for ward use, for the protection of both staff and patients.

Non-s	terile			
REF	Colour	Pack contents	Shipping unit	PZN
80-470M	Blue	50 pcs	12 packs	00559457

### Usage

For highly infectious situations, protects user against splashes of blood/wound secretions, also for rinsing operations. For all surgical applications and for ward use.

Non-s	terile			
REF	Colour	Pack contents	Shipping unit	PZN
80-455M	Blue	50 pcs	12 packs	00559517

# **Suavel**<sup>®</sup> Sensima

### Our surgical mask for allergy sufferers

### Features

- hypoallergenic
- Type II surgical mask (EN 14683 standard) for users with sensitive skin
- Bacterial filter efficiency  $\geq$  98%
- 3-ply, white exterior
- tie-on mask
- long ties
- disposable





#### Usage

For allergy sufferers and users with sensitive skin. For all surgical applications and for ward use, for the protection of both staff and patients.

Non-s	terile			
REF	Colour	Pack contents	Shipping unit	PZN
80-440M	White	50 pcs	12 packs	11362700

### **Suavel**<sup>®</sup> Apollo Astronaut-shaped surgical cap

### Features

- surgical cap in astronaut helmet style, made from a tear-resistant polypropylene non-woven
- with long ties and neck protection
- breathable
- with integrated sweatband





### Usage

For operations, for covering the hair and preventing excessive perspiration. Suitable for beard wearers.

Non-sterile						
REF	Colour	Pack contents	Shipping unit	PZN		
70-131G	Green	100 pcs	10 packs	11362597		

# Suavel<sup>®</sup> Astrid

Our universal bouffant cap

### Features

- beret cap made from soft, lightweight non-woven polypropylene
- with double-sewn rubber drawstring
- breathable





### Usage

For use in hospitals and in the healthcare sector, for both nursing staff and patients. Usable anywhere were the hair has to be covered, also including applications in the food and catering sector.

### Suavel<sup>®</sup> Astrid, non-sterile

REF	Size	Colour	Pack contents	Shipping unit	PZN
70-121B	51 cm	Blue	100 pcs	10 packs	11362611
70-121G	51 cm	Green	100 pcs	10 packs	11362628
70-121W	51 cm	White	100 pcs	10 packs	11362605

### Suavel<sup>®</sup> Astrid XL, non-sterile

70-124G	60 cm	Green	100 pcs	10 packs	11362634
70-124W	60 cm	White	100 pcs	10 packs	11362640

### BeeSana<sup>®</sup> PP gown 35 g

### The high-performance

scrubs for personnel and visitors

### Features

- scrubs made from 35 g/m<sup>2</sup> PP non-woven
- fastening ties at the neck and waist
- cuffs at the end of the sleeves
- liquid-repellent
- non-sterile
- disposable

### Usage

Ward and visitor gowns to prevent the spread of germs in hospitals and in geriatric nursing; for personnel, patients and visitors.



Non-sterile						
REF	Size	Colour	Pack contents	Shipping unit	PZN	
3700	140 x 140 cm	Green	10 pcs	10 packs	09747462	

### BeeSana<sup>®</sup> PP gown 23 g

### The industry-standard

### protective gown for personnel and visitors

### Features

- protective gown made from 23 g/m<sup>2</sup> PP non-woven
- fastening ties at the neck and waist
- sleeves come with either sewn-in, elastic rubber ties or with cuffs
- lightweight material
- moisture-repellent
- non-sterile
- sisposable

### Usage

Ward and visitor gowns to prevent the spread of germs in hospitals and in geriatric nursing; for personnel, patients and visitors.





### Non-sterile, with elastic rubber ties

REF	Size	Colour	Pack contents	Shipping unit	PZN
3684	110 x 140 cm	White	10 pcs	10 packs	11509391
3685	110 x 140 cm	Green	10 pcs	10 packs	11509416
3699	120 x 140 cm	Yellow	10 pcs	10 packs	09747479
3729	130 x 150 cm	White	10 pcs	10 packs	11527609

### Non-sterile, with elastic cuffs

	-				
4685 1	110 x 140 cm	Green	10 pcs	10 packs	11363036
4699 1	120 x 140 cm	Yellow	10 pcs	10 packs	11363042

### BeeSana<sup>®</sup> PP/PE gown 26 g

### The waterproof

### semi-laminated protective gown for personnel and visitors

### Features

- protective gown made from 26 g/m<sup>2</sup> PP non-woven with PE coating covering front/arm areas
- fastening ties at the neck and waist
- cuffs at the end of the sleeves
- waterproof across the coated areas
- non-sterile
- disposable

### Usage

Ward and visitor gowns for applications requiring a high level of protection to prevent the spread of germs in hospitals and in geriatric nursing; for personnel, patients and visitors.



### Non-sterile

REF	Size	Colour	Pack contents	Shipping unit	PZN
4698	136 x 136 cm	Yellow	10 pcs	10 packs	11642842
4698-XL	150 x 170 cm	Yellow	10 pcs	10 packs	17376146

### BeeSana<sup>®</sup> SMS gown 40 g

### **High protection**

Protective gown offering a high level of protection against liquids and particulates

### Features

- protective gown made from soft, 3-ply SMS composite material
- combined layers provide a highly effective barrier against fluids, germs and light splashes of chemicals
- fastening ties at the neck and waist and cuffs at the end of the sleeves
- non-sterile
- disposable

### Usage

Suitable for highly infectious situations. To prevent the spread of germs on emergency and isolation wards, in care homes, and in hospitals, for personnel, patients and visitors.



Nor	n-sterile				
REF	Size	Colour	Pack contents	Shipping unit	PZN
3597	115 x 150 cm	Blue	10 pcs	5 packs	11652154



### BeeSana<sup>®</sup> SMS gown 15 g

### Suavel<sup>®</sup>Yeti Coverall

### The lightweight

protective gown for personnel and visitors

### Features

- protective gown made from soft, 3-ply 15 g/m<sup>2</sup> SMS composite material
- fastening ties at the neck and waist and cuffs at the end of the sleeves
- non-sterile
- disposable

### Usage

For highly infectious situations. To prevent the spread of germs on emergency and isolation wards, in care homes, and in hospitals, for personnel, patients and visitors.





	Julia				
REF	Size	Colour	Pack contents	Shipping unit	PZN
3598	140 x 140 cm	Blue	10 pcs	10 packs	11518680

### Protective suit

#### Features

- protective suit made from polypropylene with an air-permeable coating, offering protection against harmful dusts (type 5) and splashes of liquid (type 6)
- elastic rubber ties at wrists and ankles
- concealed zip closure at front
- non-sterile
- disposable

### Usage

Category III protective suit for use in hospitals, e.g. in the ongoing pandemic situation. Offers protection against harmful dusts and splashes of liquid.

### Non-sterile

REF	Size	Height	Chest	Packing unit
4798	М	168 - 176 cm	92 - 100 cm	1 pc./box 50 packs
4799	L	174 - 182 cm	100 - 108 cm	1 pc./box 50 packs
4800	XL	182 - 188	108 - 116 cm	1 pc./box 50 packs
4801	XXL	188 - 194	116 - 124 cm	1 pc./box 50 packs

Non-sterile

# **PE disposable apron**

Our PE apron in various models

#### Features

- disposable LDPE apron
- tested to the EN 14126 standard
- with closing ties at the waist to ensure a good fit
- lightweight material
- hard-wearing and moisture-repellent
- in a handy pouch or on hanger block



#### Usage

Protection for clothing in nursing/care activities, and in all areas with stringent hygiene requirements, including large caterers and hospitals.

### Non-sterile, in pouch, for individual removal

REF	Size	Thickness	Colour	Pack contents	Shipping unit	PZN
3591	140 x 75 cm	30 µm	Blue	50 pcs	10 packs	11857643
3690	125 x 81 cm	16 µm	White	100 pcs	10 packs	04353379
3691	140 x 81 cm	16 µm	White	100 pcs	10 packs	04353422
3692	160 x 81 cm	16 µm	White	100 pcs	10 packs	04353474
3696	117 x 69 cm	16 µm	White	100 pcs	10 packs	-

### Non-sterile, on hanger block, for tearing-off individually

	3687	125 x 81 cm	16 µm	White	50 pcs	10 packs	04353385
	3688	140 x 81 cm	16 µm	White	50 pcs	10 packs	04353468
	3694	160 x 81 cm	16 µm	White	50 pcs	10 packs	04353480
-	3697	117 x 69 cm	16 µm	White	200 pcs	5 rolls	-

### **Suave** OP PRO Standard SMMS surgical scrubs, sterile

#### Features

- surgical scrubs made from 35 g/m<sup>2</sup> SMMS composite material (DIN EN 13795 standard)
- combined layers provide a highly effective barrier against fluids, germs and light chemical spills
- fastening ties with sterile card at the waist and cuffs at the end of the sleeves; hook-and-pile fastener at neck
- SMS wrapping cloth plus two paper towels
- counting cards on the peel pack
- sterile and disposable

#### Usage

During surgical interventions, wearing the Suavel<sup>®</sup> OP PRO Standard prevents the transfer of particulates and infectious material from the surgical team to the surgical wound, and vice versa.

### Sterile

REF	Size	Colour	Pack contents	Shipping unit	PZN
5541	M: 120 x 140 cm	Blue	1 pc.	30 packs	15295662
5542	L: 130 x 150 cm	Blue	1 pc.	30 packs	15295679
5543	XL: 150 x 160 cm	Blue	1 рс.	30 packs	15295685
5544	XXL: 170 x 160 cm	Blue	1 pc.	30 packs	15295691



### Suave OP PRO Expert SMMS surgical scrubs, sterile, reinforced

### Features

- surgical scrubs made from 43 g/m<sup>2</sup> SMMS composite material (DIN EN 13795 standard)
- combined layers provide a highly effective barrier against fluids, germs and light chemical spills
- fastening ties with sterile card at the waist and cuffs at the end of the sleeves; hook-and-pile fastener at neck
- reinforced chest, abdominal area and sleeves
- SMS wrapping cloth plus two paper towels
- counting cards on the peel pack
- sterile and disposable

### Usage

During surgical interventions, wearing the Suavel<sup>®</sup> OP PRO Expert prevents the transfer of particulates and infectious material from the surgical team to the surgical wound, and vice versa.

### Sterile

REF	Size	Colour	Pack contents	Shipping unit	PZN
5551	M: 120 x 140 cm	Blue	1 pc.	30 packs	15295716
5552	L: 130 x 150 cm	Blue	1 рс.	30 packs	15295722
5553	XL: 150 x 160 cm	Blue	1 pc.	30 packs	15295739
5554	XXL: 170 x 160 cm	Blue	1 pc.	30 packs	15295745





# Wound treatment and dressings – the pillars that shoulder the weight of medical care

Alongside our wound dressings, we aim to offer the ideal range of products to treat every type of wound at each stage of recovery. After all, a wound dressing alone cannot hope to satisfy every aspect or goal for treatment.

Our product range includes all non-sterile and sterile gauze/non- woven swabs as well as absorbent dressing pads. In the field of therapeutic dressings, we aim to offer a comprehensive set of products for an array of applications that include dressing, supporting, and relieving the wound.

Patient safety is at the forefront of everything we do, and all our sterile products are fully compliant with the relevant safety standards. Features here include woven-in x-ray contrast threads and the counting check card.

# BeeSana<sup>®</sup> Gauze cutting

### For individual wound coverings

### Features

- pre-cut gauze (EN 14079 standard)
- 20-thread
- absorbent



### Usage

For individual wound covering applications. For the external wound care of exuding and/or contaminated wounds. For blood and exudate absorption.

Non-s	terile			
REF	Size (width x length)	Pack contents	Shipping unit	PZN
1330	20 x 20 cm (approx.)	500 pcs	10 packs	-
1335	30 x 40 cm (approx.)	500 pcs	6 packs	_

# BeeSana® Medical gauze

### For individual wound coverings

### Features

- medical gauze (EN 14079 standard)
- 20-thread

### Usage

For preparing individual bandages and wound coverings.

sterile			
e, 1-ply and 4-ply			
Size (width x length)	Pack contents	Shipping unit	PZN
80 cm x 40 m (approx.)	5 pcs	2 packs	04852674
e sterile with self-a	dhesive counting c	ards	
	e, 1-ply and 4-ply Size (width x length) 80 cm x 40 m (approx.)	e, 1-ply and 4-ply Size (width x length) Pack contents 80 cm x 40 m (approx.) 5 pcs	e, 1-ply and 4-ply Size Pack contents Shipping unit (width x length)

80 cm wide, 5-ply, folded to 16 x 20 cm, 3 sterile pcs per pack

REF	Length	Pack contents	Shipping unit	PZN
4179	100 cm	8 x 3 pcs	2 packs	13581004

# **BeeSana**<sup>®</sup> Gauze swab

Our standard gauze swab with selvedge cut edges

### Features

- gauze swab, medical gauze (EN 14079 standard)
- selvedge cut edges, 100% cotton
- folded, 17-thread
- absorbent



### Usage

For general-purpose wound care of exuding wounds and for absorbing bodily fluids. For applying ointments. For smaller-scale outpatient surgical interventions, in wards and for first aid applications.

### Non-sterile

		REF	Size (width x length)	Pack contents	Shipping unit	PZN
8-ply		1000	5 x 5 cm (approx.)	100 pcs	10 packs	04852705
		1006	7.5 x 7.5 cm (approx.)	100 pcs	10 packs	04852711
		1012	10 x 10 cm (approx.)	100 pcs	10 packs	04852728
		1053	10 x 20 cm (approx.)	100 pcs	5 packs	01330633
	10 pcs, pre-counted For set packers	2712	10 x 10 cm (approx.)	200 x 10 pcs	1 box	-
12-ply		1002	5 x 5 cm (approx.)	100 pcs	10 packs	04852734
		1008	7.5 x 7.5 cm (approx.)	100 pcs	10 packs	04852740
		1014	10 x 10 cm (approx.)	100 pcs	10 packs	04852757
		1038	10 x 13 cm (approx.)	100 pcs	10 packs	04852763
		1023	10 x 20 cm (approx.)	100 pcs	10 packs	04852786
	10 pcs, pre-counted	2708	7.5 x 7.5 cm (approx.)	200 x 10 pcs	1 box	-
	For set packers	2714	10 x 10 cm (approx.)	200 x 10 pcs	1 box	-
		2723	10 x 20 cm (approx.)	100 x 10 pcs	1 box	-
16-ply		1004	5 x 5 cm (approx.)	100 pcs	10 packs	04852792
		1010	7.5 x 7.5 cm (approx.)	100 pcs	10 packs	04852800
		1016	10 x 10 cm (approx.)	100 pcs	10 packs	04852817
		1052	10 x 13 cm (approx.)	100 pcs	10 packs	04852823
		1025	10 x 20 cm (approx.)	100 pcs	5 packs	11849589
32-ply		1018	10 x 10 cm (approx.)	50 pcs	10 packs	_
64-ply	20-thread	1235	10 x 13 cm (approx.)	25 pcs	8 packs	11849626

# **BeeSana**<sup>®</sup> Gauze swab

Our standard gauze swab with selvedge cut edges

### Features

- gauze swab, medical gauze (EN 14079 standard)
- selvedge cut edges, 100% cotton
- folded, 17-thread
- absorbent





### Single sterile

		REF	Size (width x length)	Pack contents	Shipping unit	PZN
8-ply	Pack of 2 pcs	1043	5 x 5 cm (approx.)	25 x 2 pcs	27 packs	04852846
	-	1044	7.5 x 7.5 cm (approx.)	25 x 2 pcs	16 packs	04852852
	-	1045	10 x 10 cm (approx.)	25 x 2 pcs	16 packs	04852869
	-	1046	10 x 20 cm (approx.)	25 x 2 pcs	8 packs	04852875
	Pack of 5 pcs	4600	5 x 5 cm (approx.)	80 x 5 pcs	8 packs	13458367
	_	4601	7.5 x 7.5 cm (approx.)	80 x 5 pcs	8 packs	01330662
	-	1062	10 x 10 cm (approx.)	35 x 5 pcs	6 packs	01330679
	Pack of 10 pcs	4606	10 x 10 cm (approx.)	50 x 10 pcs	4 packs	12348949
12-ply	Pack of 2 pcs	4130	5 x 5 cm (approx.)	25 x 2 pcs	24 packs	07216329
	-	4131	7.5 x 7.5 cm (approx.)	25 x 2 pcs	18 packs	07216335
	-	4132	10 x 10 cm (approx.)	25 x 2 pcs	12 packs	07216341
	-	4133	10 x 20 cm (approx.)	25 x 2 pcs	6 packs	07216358
	Pack of 5 pcs	1064	5 x 5 cm (approx.)	30 x 5 pcs	12 packs	01330722
	_	1065	7.5 x 7.5 cm (approx.)	30 x 5 pcs	6 packs	01330739
	_	1063	10 x 10 cm (approx.)	30 x 5 pcs	6 packs	01330745
	_	1066	10 x 20 cm (approx.)	15 x 5 pcs	6 packs	01330751
	Pack of 10 pcs	4616	5 x 5 cm (approx.)	50 x 10 pcs	8 packs	12348955
	-	4617	7.5 x 7.5 cm (approx.)	30 x 10 pcs	8 packs	12349044
	-	1026	10 x 10 cm (approx.)	40 x 10 pcs	4 packs	12348694
	-	4619	10 x 20 cm (approx.)	20 x 10 pcs	4 packs	12349015
	Pack of 20 pcs	4622	10 x 10 cm (approx.)	20 x 20 pcs	4 packs	12349021
	Pack of 30 pcs	4640	7.5 x 7.5 cm (approx.)	25 x 30 pcs	4 packs	12349185
16-ply	Pack of 4 pcs	1092	7.5 x 7.5 cm (approx.)	64 x 4 pcs	8 packs	12348783
32-ply	Pack of 2 pcs	4628	10 x 10 cm (approx.)	30 x 2 pcs	8 packs	12349038



### Single sterile, rolled gauze dressing pad

		REF	Size (width x length)	Pack contents	Shipping unit	PZN
16-ply	Pack of 2 pcs	1146	20 x 20 cm (approx.)	20 x 2 pcs	4 packs	13582587



### Double sterile, with self-adhesive counting cards

16-ply	Pack of 10 pcs	1058	10 x 10 cm (approx.)	25 x 10 pcs	4 packs	13582601
	Pack of 60 pcs	4644	10 x 10 cm (approx.)	8 x 60 pcs	2 packs	12349216
32-ply	Pack of 5 pcs	4636	12.5 x 20 cm (approx.)	12 x 5 pcs	4 packs	12349239

# BeeSana<sup>®</sup> Slit gauze swab

### Our gauze swabs with Y-slit

### Features

- gauze swab, medical gauze (EN 14079 standard)
- with Y-slit
- selvedge cut edges
- folded, 17-thread





### Usage

For wound care for drainage tubes, cannulas, tracheotomies and tractions.

### Non-sterile

	REF	Size (width x length)	Pack contents	Shipping unit	PZN
8-ply	4291	7.5 x 7.5 cm (approx.)	100 pcs	10 packs	12491084
12-ply	4192	5 x 5 cm (approx.)	100 pcs	10 packs	12491090
	4122	7.5 x 7.5 cm (approx.)	100 pcs	10 packs	12491109
	4160	10 x 10 cm (approx.)	100 pcs	10 packs	12491115
16-ply	4159	5 x 5 cm (approx.)	100 pcs	10 packs	12491121





### Single sterile, pack of 2 pcs

8-ply	4158	10 x 10 cm (approx.)	100 x 2 pcs	8 packs	03870954
12-ply	4115	5 x 5 cm (approx.)	25 x 2 pcs	16 packs	01916669
	4116	7.5 x 7.5 cm (approx.)	25 x 2 pcs	16 packs	01916652
	4117	10 x 10 cm (approx.)	25 x 2 pcs	16 packs	01916646

# BeeSana® Non-woven swab

### Our non-woven swab for external wound care

### Features

- viscose and polyester non-woven swab
- absorbent
- virtually no fraying



### Usage

For general-purpose wound care of exuding wounds and for absorbing bodily fluids. For applying ointment in smaller-scale outpatient surgical interventions, in wards and for first aid applications.

### Non-sterile

		REF	Size (width x length)	Pack contents	Shipping unit	PZN
<b>4-ply</b> 30 g		1789	5 x 5 cm (approx.)	100 pcs	10 packs	01450254
	_	1790	7.5 x 7.5 cm (approx.)	100 pcs	10 packs	01450225
	-	1791	10 x 10 cm (approx.)	100 pcs	10 packs	01450231
	-	1792	10 x 20 cm (approx.)	100 pcs	5 packs	01445773
	10 pcs, pre-counted For set packers	2791	10 x 10 cm (approx.)	200 x 10 pcs	1 box	_
40 g		1124	5 x 5 cm (approx.)	100 pcs	10 packs	04852881
	_	1125	7.5 x 7.5 cm (approx.)	100 pcs	10 packs	04852898
	_	1126	10 x 10 cm (approx.)	100 pcs	10 packs	04852906
	Large box	1127	10 x 10 cm (approx.)	100 pcs	76 packs	_
	-	11261	10 x 20 cm (approx.)	100 pcs	5 packs	04852912
<b>6-ply</b> 30 g		1796	5 x 5 cm (approx.)	100 pcs	10 packs	01330886
	-	1797	7.5 x 7.5 cm (approx.)	100 pcs	10 packs	01330900
	-	1798	10 x 10 cm (approx.)	100 pcs	10 packs	01330917
	_	1799	10 x 20 cm (approx.)	100 pcs	5 packs	01330923
	10 pcs, pre-counted	2798	10 x 10 cm (approx.)	200 x 10 pcs	1 box	_
	For set packers	2799	10 x 20 cm (approx.)	100 x 10 pcs	1 box	_

### Single sterile

<b>4-ply</b> 30 g	Pack of 2 pcs	1137	5 x 5 cm (approx.)	25 x 2 pcs	27 packs	03036843
	-	1138	7.5 x 7.5 cm (approx.)	25 x 2 pcs	16 packs	03036895
	-	1139	10 x 10 cm (approx.)	25 x 2 pcs	16 packs	03036903
		1140	10 x 20 cm (approx.)	25 x 2 pcs	8 packs	03036949
	Pack of 5 pcs	4376	7.5 x 7.5 cm (approx.)	100 x 5 pcs	8 packs	12350768
		4377	10 x 10 cm (approx.)	60 x 5 pcs	8 packs	12350774
	Pack of 10 pcs	1219	10 x 10 cm (approx.)	25 x 10 pcs	8 packs	12350691
<b>6-ply</b> 30 g	Pack of 2 pcs	1800	5 x 5 cm (approx.)	25 x 2 pcs	27 packs	01330834
		1801	7.5 x 7.5 cm (approx.)	25 x 2 pcs	16 packs	01330840
		1802	10 x 10 cm (approx.)	25 x 2 pcs	16 packs	01330857
		1803	10 x 20 cm (approx.)	25 x 2 pcs	8 packs	03036955
	Pack of 5 pcs	4225	5 x 5 cm (approx.)	120 x 5 pcs	8 packs	12350662
		4226	7.5 x 7.5 cm (approx.)	70 x 5 pcs	8 packs	12350679
	-	4227	10 x 10 cm (approx.)	40 x 5 pcs	8 packs	12350685

# BeeSana<sup>®</sup> Slit non-woven swab

### Our non-woven swab with Y-slit

#### Features

- viscose and polyester non-woven dressing swab
- with Y-slit
- absorbent
- virtually no fraying



#### Usage

For wound care for drainage tubes, cannulas, tracheotomies and tractions.

### Non-sterile

	REF	Size (width x length)	Pack contents	Shipping unit	PZN
4-ply, 30 g	1793	5 x 5 cm (approx.)	100 pcs	10 packs	01330768
	1794	7.5 x 7.5 cm (approx.)	100 pcs	10 packs	01330774
	1795	10 x 10 cm (approx.)	100 pcs	10 packs	01330780



### Single sterile, pack of 2 pcs

4-ply, 30 g	1232	5 x 5 cm (approx.)	25 x 2 pcs	27 packs	01364684
	1234	7.5 x 7.5 cm (approx.)	25 x 2 pcs	16 packs	01364678
	1233	10 x 10 cm (approx.)	25 x 2 pcs	16 packs	01364661
6-ply, 30 g	4105	10 x 10 cm (approx.)	25 x 2 pcs	28 packs	13507290

## **ABE**<sup>®</sup> drainal

### Our slit/hole dressing pad with metallic coating

#### Features

- non-woven slit dressing pad with hole
- metallic coating lessens the risk of the material adhering to the wound
- hole diameter: 5.5–11 mm (approx.)



#### Usage

For external wound care with drainage tubes and infusion catheters. Place the metallic coating against the skin.

Single sterile, pack of 1 pc.						
REF	Size (width x length)	Pack contents	Shipping unit	PZN		
4118	6 x 7 cm (approx.)	50 x 1 pc.	6 packs	01333867		

## **ABE**<sup>®</sup> tracheal

## Our slit/hole dressing pads with metal coating

#### Features

- non-woven slit dressing pad with hole
- metallic coating lessens the risk of the material adhering to the wound
- sterile packed
- hole diameter: 13-30 mm (approx.)



#### Usage

Slit/hole dressing pad for would care with catheters, stomata and drainage tubes. The hole in the pad also has a slit, so the ABE® Tracheal can create a perfect seal around a drainage tube or stoma, and the exudate can be quickly absorbed by the the non-woven dressing pad. Place the metallic coating against the skin.

Single	sterile, pack of 1 p			
REF	Size (width x length)	Pack contents	Shipping unit	PZN
4119	8 x 10 cm (approx.)	50 x 1 pc.	6 packs	_

## **BeeSana**<sup>®</sup> Absorbent dressing pad

Our absorbent pad for superior fluid absorption

#### Features

- wound-friendly absorbent dressing pad
- absorbent core made from 100% cellulose
- fully enveloped in non-woven material
- blue top surface for identifying the reverse side and as a hydrophobic barrier layer to avoid soaking through, and protect underwear and clothing from contamination

#### Usage

For wound care of exuding wounds in all areas of healthcare. Absorbent dressing pads are also ideal for padding wounds.







### Non-sterile, in polybag

REF	Size (width x length)	Pack contents	Shipping unit	PZN
1227	10 x 10 cm (approx.)	25 pcs	16 packs	01363816
1228	10 x 20 cm (approx.)	25 pcs	8 packs	07607202
1224	10 x 30 cm (approx.)	25 pcs	15 packs	12146619
1225	15 x 20 cm (approx.)	25 pcs	12 packs	12146625
1231	15 x 25 cm (approx.)	25 pcs	7 packs	00137727
1229	20 x 20 cm (approx.)	25 pcs	12 packs	08696337
1244	20 x 30 cm (approx.)	25 pcs	5 packs	12146631
1230	20 x 40 cm (approx.)	25 pcs	6 packs	08442499

## Single sterile, shrink-wrapped pack of 1 pc.

1354	10 x 10 cm (approx.)	60 x 1 pc.	10 packs	08416059
1355	10 x 20 cm (approx.)	30 x 1 pc.	12 packs	06193514
1388	15 x 25 cm (approx.)	30 x 1 pc.	8 packs	12146602
1356	20 x 20 cm (approx.)	30 x 1 pc.	6 packs	08416042
1357	20 x 40 cm (approx.)	25 x 1 pc.	4 packs	07737904

## **BeeSana**<sup>®</sup> Eye pad

Our eye dressing pad with cushioning effect

### Features

- oval dressing pad with all-round seal
- absorbent inner filling



#### Usage

For absorbing fluids from the eye. For padding and for eye bandages. Important: For injuries to the external eye and for postoperative applications, a sterile eye dressing pad/wound covering must be used. The oval design is ideal for fitting the eye socket.

## Non-sterile

REF	Size (width x length)	Pack contents	Shipping unit	PZN
4195	5.6 x 7.2 cm (approx.)	50 pcs	10 packs	02214099
4155	5.6 x 7.2 cm (approx.)	500 pcs	1 pack	12348872
4194	7 x 8.5 cm (approx.)	50 pcs	10 packs	02435297

### Single sterile, shrink-wrapped pack of 1 pc

4196	5 x 7 cm (approx.)	25 x 1 pc.	9 packs	02214107
4156	7 x 8.5 cm (approx.)	25 x 1 pc.	6 packs	02214136

## **ABE**<sup>®</sup>nose bandage

## Our flexible nasal sling

#### Features

- bi-elastic nose bandage
- with elastic ear loops or an elastic strap
- with integrated dressing pad in centre





#### Usage

For post-operative ENT applications with exudate. Worn with ear loops or by being tied at the back of the head.

Non-sterile, with elastic strap						
REF	Pack contents	Shipping unit	PZN			
4039	20 pcs	10 packs	01333749			

Non-sterile, with ear loops						
REF	Pack contents	Shipping unit	PZN			
9814	20 pcs	10 packs	11885562			
9818	500 pcs	1 pack	13579533			

## **BeeSana**<sup>®</sup> Gauze balls

## Our super-absorbent gauze balls

#### Features

- gauze balls, medical gauze (EN 14079 standard)
- innovative wrapping system from just one piece of pre-cut gauze
- inner-facing cut edges
- 20-thread
- absorbent

#### Usage

For cleaning wounds and staunching bleeding. Available in various swab sizes to suit each wound. Non-sterile,

## Non-sterile, in polybag

	REF	Size (width x length)	Shape	Pack contents	Shipping unit	PZN
Hazelnut-sized	1102	10 x 10 cm (approx.)	Round	500 pcs	8 packs	01364017
Walnut-sized	4460	15 x 16 cm (approx.)	Round	500 pcs	2 packs	00822216
Plum-sized	4461	20 x 20 cm (approx.)	Round	250 pcs	16 packs	00822239
Plum-sized	1106	20 x 20 cm (approx.)	Thimble	250 pcs	16 packs	00822222
Egg-sized	4462	25 x 26 cm (approx.)	Round	250 pcs	10 packs	00822251
Egg-sized	1109	25 x 26 cm (approx.)	Thimble	250 pcs	10 packs	00822245
X-sized	1348	30 x 30 cm (approx.)	Round	100 pcs	20 packs	11862578
XX-sized	1113	30 x 40 cm (approx.)	Round	100 pcs	15 packs	01363940
Fist-sized	1116	40 x 50 cm	Round	100 pcs	10 packs	00822268

## Single sterile, in peel pack

Pack of 2 pcs	Plum-sized	4356	20 x 20 cm (approx.)	Round	100 x 2 pcs	8 packs	13506801
Pack of 3 pcs	Plum-sized	4357	20 x 20 cm (approx.)	Round	80 x 3 pcs	8 packs	00137696
	Egg-sized	4308	25 x 26 cm (approx.)	Round	60 x 3 pcs	8 packs	13506818
	Fist-sized	4085	45 x 45 cm (approx.)	Round	35 x 3 pcs	4 packs	13506824
2+2 pcs per pack	Plum-sized	1084	20 x 20 cm (approx.)	Round	50x (2+2) pcs	8 packs	01364129
Packed in soft peel pack with 2 separate compartments (tear-off perforation)	Plum-sized	1164	20 x 20 cm (approx.)	Thimble	50x (2+2) pcs	8 packs	01364106
2+3 pcs per pack	Walnut-sized	1074	15 x 16 cm (approx.)	Round	45x (2+3) pcs	8 packs	13506830
Packed in soft peel pack with — 2 separate compartments	Plum-sized	1131	20 x 20 cm (approx.)	Round	40x (2+3) pcs	8 packs	01364135
	Plum-sized	1189	20 x 20 cm (approx.)	Thimble	40x (2+3) pcs	8 packs	01364075
(tear-off perforation)	Egg-sized	1141	25 x 26 cm (approx.)	Round	35x (2+3) pcs	8 packs	-
	XX-sized	4216	30 x 40 cm (approx.)	Round	35x (2+3) pcs	4 packs	13506847
Pack of 5 pcs	Plum-sized	1184	20 x 20 cm (approx.)	Round	50 x 5 pcs	8 packs	13506853
	Egg-sized	4209	25 x 26 cm (approx.)	Round	40 x 5 pcs	8 packs	13506876
-	Egg-sized	1165	25 x 26 cm (approx.)	Thimble	40 x 5 pcs	8 packs	01364069
	XX-sized	4211	30 x 40 cm (approx.)	Round	40 x 5 pcs	4 packs	13506882
Pack of 10 pcs	Plum-sized	4215	20 x 20 cm (approx.)	Round	30 x 10 pcs	8 packs	07393706

## **BeeSana**<sup>®</sup> Pointed swabs

Our pointed swabs made from pure cotton

#### Features

- gauze swab, 100% cotton (EN 14079 standard)
- 20-thread
- tapered, formed from just one piece of gauze



#### Usage

For blood and exudate absorption, especially in ENT applications. For staunching bleeding or as a tamponade.

Double ste	rile, with a	counting cards				
	REF	Size (unrolled) (width x length)	Length	Pack contents	Shipping unit	PZN
Pack of 10 pcs	4170	Large, 15 x 16 cm (approx.)	13.5 cm (approx.)	40 x 10 pcs	8 packs	_

## **ABE**<sup>®</sup> ribbon gauze

### Our tamponade made from pure cotton

#### Features

- medical gauze (EN 14079 standard)
- 100% cotton
- woven edges
- white
- rolled
- woven-in x-ray contrast threads (BaSO<sub>4</sub>)

#### Usage

As a tamponade in natural bodily orifices.

## Single sterile, pack of 1 pc.

	REF	Size (width x length)	Pack contents	Shipping unit	PZN
24-thread	4300	1 cm x 5 m (approx.)	100 x 1 pc.	8 packs	13574814
	4301	2 cm x 5 m (approx.)	70 x 1 pc.	8 packs	13574820



## BeeSana<sup>®</sup> Cellulose swabs

## Our cellulose swab on a roll

#### Features

- swab made from highly bleached cellulose (German pharmacopoeia (DAB) standard)
- ready-to-use
- pre-cut





#### Usage

For local disinfection before injections and infusions.

### Non-sterile

	REF	Size (width x length)	Pack contents	Shipping unit	PZN
	1190	4 x 5 cm (approx.)	2 rolls x 500 pcs	16 packs	07210048
	1278	8 x 10 cm (approx.)	2 rolls x 125 pcs	10 packs	01364253
Reduced bioburden	1177	4 x 5 cm (approx.)	1 roll x 500 pcs	20 packs	01331070
Cellulose swab box	9871		1 pc.	1 pack	01363756

**Alcohol swabs** 

### Our ready-to-use swabs for skin disinfection

#### Features

- non-woven alcohol swabs
- impregnated with 70% isopropyl alcohol
- ready-to-use

#### Usage

For disinfecting skin before subcutaneous injections and capillary blood sampling.



REF	Size (width x length)	Pack contents	Shipping unit	PZN
4465	Approx. 3 x 3 cm (folded) Approx. 6 x 3 cm (unfolded)	100 pcs	20 packs	09391893

## Cotton buds

#### Features

- various designs and variations for different types of applications
- sterile and non-sterile formats
- various head sizes
- wood or flexible plastic variants



#### Usage

For a wide range of uses in medical applications, such as cleaning, applying liquid substances and ointments as well as for laboratory work.

## Non-sterile, in Ziploc bag

REF	Length	Head	Handle	Pack contents	Shipping unit	PZN
4046	7.5 cm	Small head	Plastic	10 x 200 pcs	10 packs	11857620
1260	15 cm	Small head	Wood	10 x 100 pcs	20 packs	06177225
1261	15 cm	Medium head	Wood	10 x 100 pcs	10 packs	06177231
1263	15 cm	Large head	Wood	10 x 50 pcs	10 packs	01363762
1268	15 cm	Large head	Plastic	10 x 50 pcs	10 packs	11857494

## Single sterile, 2 pcs per pack

1262	15 cm	Small head	Wood	100 x 2 pcs	6 packs	06177254
1264	15 cm	Large head	Wood	100 x 2 pcs	3 packs	01363779
1270	15 cm	Large head	Plastic	100 x 2 pcs	3 packs	13711168

## Wooden tongue depressor

Our tongue depressor for examinations and for applying ointment



Non-s	terile					
REF	Size (width x length)	Pack contents	Shipping unit	PZN		
9802	150 x 18 mm (approx.)	100 pcs	20 packs	00822334		
Single	Single sterile, 1 pc. per pack					
9665	150 x 18 mm (approx.)	100 х 1 рс.	10 packs	01333821		

## **ABE**<sup>®</sup> injection plaster

## Our fabric injection plaster

#### Features

- ready-to-use injection plaster
- with wound pad
- polyester net lessens risk of dressing pad becoming stuck to the wound
- with polyacrylate adhesive
- woven fabric structure



#### Usage

For minor injuries such as scrapes/grazes and shallow cuts, and for wound care after taking blood.

## Non-sterile

REF	Size (width x length)	Pack contents	Shipping unit	PZN
9525	2 x 4 cm (approx.)	500 pcs	20 packs	01664860

## **ABE**<sup>®</sup> light

## Our elastic wound bandage with wound dressing pad

#### Features

- white, slightly elastic wound bandage
- pad placed in the middle
- polyester net lessens risk of dressing pad becoming stuck to the wound
- with polyacrylate adhesive



#### Usage

For first aid/follow-up treatment of minor injuries, also on sensitive parts of the body. Good fit for wounds to joints.

REF	Size (width x length)	Pack contents	Shipping unit	PZN
9547	4 cm x 5 m (approx.)	1 pc.	14 packs	06790233
9548	6 cm x 5 m (approx.)	1 pc.	10 packs	06790256
9549	8 cm x 5 m (approx.)	1 pc.	8 packs	06790262

## **ABE**<sup>®</sup> Absorbent adhesive dressing

## Our absorbent and hypoallergenic instant wound dressing

### Features

- self-adhesive wound dressing made from white non-woven material
- with polyacrylate adhesive
- with central wound pad
- less risk of sticking to the wound



### Usage

For post-operative wound care or in the ward as a sterile instant wound bandage for absorption of blood and exudate, and to protect wounds from (cross-)infection.

## Single sterile, shrink-wrapped pack of 1 pc.

REF	Size (width x length)	Pack contents	Shipping unit	PZN
4750	5 x 7 cm (approx.)	100 x 1 pc.	12 packs	01452997
4752	6 x 10 cm (approx.)	50 x 1 pc.	6 packs	01126022
4754	8 x 10 cm (approx.)	50 x 1 pc.	6 packs	01446011
4760	15 x 10 cm (approx.)	50 x 1 pc.	6 packs	01445980
4762	20 x 10 cm (approx.)	50 x 1 pc.	6 packs	01445968
4764	25 x 10 cm (approx.)	25 x 1 pc.	6 packs	01445951
4768	35 x 10 cm (approx.)	25 x 1 pc.	6 packs	01446057

## **ABE**<sup>®</sup> derm

### Transparent film without wound dressing

#### Features

- elastic, semi-permeable polyurethane film
- coated with polyacrylate adhesive
- wound environment visible through the transparent film



#### Usage

Suitable for postoperative wound care and for monitoring the wound environment, wound care and catheter fixation. Can also be placed directly on primary wounds for wound monitoring.

### Single sterile, shrink-wrapped pack of 1 pc.

REF	Size (width x length)	Pack contents	Shipping unit	PZN
9571	6 x 7 cm (approx.)	100 х 1 рс.	40 packs	06790204
9572	10 x 12 cm (approx.)	50 x 1 pc.	12 packs	06790210
9573	15 x 20 cm (approx.)	50 x 1 pc.	15 packs	06790227

## **ABE**<sup>®</sup> derm plus Transparent film wound dressing

#### Features

- elastic, semi-permeable polyurethane film
- coated with polyacrylate adhesive
- integrated wound pad
- wound environment visible through the transparent film



#### Usage

Suitable for postoperative wound care and for monitoring the wound environment, wound care and catheter fixation. Can also be placed directly on primary wounds for wound monitoring.

Single st				
REF	Size (width x length)	Pack contents	Shipping unit	PZN
9580	5 x 7 cm (approx.)	50 x 1 pc.	24 packs	04589917
9581	6 x 10 cm (approx.)	50 x 1 pc.	12 packs	04589923
9584	15 x 10 cm (approx.)	50 x 1 pc.	12 packs	04656720
9585	20 x 10 cm (approx.)	50 x 1 pc.	24 packs	04656737

## **ABE**<sup>®</sup> cannula retention dressing

## Our cannula fixing plaster

### Features

- self-adhesive fixing dressing
- with central slit
- less risk of the integrated wound pad sticking to the wound
- coated with polyacrylate adhesive

### Usage

For sterile wound dressing of the puncture site and the fixing of venous access ports.

Single st				
REF	Size (width x length)	Pack contents	Shipping unit	PZN
4457	8 x 6 cm (approx.)	50 x 1 pc.	10 packs	12589854

## **ABE**<sup>®</sup> non-woven fabric

## Our self-adhesive non-woven fixing material

### Features

- self-adhesive non-woven fixing material, white
- stretchable widthways (min. 100%)
- sess risk of the integrated wound pad sticking to the wound
- coated with polyacrylate adhesive



#### Usage

For covering and fixing large bandaged areas in medical applications. For securing cannulas, drainage tubes, probes, catheters, measuring instruments, etc. For securing bandages to joints, and other round/tapered parts of the body.

Non-s	terile			
REF	Size (width x length)	Pack contents	Shipping unit	PZN
4470	5 cm x 10 m (approx.)	1 pc.	50 packs	06168663
4471	10 cm x 10 m (approx.)	1 pc.	24 packs	06168686
4472	15 cm x 10 m (approx.)	1 pc.	24 packs	06168692
4473	20 cm x 10 m (approx.)	1 pc.	52 packs	06168700



ditrade

## **ABE**<sup>®</sup> silk

## Our synthetic silk plaster on a roll

#### Features

- fixing plaster made from white cellulose acetate
- coated with polyacrylate adhesive
- zigzag edge, easily torn off widthways or lengthways



#### Usage

For fixing wound coverings, bandages, drainage tubes, catheters and cannulas in place.

## Non-sterile, with protective ring

REF	Size (width x length)	Pack contents	Shipping unit	PZN
9531	2.5 cm x 5 m (approx.)	24 pcs	10 packs	11732887
Non-sterile, without protective ring				
9535	1.25 cm x 5 m (approx.)	36 pcs	10 packs	11734544
9536	2.5 cm x 5 m (approx.)	24 pcs	10 packs	11745536
9540	1.25 cm x 9.14 m (approx.)	36 pcs	10 packs	03916834
9541	2.5 cm x 9.14 m (approx.)	24 pcs	10 packs	03916840

## BeeSana<sup>®</sup> SafeCut

## Our self-adhering finger plaster bandage

### Features

- self-adhering, elastic fixation bandage
- fine latex coating prevents slippage
- stretches by approx. 100%







### Usage

For personal first aid for finger and joint injuries, and for torn fingernails. For fixing wound coverings, and for use as a light support and compression bandage. For immobilising injured parts of the body and limbs.

### Non-sterile

REF	Size (width x length, stretched)	Colour	Pack contents	Shipping unit	PZN
1844	2.5 cm x 4.5 m (approx.)	Blue	1 pc.	252 packs	10544882

### Fold the support pad (fixation possible with or without pad)



**1.** Fold the end of the bandage over itself a few times until the support pad created is adequate.



**2.** Fix the support pad onto the intended location with the rest of the bandage until the bandage is stable enough.



**3.** Apply the bandage firmly and then tear off.

## **ABE**<sup>®</sup> gauze bandage Our traditional fixing bandage

#### Features

- rigid fixing bandage (DIN 61631 standard)
- smooth structure
- woven edges
- 65% viscose/35% cotton



#### Usage

Broad field of use for medical applications, e.g. as an emergency bandage for first aid or as a carrier material for zinc paste bandages.

REF	Size (width x length)	Pack contents	Shipping unit	PZN
1400	4 cm x 4 m (approx.)	20 pcs	16 packs	01331029
1401	6 cm x 4 m (approx.)	20 pcs	16 packs	01331035
1402	8 cm x 4 m (approx.)	20 pcs	12 packs	01331041
1403	10 cm x 4 m (approx.)	20 pcs	8 packs	01331058

## **ABE<sup>®</sup> Last**

## Our elastic fixation bandage with a smooth structure

### Features

- elastic fixing bandage (DIN 61634 standard)
- fine, smooth structure
- fully stretched length 4 m (approx.)
- 54% polyester/46% viscose



### Usage

For fixing wound coverings, for fixing support splints, and for use as a light support and compression bandage. Suitable for use on joints, tapered and rounded body parts. For immobilising injured parts of the body and limbs.

Non-sterile, pack of 20 pcs.				
REF	Size (width x length, stretched)	Pack contents	Shipping unit	PZN
1409	4 cm x 4 m (approx.)	20 pcs	28 packs	04853248
1410	6 cm x 4 m (approx.)	20 pcs	24 packs	01364282
1411	8 cm x 4 m (approx.)	20 pcs	20 packs	04853260
1412	10 cm x 4 m (approx.)	20 pcs	15 packs	04853277
1413	12 cm x 4 m (approx.)	20 pcs	6 packs	04853283

## **ABE<sup>®</sup> Last** crepe

## Our elastic fixation bandage with a crepe structure

#### Features

- elastic fixation bandage (DIN 61634 standard)
- crepe structure
- adhesion of individual layers
- fully-stretched length 4 m (approx.)
- 70% viscose/30% polyamide



#### Usage

For fixing wound coverings, for fixing support splints, and for use as a light support and compression bandage. Suitable for use on joints, tapered and rounded body parts.

Non-st	terile			
REF	Size (width x length, stretched)	Pack contents	Shipping unit	PZN
1749	4 cm x 4 m (approx.)	40 pcs	12 packs	14439863
1750	6 cm x 4 m (approx.)	20 pcs	16 packs	07198324
1751	8 cm x 4 m (approx.)	20 pcs	12 packs	07198330
1752	10 cm x 4 m (approx.)	20 pcs	10 packs	07198347
1753	12 cm x 4 m	20 pcs	8 packs	07198353

## Double sterile, pack of 1 pc. with 2 counting cards

4741	6 cm x 4 m (approx.)	30 x 1 pc.	4 packs	13655731
4742	8 cm x 4 m (approx.)	30 x 1 pc.	4 packs	13655642
4588	12 cm x 4 m (approx.)	30 x 1 pc.	4 packs	13655725

## **ABE**<sup>®</sup> fix color

## Our elastic and self-adhering coloured fixing bandage, smooth

#### Features

- self-adhering and elastic fixation bandage
- smooth structure
- distensibility (DIN 61632 standard) min. 40%
- 61% polyamide/39% cotton



#### Usage

For fixing wound coverings and bandages on body parts that are moved frequently. Such as for isolating wrist or ankle joints, for example.

REF	Size (width x length, stretched)	Colour	Pack contents	Shipping unit	PZN
19461	6 cm x 20 m (approx.)	Blue	1 pc.	15 packs	08520615
19463	6 cm x 20 m (approx.)	Green	1 pc.	36 packs	08520638
19465	6 cm x 20 m (approx.)	Red	1 pc.	15 packs	08520650

## **ABE**<sup>®</sup> Last adhesive

## Our elastic and self-adhering fixing bandage

#### Features

- self-adhering, elastic fixation bandage
- crepe structure
- no additional fixation required
- microdot latex impregnation
- large-pored
- stretches by approx. 60%
- 70% polyester/30% polyamide





#### Usage

For fixing wound coverings, for fixing support splints, and for use as a light support and compression bandage. Suitable for use on joints, tapered and rounded body parts. For immobilising injured parts of the body and limbs.

## Non-sterile

REF	Size (width x length, stretched)	Pack contents	Shipping unit	PZN
1289	4 cm x 4 m (approx.)	1 pc.	30 packs	04019746
1391	6 cm x 4 m (approx.)	1 pc.	36 packs	07198436
1392	8 cm x 4 m (approx.)	1 pc.	36 packs	07198442
1393	10 cm x 4 m (approx.)	1 pc.	36 packs	07198459
1394	12 cm x 4 m (approx.)	1 pc.	36 packs	07198465
1288	4 cm x 20 m (approx.)	1 pc.	22 packs	02435498
1383	6 cm x 20 m (approx.)	1 pc.	16 packs	04853308
1384	8 cm x 20 m (approx.)	1 pc.	12 packs	04853314
1385	10 cm x 20 m (approx.)	1 pc.	10 packs	04853320
1386	12 cm x 20 m (approx.)	1 рс.	10 packs	04853337

### Double sterile, pack of 1 pc. with counting card

4793 12 cm x 4 m 15 x 1 pc. 4 packs –

## ABE<sup>®</sup> Ideal 100

## Our short-stretch elastic fabric Ideal bandage

#### Features

- fabric elastic bandage (DIN 61632 standard)
- natural white
- 100% cotton
- fully stretched length 5 m (approx.)
- selvedge
- washable



### Usage

For relief bandages and follow-up treatment of dislocations or fractures.

## Non-sterile

REF	Size (width x length, stretched)	Pack contents	Shipping unit	PZN
1589	4 cm x 5 m (approx.)	10 pcs	10 packs	01364299
1423	6 cm x 5 m (approx.)	10 pcs	10 packs	01364477
1424	8 cm x 5 m (approx.)	10 pcs	10 packs	01364448
1425	10 cm x 5 m (approx.)	10 pcs	5 packs	04853484
1426	12 cm x 5 m (approx.)	10 pcs	5 packs	04853490
1576	15 cm x 5 m (approx.)	10 pcs	5 packs	01364460
1577	20 cm x 5 m (approx.)	10 pcs	5 packs	01364454

## Double sterile, pack of 1 pc. with counting card

4771	6 cm x 5 m (approx.)	25 x 1 pc.	4 packs	-
4772	8 cm x 5 m (approx.)	25 x 1 pc.	4 packs	_
4773	10 cm x 5 m (approx.)	20 x 1 pc.	4 packs	_
4774	12 cm x 5 m (approx.)	15 x 1 pc.	4 packs	_

## **ABE**<sup>®</sup> universal

## Our permanently-elastic, medium-stretch universal bandage

#### Features

- permanently elastic bandage for medium compression
- open, air-permeable structure
- fully stretched length: approx. 5 m
- latex-free
- 67% cotton/30% polyester/3% elastane



#### Usage

For bandages where medium compression is required, and especially for compressive wound bandages. For fixing wound coverings. For immobilising injured parts of the body and limbs, and for supporting and relieving pressure on limbs.

### Non-sterile

REF	Size (width x length, stretched)	Pack contents	Shipping unit	PZN
1414	6 cm x 5 m (approx.)	10 pcs	16 packs	00761561
1415	8 cm x 5 m (approx.)	10 pcs	12 packs	00761590
1416	10 cm x 5 m (approx.)	10 pcs	10 packs	00761526
1432	12 cm x 5 m (approx.)	10 pcs	8 packs	00761532
1433	15 cm x 5 m (approx.)	10 pcs	6 packs	00761549
1434	20 cm x 5 m (approx.)	10 pcs	5 packs	00761555

## Double sterile, pack of 1 pc. with counting card

4722	6 cm x 5 m (approx.)	25 x 1 pc.	4 packs	13655636
4723	8 cm x 5 m (approx.)	25 x 1 pc.	4 packs	13655694
4572	10 cm x 5 m (approx.)	20 x 1 pc.	4 packs	13655748
4560	12 cm x 5 m (approx.)	15 x 1 pc.	4 packs	13655719

## **ABE**<sup>®</sup> lastic

## Our permanently elastic, universal short-stretch bandage

#### Features

- permanently elastic bandage for medium compression
- open, air-permeable structure
- fully stretched length: approx. 5 m
- latex-free
- 60% cotton/38% polyester/2% elastane



#### Usage

For bandages where light compression is required, and especially for compressive wound support and pressure relief bandages. For immobilising injured parts of the body and limbs. For fixing wound coverings.

### Non-sterile

REF	Size (width x length, stretched)	Pack contents	Shipping unit	PZN
1780	4 cm x 5 m (approx.)	10 pcs	10 packs	02423265
1781	6 cm x 5 m (approx.)	10 pcs	10 packs	02423259
1782	8 cm x 5 m (approx.)	10 pcs	10 packs	02423242
1783	10 cm x 5 m (approx.)	10 pcs	10 packs	02423236
1784	12 cm x 5 m (approx.)	10 pcs	5 packs	02423213
1785	15 cm x 5 m (approx.)	10 pcs	5 packs	02423207
1786	20 cm x 5 m (approx.)	10 pcs	5 packs	02423182

## Double sterile, pack of 1 pc. with counting card

4779	6 cm x 5 m (approx.)	30 x 1 pc.	4 packs	13655671
4780	8 cm x 5 m (approx.)	25 x 1 pc.	4 packs	13655659
4781	10 cm x 5 m (approx.)	20 x 1 pc.	4 packs	13655702
4782	12 cm x 5 m (approx.)	15 x 1 pc.	4 packs	13655665

## **ABE**<sup>®</sup> padded bandage

## Our synthetic wadding for padding

#### **Features**

- synthetic wadding bandage made from polyester fibre
- latex-free
- can be sterilised

#### Usage

As a padding material under plaster casts, and support/compression bandages.

### Non-sterile

REF	Size (width x length)	Pack contents	Shipping unit	PZN	
1701	6 cm x 3 m (approx.)	1 pc.	108 packs	06627958	
1806	8 cm x 3 m (approx.)	1 pc.	72 packs	04351788	
1702	10 cm x 3 m (approx.)	1 pc.	72 packs	06627964	
1703	15 cm x 3 m (approx.)	1 pc.	36 packs	06627970	
Single sterile, pack of 1 pc.					

1705	10 cm x 3 m (approx.)	15 x 1 pc.	4 packs	13715901
------	-----------------------	------------	---------	----------

## **ABE**<sup>®</sup> extra

## Our short-stretch bandage for support and compression bandages

#### **Features**

- powerful, fabric elastic compression bandage
- short-stretch
- open, air-permeable structure
- fully stretched length 5 m (approx.)
- 100% cotton
- with bandage clips

#### Usage

For strong compression. Suitable for opposing bandaging techniques, according to Pütter. For stabilising, relief and fixation during the treatment of injuries to the musculoskeletal system.

REF	Size (width x length, stretched)	Pack contents	Shipping unit	PZN
1421	8 cm x 5 m (approx.)	1 pc.	72 packs	01364402
1422	10 cm x 5 m (approx.)	1 pc.	36 packs	01364394
1735	6 cm x 5 m (approx.)	10 pcs	10 packs	08490759
1736	8 cm x 5 m (approx.)	10 pcs	10 packs	08490765
1737	10 cm x 5 m (approx.)	10 pcs	5 packs	08490771
1738	12 cm x 5 m (approx.)	10 pcs	5 packs	08490788





## **ABE**<sup>®</sup> extra adhesive

## Our self-adhering, permanently elastic short-stretch compression bandage

### Features

- strong, permanently elastic compression bandage
- short-stretch
- self-adhering
- selvedge
- can be sterilised with gas and steam at 121 °C, x-ray and gamma
- fully-stretched length 5 m (approx.)
- 88% cotton/11% polyamide/1% polyurethane

### Usage

For stabilising, relief and fixation during the treatment of injuries to the musculoskeletal system.

### Non-sterile

REF	Size (width x length, stretched)	Pack contents	Shipping unit	PZN
1777	8 cm x 5 m (approx.)	1 pc.	64 packs	01340933
1778	10 cm x 5 m (approx.)	1 pc.	32 packs	01340956
1779	12 cm x 5 m (approx.)	1 pc.	32 packs	01340979

## **ABE<sup>®</sup> Ven**

## Our permanently-elastic, long-stretch compression bandage

#### Features

- permanently elastic
- with latex-free end fixation
- 93% cotton/7% polyurethane
- fully stretched length 7 m (approx.)
- washable
- with bandage clips

### Usage

For stabilising, relief and fixation during the treatment of injuries to the musculoskeletal system.

REF	Size (width x length, stretched)	Pack contents	Shipping unit	PZN
1563	8 cm x 7 m (approx.)	1 pc.	64 packs	04854147
1564	10 cm x 7 m (approx.)	1 pc.	32 packs	04854153
1565	12 cm x 7 m (approx.)	1 pc.	32 packs	04854176





## **ABE**<sup>®</sup> tubular bandage

## Our fine mesh knitted fabric tube for universal use

#### Features

- tube bandage made from 100% cotton
- stretches lengthways and widthways
- easy to use
- can be cut anywhere to fit
- seamless to avoid creating pressure points



#### Usage

As a pressure-free bandage and an underlay or overlay to protect the skin for zinc paste/plaster bandages. For immobilising and traction bandages. As an overlay for splints.

REF	Size (width x length)	Usage	Pack contents	Shipping unit	PZN
5181	Size 1 = 1.5 cm x 20 m (approx.)	Finger/toe	1 рс.	24 packs	01331354
5182	Size 2 = 2.5 cm x 20 m (approx.)	Several fingers/toes	1 рс.	16 packs	01331360
5183	Size $3 = 4$ cm x 20 m (approx.)	Hand/forearm/child's arm	1 рс.	12 packs	01331377
5184	Size 4 = 6 cm x 20 m (approx.)	Upper arm/lower leg	1 рс.	28 packs	01331383
5185	Size 5 = 8 cm x 20 m (approx.)	Thigh	1 рс.	22 packs	01331420
5186	Size $6 = 10 \text{ cm x } 20 \text{ m (approx.)}$	Head/child's torso	1 рс.	18 packs	01331437
5187	Size 7 = 12 cm x 20 m (approx.)	Torso	1 pc.	16 packs	01331443

## **ABE**<sup>®</sup> tubular net bandage

## Our highly elastic tubular net bandage

#### Features

- elastic tube bandage
- simple to use
- can be cut anywhere to fit
- binding prevents bandage unravelling further
- seamless to avoid creating pressure points
- follows anatomical contours to ensure a good fit



#### Usage

For fixing bandages and padding materials. Scissors can be used to cut a hole in the mesh material, so as to permit patients to exclude parts of the body (such as the thumb with a hand bandage) from being wrapped in the bandage. Suitable for use on joints, tapered and rounded body parts.

REF	Size (width x length)	Pack contents	Shipping unit	PZN
5164	Size $3 = \text{forearm}$ (approx., 2.2 cm x 25 m)	1 pc.	10 packs	03065804
5166	Size 4 = upper arm (approx., 2.5 cm x 25 m)	1 pc.	8 packs	03065810
5170	Size $5 = \text{leg}$ (approx., 4 cm x 25 m)	1 pc.	5 packs	03065827
5172	Size $6 = head$ (approx., 5.5 cm x 25 m)	1 pc.	4 packs	03065833
5176	Size 7 = torso (approx., 6.2 cm x 25 m)	1 pc.	4 packs	03065856
5180	Size 8 = torso XXL (approx., 8 cm x 25 m)	1 pc.	3 packs	03065862

## **ABE**<sup>®</sup> Gilchrist's bandage

Our universal shoulder bandage

#### Features

• universal shoulder bandage



### Usage

For injuries to the shoulder and upper arm area, for immobilising or fixing the shoulder joint.

Non-ste	rile			
REF	Size	Pack contents	Shipping unit	PZN
1914	Universal	1 pc.	12 packs	01363437

## **ABE**<sup>®</sup> body bandage Our post-operative support bandage

#### Features

- support and relief bandage
- wide application padding
- broad, adjustable hook-and-pile fastener
- compression can be adjusted to therapeutic requirements





#### Usage

For post-operative support and treatment of fractured ribs.

REF	Size	(circumference)	Pack contents	Shipping unit	PZN
1890	Size 1	65–75 cm (approx.)	1 рс.	15 packs	01363727
1891	Size 2	75–85 cm (approx.)	1 рс.	15 packs	01363710
1892	Size 3	85–95 cm (approx.)	1 рс.	15 packs	01363704
1893	Size 4	95–105 cm (approx.)	1 рс.	15 packs	01363696
1894	Size 5	105–115 cm (approx.)	1 pc.	15 packs	01363609
1895	Size 6	115–125 cm (approx.)	1 pc.	15 packs	01363555
1904	Size 7	125–135 cm (approx.)	1 pc.	15 packs	02435357

## **ABE**<sup>®</sup> heart catheter pressure bandage

## Our cardiac catheter pressure bandage

#### Features

- adjustable hook-and-pile fastener
- compression can be adjusted to therapeutic requirements
- colour-coded





#### Usage

For compression after femoral artery/vein punctures for invasive diagnostic work/treatment.

### Non-sterile

REF	Size (universal)	Pack contents	Shipping unit	PZN
2907	Left (colour code: red)	1 pc.	12 packs	04351742
2908	Right (colour code: blue)	1 pc.	12 packs	04351759

## **ABE**<sup>®</sup> limb holder

### Our cuff for adults

#### Features

- cuff from special, soft material
- hypoallergenic, thanks to cotton contact area
- hook-and-pile fastener for simple fixation
- no additional fixing strap required, since already integrated into the cuff itself

#### Usage

For short-term immobilisation.

# Non-sterileREFSizePack contentsShipping unitPZN1879For adults10 pcs12 packs12553046



## **ABE**<sup>®</sup> pressure bandage roll Special bandage

### Features

- pressure bandage roll made from polyurethane, in knitted fabric tube
- diameter 7 cm (approx.)
- latex-, PVC- and DEHP-free
- disposable



### Usage

As a pressure padding, for example after heart catheterisation.

# Non-sterile REF Size Pack contents Shipping unit PZN 4037 7 cm Ø x 10 cm (approx.) 12 pcs 9 packs 06620850

## Scrotum roll

#### Features

- scrotum roll from foamed material
- with a knitted fabric tube covering

### Usage

For elevating the scrotum after operations or with oedemas.

Non-steril	е

REF	Size	Pack contents	Shipping unit	PZN
4035	Ø 8 cm x 10 cm (approx.)	6 pcs	8 packs	01333755





## BeeSana<sup>®</sup> Lap sponge X-ray

Our abdominal swabs in a wide range of variants

#### Features

- pre-washed abdominal swabs from medical gauze (EN 14079 standard)
- with X-ray contrast chip
- 20-thread
- absorbent
- specified formats: Dimensions PRIOR to washing in production

#### Usage

For absorbing and staunching blood and bodily fluids during surgical procedures. The abdominal swabs come in various sizes, layer counts and colours to suit a wide range of purposes.

Size

### Non-sterile, pack of 5 pcs, bundled with wrapper



		REF	Size (width x length)	Loop	Pack contents	Shipping unit	PZN
4-ply	White	1613	20 x 30 cm (approx.)	Yes	4 x 5 pcs	36 packs	13506907
		1547	40 x 40 cm (approx.)	Yes	4 x 5 pcs	18 packs	13506913
		1616	45 x 45 cm (approx.)	Yes	4 x 5 pcs	12 packs	13506936
		1648	45 x 45 cm (approx.)	No	4 x 5 pcs	12 packs	13506942
		1670	50 x 60 cm (approx.)	Yes	4 x 5 pcs	8 packs	13506959
	Unwashed	1602	50 x 60 cm (approx.)	No	4 x 5 pcs	10 packs	13506965
		1465	8 x 90 cm (approx.)	Yes	4 x 5 pcs	30 packs	13506994
	Green	1614	20 x 30 cm (approx.)	Yes	4 x 5 pcs	36 packs	13507002
		1617	45 x 45 cm (approx.)	Yes	4 x 5 pcs	12 packs	13507019
		1607	45 x 45 cm (approx.)	No	4 x 5 pcs	12 packs	13507025
		1696	8 x 90 cm (approx.)	Yes	4 x 5 pcs	30 packs	13507031
	Blue	4593	45 x 45 cm (approx.)	Yes	4 x 5 pcs	12 packs	13507048
6-ply	White	1612	20 x 30 cm (approx.)	Yes	4 x 5 pcs	30 packs	13507054
		1502	40 x 40 cm (approx.)	Yes	4 x 5 pcs	15 packs	13507083
		1618	45 x 45 cm (approx.)	Yes	4 x 5 pcs	10 packs	13507108
	Unwashed	1520	45 x 45 cm (approx.)	Yes	4 x 5 pcs	10 packs	13507120
	-	1659	45 x 45 cm (approx.)	No	4 x 5 pcs	10 packs	13507137
		1671	50 x 60 cm (approx.)	Yes	4 x 5 pcs	6 packs	13507143
		1449	8 x 90 cm (approx.)	Yes	4 x 5 pcs	27 packs	13507166
	Green	1615	20 x 30 cm (approx.)	Yes	4 x 5 pcs	30 packs	13507172
	24-thread	1496	20 x 30 cm (approx.)	Yes	4 x 5 pcs	10 packs	13507189
		1608	45 x 45 cm (approx.)	Yes	4 x 5 pcs	10 packs	13507203
		1651	45 x 45 cm (approx.)	No	4 x 5 pcs	10 packs	13507226
	24-thread	1583	50 x 60 cm (approx.)	Yes	2 x 5 pcs	10 packs	13507232
		1603	50 x 60 cm (approx.)	No	4 x 5 pcs	6 packs	13507249
		4058	8 x 90 cm (approx.)	Yes	4 x 5 pcs	27 packs	13507255
	Blue	4592	45 x 45 cm (approx.)	Yes	4 x 5 pcs	10 packs	13507261
		4594	45 x 45 cm (approx.)	No	4 x 5 pcs	10 packs	13507278
8-ply	White	1506	45 x 45 cm (approx.)	Yes	2 x 5 pcs	10 packs	13507284



## Double sterile, with self-adhesive counting cards

_				REF	Size (width x length)	Loop	Pack contents	Shipping unit	PZN
4-ply	Pack of 2 pcs	White	Double-peel	1664	45 x 45 cm (approx.)	Yes	15 x 2 pcs	4 packs	13166417
			Double-peel	4515	50 x 60 cm (approx.)	Yes	24 x 2 pcs	2 packs	13166423
			Double-peel	1682	8 x 90 cm (approx.)	Yes	30 x 2 pcs	4 packs	13166446
		Green	Double-peel	1665	45 x 45 cm (approx.)	Yes	15 x 2 pcs	4 packs	13166452
			Double-peel	4504	50 x 60 cm (approx.)	Yes	24 x 2 pcs	2 packs	13166469
			Double-peel	1684	8 x 90 cm (approx.)	Yes	30 x 2 pcs	4 packs	13166475
	Pack of 5 pcs	White	Double-peel	1571	20 x 30 cm (approx.)	Yes	15 x 5 pcs	4 packs	13614399
			Double-peel	1593	45 x 45 cm (approx.)	Yes	15 x 5 pcs	2 packs	13166498
	Unwashed, 24-thread	Green	Double-peel	1597	30 x 30 cm (approx.)	No	18 x 5 pcs	2 packs	13166297
			Double-peel	1572	20 x 30 cm (approx.)	Yes	15 x 5 pcs	4 packs	13166280
			Double-peel	1594	45 x 45 cm (approx.)	Yes	15 x 5 pcs	2 packs	13166334
6-ply	Pack of 2 pcs	White	Double-peel	1666	45 x 45 cm (approx.)	Yes	15 x 2 pcs	4 packs	13166357
			Double-peel	1677	45 x 45 cm (approx.)	No	15 x 2 pcs	4 packs	13166363
			Double-peel	1685	8 x 90 cm (approx.)	Yes	30 x 2 pcs	4 packs	14218665
	Pack of 2 pcs	Green	Double-peel	1667	45 x 45 cm (approx.)	Yes	15 x 2 pcs	4 packs	13166386
			Double-peel	4501	50 x 60 cm (approx.)	Yes	20 x 2 pcs	2 packs	13166185
			Double-peel	1680	8 x 90 cm (approx.)	Yes	30 x 2 pcs	4 packs	13166191
	Pack of 5 pcs	White	Double-peel	1595	45 x 45 cm (approx.)	Yes	12 x 5 pcs	2 packs	13166222
			Double-peel	1591	45 x 45 cm (approx.)	No	12 x 5 pcs	2 packs	13166239
			Double-peel	1669	8 x 90 cm (approx.)	Yes	24 x 5 pcs	2 packs	13166245
	Pack of 5 pcs	Green	Double-peel	1574	20 x 30 cm (approx.)	Yes	9 x 5 pcs	4 packs	13166251
			Double-peel	1587	40 x 40 cm (approx.)	Yes	12 x 5 pcs	2 packs	13166268
			Double-peel	1596	45 x 45 cm (approx.)	Yes	12 x 5 pcs	2 packs	13166400

## BeeSana<sup>®</sup> Gauze swab X-ray

## Our classic gauze swab with X-ray contrast thread

#### Features

- gauze swabs (EN 14079 standard)
- woven-in X-ray contrast threads (BaSO<sub>4</sub>)
- selvedge cut edges
- folded, 17-thread
- absorbent
- non-sterile in sterilisable paper packaging
- double sterile packaged with self-adhesive counting cards for secure documentation during operations





#### Usage

For (invasive) surgical interventions in the operating theatre. For staunching bleeding and wound care.

		REF	Size (width x length)	Pack contents	Shipping unit	PZN
8-ply		1013	10 x 10 cm (approx.)	100 pcs	10 packs	12348599
		1037	10 x 20 cm (approx.)	100 pcs	5 packs	12348719
12-ply		1003	5 x 5 cm (approx.)	100 pcs	10 packs	12348582
		1009	7.5 x 7.5 cm (approx.)	100 pcs	10 packs	06613703
		1015	10 x 10 cm (approx.)	100 pcs	10 packs	06613726
_		1039	10 x 13 cm (approx.)	100 pcs	10 packs	12348731
		1024	10 x 20 cm (approx.)	100 pcs	5 packs	06613732
	10 pcs, pre-counted	2715	10 x 10 cm (approx.)	200 x 10 pcs	1 box	_
	For set packers	2724	10 x 20 cm (approx.)	100 x 10 pcs	1 box	_
16-ply		1005	5 x 5 cm (approx.)	100 pcs	10 packs	06613749
		1011	7.5 x 7.5 cm (approx.)	100 pcs	10 packs	06613755
		1017	10 x 10 cm (approx.)	100 pcs	10 packs	06613761
		1054	10 x 13 cm (approx.)	100 pcs	10 packs	12348665
		1034	10 x 20 cm (approx.)	100 pcs	5 packs	12348702
	10 pcs, pre-counted	2717	10 x 10 cm (approx.)	200 x 10 pcs	1 box	-
	For set packers	2727	20 x 20 cm (approx.)	60 x 10 pcs	1 box	_
24-ply		1090	10 x 20 cm (approx.)	50 pcs	5 packs	12348642
32-ply		1019	10 x 10 cm (approx.)	50 pcs	10 packs	12348607
	20-thread	1021	10 x 10 cm (approx.)	50 pcs	10 packs	12348748



## Double sterile, with self-adhesive counting cards

		REF	Size (width x length)	Pack contents	Shipping unit	PZN
12-ply	Pack of 5 pcs	1149	10 x 10 cm (approx.)	25 x 5 pcs	8 packs	12348754
	_	1150	10 x 20 cm (approx.)	30 x 5 pcs	4 packs	12348760
	Pack of 10 pcs	4392	5 x 5 cm (approx.)	30 x 10 pcs	8 packs	12348932
	_	1174	7.5 x 7.5 cm (approx.)	40 x 10 pcs	4 packs	12348843
	-	1049	10 x 10 cm (approx.)	30 x 10 pcs	4 packs	11703472
	-	1168	10 x 20 cm (approx.)	15 x 10 pcs	4 packs	12348777
	Pack of 20 pcs	1176	7.5 x 7.5 cm (approx.)	20 x 20 pcs	4 packs	12348866
	-	4314	10 x 10 cm (approx.)	15 x 20 pcs	4 packs	12348820
	_	1041	10 x 20 cm (approx.)	20 x 20 pcs	2 packs	12348659
16-ply	Pack of 10 pcs	4368	5 x 5 cm (approx.)	30 x 10 pcs	8 packs	12348990
	-	1071	7.5 x 7.5 cm (approx.)	30 x 10 pcs	4 packs	12348613
	-	1072	10 x 10 cm (approx.)	25 x 10 pcs	4 packs	12348636
	-	1142	10 x 20 cm (approx.)	15 x 10 pcs	4 packs	12348808
-	Pack of 20 pcs	4351	10 x 10 cm (approx.)	12 x 20 pcs	4 packs	12348984
32-ply	Pack of 10 pcs	4318	10 x 10 cm (approx.)	15 x 10 pcs	4 packs	12348837

## **BeeSana**<sup>®</sup> Gauze balls X-ray

#### Features

- gauze ball, medical gauze (EN 14079 standard)
- woven-in X-ray contrast thread (BaSO<sub>4</sub>)
- innovative wrapping system from just one piece of pre-cut gauze
- inner-facing cut edges
- 20-thread
- absorbent
- for documentation during the operation with counting cards (for sterile gauze swabs with X-ray contrast only)

## Non-sterile, in polybag

	REF	Size (width x length)	Shape	Pack contents	Shipping unit	PZN
Plum-sized	1154	20 x 20 cm (approx.)	Round	250 pcs	16 packs	00822297
Egg-sized	1155	25 x 26 cm (approx.)	Round	250 pcs	10 packs	00822311
Egg-sized	1110	25 x 26 cm (approx.)	Thimble	250 pcs	10 packs	00822305
X-sized	1349	30 x 30 cm (approx.)	Round	100 pcs	20 packs	13508243
XX-sized	1114	30 x 40 cm (approx.)	Round	100 pcs	15 packs	01363868
XXX-sized	1115	40 x 40 cm (approx.)	Round	100 pcs	10 packs	13508266
Fist-sized	1117	40 x 50 cm (approx.)	Round	100 pcs	10 packs	00822328
Abdominal swab	1119	60 x 80 cm (approx.)	Round	50 pcs	8 packs	13508295



## Double sterile, with self-adhesive counting cards

Pack of 5 pcs	Egg-sized	4324	25 x 26 cm (approx.)	Round	25 x 5 pcs	8 packs	13165145
_	XX-sized	4326	30 x 40 cm (approx.)	Round	30 x 5 pcs	4 packs	13165151
	Fist-sized	4325	40 x 50 cm (approx.)	Round	20 x 5 pcs	4 packs	13165168
Pack of 10 pcs	Walnut-sized	1180	15 x 16 cm (approx.)	Round	30 x 10 pcs	8 packs	13165174
	Plum-sized	1181	20 x 20 cm (approx.)	Round	18 x 10 pcs	8 packs	13165180
	Egg-sized	1075	25 x 26 cm (approx.)	Round	30 x 10 pcs	4 packs	11703503
_	Egg-sized	1186	25 x 26 cm (approx.)	Thimble	30 x 10 pcs	4 packs	13165197
	X-sized	1182	30 x 30 cm (approx.)	Round	20 x 10 pcs	4 packs	13165211
	XX-sized	4319	30 x 40 cm (approx.)	Round	18 x 10 pcs	4 packs	13165292
-	Fist-sized	4333	40 x 50 cm (approx.)	Round	12 x 10 pcs	4 packs	13165228
Pack of 20 pcs	Plum-sized	4329	20 x 20 cm (approx.)	Round	25 x 20 pcs	4 packs	13165234
_	Egg-sized	4337	25 x 26 cm (approx.)	Round	18 x 20 pcs	4 packs	13165240
-	Egg-sized	4012	25 x 26 cm (approx.)	Thimble	18 x 20 pcs	4 packs	13165257
-	XX-sized	4320	30 x 40 cm (approx.)	Round	20 x 20 pcs	2 packs	13165263



## **BeeSana**<sup>®</sup> Preparation balls X-ray

## Features

- gauze swab, 100% cotton (EN 14079 standard)
- woven-in X-ray contrast thread (BaSO<sub>4</sub>)
- 24-thread
- innovative wrapping system from just one piece of pre-cut gauze
- highly absorbent

## Usage

For free preparation of tissue layers and anatomical structures.

## Non-sterile

	REF	Size (width x length)	Pack contents	Shipping unit	PZN
Packed in polybag	1169	4 x 4 cm (approx.)	100 pcs	50 packs	13583440
	1170	6 x 6 cm (approx.)	100 pcs	40 packs	01363839
	1171	8 x 8 cm (approx.)	100 pcs	30 packs	01363822
	1173	10 x 10 cm (approx.)	100 pcs	20 packs	12386275





## **Double sterile,** with self-adhesive counting cards

Pack of 10 pcs	4089	6 x 6 cm (approx.)	60 x 10 pcs	2 packs	13165731
In paper bag	4311	8 x 8 cm (approx.)	60 x 10 pcs	2 packs	13165754
	4339	12 x 12 cm (approx.)	50 x 10 pcs	2 packs	13165777
Pack of 10 pcs	4343	4 x 4 cm (approx.)	20 x 10 pcs	9 packs	13585373
In counting box	4340	6 x 6 cm (approx.)	20 x 10 pcs	9 packs	12731105
	4341	8 x 8 cm (approx.)	20 x 10 pcs	9 packs	12731111
	4342	12 x 12 cm (approx.)	10 x 10 pcs	4 packs	12731128

## **BeeSana**<sup>®</sup> Pointed swabs X-ray

## Our pointed swabs made from pure cotton

## Features

- gauze swab, 100% cotton (EN 14079 standard)
- 20-thread
- tapered, formed from just one piece of gauze
- woven-in X-ray contrast thread (BaSO<sub>4</sub>)



## Usage

For blood and exudate absorption, especially in ENT applications. For staunching bleeding or as a tamponade.

## Non-sterile

REF	Size (width x length, unrolled)	Length	Pack contents	Shipping unit	PZN
1238	Large, 15 x 16 cm (approx.)	13.5 cm (approx.)	500 pcs	20 packs	13509403
			BooSar	ditrade	

## Double sterile, with self-adhesive counting cards

	REF	Size (width x length, unrolled)	Length	Pack contents	Shipping unit	PZN
Pack of 10 pcs	4004	Small, 6 x 8 cm (approx.)	8 cm (approx.)	60 x 10 pcs	8 packs	13165783
	4161	Large, 15 x 16 cm (approx.)	13.5 cm (approx.)	40 x 10 pcs	8 packs	13165808

## BeeSana<sup>®</sup> Gauze ball X-ray

## Our gauze balls made from pure cotton

## Features

- gauze swab, medical gauze (EN 14079 standard)
- with green retrieval thread
- innovative wrapping system from just one piece of pre-cut gauze
- inner-facing cut edges
- 20-thread
- absorbent

## Usage

For blood and exudate absorption. The swab can be safely retrieved using the thread.

## Non-sterile

	REF	Size (width x length)	Shape	Ball	Pack contents	Shipping unit	PZN
Plum-sized	4043	20 x 20 cm (approx.)	Round	Double	200 pcs	5 packs	12390911

## **ABE**<sup>®</sup> ribbon gauze X-ray

## Our tamponade made from pure cotton

## Features

- medical gauze (EN 14079 standard)
- 100% cotton
- woven edges
- 20-thread
- 4-ply
- rolled
- woven-in X-ray contrast thread (BaSO<sub>4</sub>)

## Usage

For tamponading wound cavities or fistula tracts. Absorption of bodily fluids in hard-to-reach wound and bodily cavities.

#### Non-sterile Size PZN REF **Pack contents** Shipping unit (width x length) 20-thread 1446 6 cm x 2 m (approx.) 40 pcs 12386298 10 packs Double sterile, pack of 1 pc. with counting card Size REF Pack contents Shipping unit PZN (width x length) 20-thread 02214231 4312 6 cm x 2 m (approx.) 25 x 1 pc. 4 packs







# **BeeSana**<sup>®</sup> Sets – extremely cost-effective yet offering maximum customisation

Our set solutions offer a wide range of combinations of components in a single pack enabling our clients to use our products safely and efficiently. Material use, preparation time, waste volume and logistics costs are all reduced with our sets that come in a wide range of individual components ensuring high quality patient care. All components are ready for use in the amount precisely allocated to the patient and procedure concerned.

Of course, all our sets for the operating theatre are fitted with double, self-adhesive counting cards to ensure reliable recording of the details in the patient files. The use of set solutions leads to an optimisation of process and great costeffectiveness for the client. Set solutions are available both as prefabricated standard sets and as bespoke sets put together according to the client's wishes.

#### **Bespoke set solutions**

The contents of each set are put together according to the specific requirements of our clients. Our trained and competent sales staff work with our clients to understand the requirements for special procedures and come up with tailored solutions to meet the highest demands. Clients can draw from our entire product range, comprising bandages, coverings, protective clothing and other materials.

### Standardised set solutions

Given our years of experience in creating bespoke set solutions for our clients, we have been able to develop a range of standard sets for various procedures both in the operating theatre and in the ward. Combinations that have been specially adapted to the requirements of various procedures make it easier for our clients to prepare and carry out standard procedures.

## **BeeSana**<sup>®</sup> Ward catheter set

## Features

The ward catheter set is used for hygienic intermittent bladder catheterisation and is sterile, to enable immediate use on the ward.



## Usage

For bladder catheterisation.

## Sterile

REF	Contents	Pack contents	Shipping unit	PZN
4265	1x rectangular dish, 650 ml 1x swab dish, 125 ml 2x gauze dressing pads, 7.5 x 7.5 cm, 8-ply 6x gauze swabs, 20 x 20 cm, plum size 1x pair of plastic tweezers 1x slit/hole cloth, 60 x 60 cm Wrapped in: 1 cover cloth, 60 x 60 cm	21 sets	2 packs	04136073
4267	1x rectangular dish, 650 ml 1x swab dish, 125 ml 2x gauze dressing pads, 7.5 x 7.5 cm, 8-ply 6x gauze swabs, 20 x 20 cm, plum size 2x pairs of examination gloves, size L, latex 1x pair of plastic tweezers 1x slit/hole cloth, 60 x 60 cm Wrapped in: 1 cover cloth, 60 x 60 cm	21 sets	2 packs	04136067

## **BeeSana**<sup>®</sup> PEG bandage set

## Features

The PEG bandage set is used to hygienically position a bandage during percutaneous endoscopic gastrostomy and is sterile for immediate use on the ward.



## Usage

For external wound care.

## Sterile

REF	Contents	Pack contents	Shipping unit	PZN
4653	2x non-woven slit dressing pads, 7.5 x 7.5 cm, 4-ply, 30 g 2x non-woven dressing pads, 7.5 x 7.5 cm, 6-ply, 30 g 1 x ABE® non-woven fixing material 12 x 15 cm	10 sets	30 packs	14399023
4577	2x slit non-woven dressing pads, 7.5 x 7.5 cm, 30 g, 4-ply 2x non-woven dressing pads, 7.5 x 7.5 cm, 30 g, 6-ply 1 x ABE <sup>®</sup> non-woven fixing material 12 x 15 cm	40 sets	8 packs	01333962
4589	2x slit non-woven dressing pads, 7.5 x 7.5 cm, 30 g, 4-ply 2x non-woven dressing pads, 7.5 x 7.5 cm, 30 g, 6-ply 1 x ABE® non-woven fixing material 12 x 15 cm	30 sets	8 packs	04152190
4836	2x gauze dressing pads, 7.5 x 7.5 cm, 8-ply 1x drainage tube dressing pad (aluminium/non-woven) 6 x 7 cm, with slit/hole 1 x ABE <sup>®</sup> non-woven fixing material 12 x 15 cm	40 sets	8 packs	15995170

1 x fixation holder, size 2, small, with hook-and-pile fastener

## **BeeSana**<sup>®</sup> Bandage change set

## Features

The bandage change set is used for hygienically changing bandages and is sterile for immediate use on the ward.



## Usage

For external wound care.

REF	Contents	Pack contents	Shipping unit	PZN
4820	2x gauze dressing pads, 10 x 10 cm, 8-ply 4x gauze swabs, 20 x 20 cm, plum size	25 sets	4 packs	00549660
5413	2x gauze dressing pads, 5 x 5 cm, 16-ply 4x gauze dressing pads, 10 x 10 cm, 12-ply	20 sets	4 packs	_

## **BeeSana**<sup>®</sup> Suture removal set

Sterile

### Features

The suture removal set is used for the hygienic removal of sutures and is sterile, so as to enable its immediate use on the ward.



#### Usage

For external wound care.

## Sterile

REF	Contents	Pack contents	Shipping unit	PZN
5450	4x gauze dressing pads, 10 x 10 cm, 12-ply	15 sets	4 packs	_
	4x gauze swabs, 25 x 26 cm, egg size			
	1x pair of suture removal scissors, 11.5 cm			
	2x disposable tweezers, 130 mm, blue			
	1x protective covering, 45 x 75 cm, 2-ply, green			

## **BeeSana**<sup>®</sup> Ward protection set

## Features

The ward protection set is designed to prevent patient-to-patient contamination and the transfer of an infection to medical and/or nursing staff.





## Non-sterile

	REF	Contents	Pack contents	Shipping unit	PZN
t germs.	4074	1x PP gown	80 sets	1 pack	
		1x nurse's cap			
		1x surgical mask			_
		2x CPE shoe covers			
		2x examination gloves, latex, size M			

## Usage

To protect against germs.



# Our incontinence products are designed to ensure quality of life for patients

A variety of conservative and surgical treatments are now available to cure or at least alleviate urinary and faecal incontinence. In cases where patients can no longer be treated, modern incontinence products provide a solution. When developing and producing Meditrade incontinence products, we strive to achieve well-designed, ergonomic solutions that provide discretion for the user.

Our comprehensive range of products enables costeffective use with a wide range of sizes, absorbency levels and qualities, always following the principle of achieving targeted and acceptable patient care.

# **Superslip**<sup>®</sup>

## Anatomically shaped incontinence slip with flexible waist

#### Features

Incontinence briefs with high liquid absorption, elastic back closure and double hook-and-pile fastener. Hydrophobic side seams and a dampness indicator to show when it's time to change.

#### Usage

For urinary and faecal incontinence, and for people requiring all-round care with reduced mobility.



## Superslip<sup>®</sup> Classic

REF	Size	Colour code	Hip measurement	Absorbency (ISO 11948-1 standard)	Pack contents	Shipping unit	PZN	Assistive technology no.
2431	Medium	White/blue	80–115 cm	1,850 ml (approx.)	20 pcs	4 packs	07757522	15.25.03.1118
2432	Large	White/green	115–150 cm	2,100 ml (approx.)	20 pcs	4 packs	07757551	15.25.03.2112

## Superslip<sup>®</sup> Comfort

2406	Small	White	50–80 cm	1,800 ml (approx.)	15 pcs	6 packs	07757485	15.25.03.0066
2407	Medium	White/blue	80–115 cm	2,450 ml (approx.)	15 pcs	8 packs	07757580	15.25.03.1119
2408	Large	White/green	115–150 cm	2,600 ml (approx.)	15 pcs	6 packs	07757605	15.25.03.2113
2414	X-Large	White/violet	>150 cm	2,700 ml (approx.)	15 pcs	4 packs	07757634	15.25.03.2312

## Superslip<sup>®</sup> Extra

2436 Large White/green 115–150 cm 3,100 ml (approx.)	15 pcs	4 packs	09747019	15.25.03.2314
------------------------------------------------------	--------	---------	----------	---------------

## Superslip<sup>®</sup> Super

2403 Medium	White/blue	80–115 cm	3,000 ml (approx.)	15 pcs	4 packs	09754143	15.25.03.1278
2404 Large	White/green	115–150 cm	3,250 ml (approx.)	15 pcs	4 packs	00351917	15.25.03.2313

## **Superslip**<sup>®</sup> Air Comfort Anatomically shaped incontinence slip with flexible waist

#### Features

Incontinence briefs with high liquid absorption, elastic back closure and double hook-and-pile fastener. Hydrophobic side seams and a dampness indicator to show when it's time to change.

#### Usage

For urinary and faecal incontinence, and for people requiring all-round care with reduced mobility.



RE	F	Size	Colour code	Hip measurement	Absorbency (ISO 11948-1 standard)	Pack contents	Shipping unit	PZN	Assistive technology no.
243	37	Medium	White/blue	80–115 cm	1,750 ml (approx.)	15 pcs	8 packs	04917043	15.25.03.1276
243	38	Large	White/green	100–150 cm	2,000 ml (approx.)	15 pcs	6 packs	04918189	15.25.03.2311

# BeeSana® Medi-Pants

## Discreet disposable briefs

#### Features

Incontinence briefs for those with a weak bladder, easy to put on and take off like normal underwear. Provides unrestricted freedom of movement and a normal level of independence. For women and men, disposable.

#### Usage

For patients suffering from a weak bladder. For mobile or dementia patients without any additional restrictions to daily life. Suitable for use with toilet training. To be worn like normal underwear.

REF	Size	Hip measurement	Pack contents	Shipping unit	PZN
2411	Medium	70–110 cm (approx.)	14 pcs	6 packs	11141726
2412	Large	100–140 cm (approx.)	14 pcs	6 packs	11141732
2413	X-Large	120–165 cm (approx.)	14 pcs	6 packs	11141749



## Liberty Light Anatomically formed incontinence pad

## Features

Anatomical fit with hydrophobic seams in the crotch area. The secure adhesive strip is used to fix this pad into underwear.

## Usage

For patients with a weak bladder, for protection after vaginal procedures or as a pad for post-natal women.



REF	Size	Colour code	Dimensions (approx.)	Absorbency (ISO 11948-1 standard)	Pack contents	Shipping unit	PZN	Assistive technology no.
2052	Normal	Rose	28 x 10.5 cm	min. 350 ml	22 pcs	10 packs	07305197	15.25.01.5175
2053	Extra	Blue	31.5 x 12 cm	min. 400 ml	18 pcs	10 packs	07305180	15.25.01.5177
2054	Super	Green	43 x 22 cm	min. 850 ml	30 pcs	5 packs	07305174	15.25.01.5176

# **Superform**<sup>®</sup>

## Anatomically formed incontinence pad with double cushion

## Features

Incontinence pad with anatomical fit. Hydrophobic side seams as a leak barrier and a dampness indicator to show when it's time to change.

## Usage

Suitable for both bed-bound and mobile incontinence patients.

**Note on product combination:**To fix the pad in place, we recommend using our washable fixation briefs, PANTY or LIGHT.



REF	Size	Colour code	Dimensions	Absorbency (ISO 11948-1 standard)	Pack contents	Shipping unit	PZN	Assistive technology no.
2047	Classic	Blue	63 x 30 cm	1,200 ml (approx.)	30 pcs	5 packs	07520263	15.25.01.0076
2048	Plus	Blue	63 x 30 cm	1,350 ml (approx.)	30 pcs	5 packs	07414302	15.25.01.1096
2049	Extra	Green	70 x 37 cm	1,450 ml (approx.)	30 pcs	4 packs	07414288	15.25.01.1105
2050	Super	Green	70 x 37 cm	1,850 ml (approx.)	20 pcs	4 packs	07414271	15.25.01.1065
2051	Maxi	White	70 x 37 cm	2,450 ml (approx.)	20 pcs	4 packs	07414070	15.25.01.2101

# Superform<sup>®</sup> Air

## Anatomically formed incontinence pad with double cushion

## Features

Incontinence pad with anatomical fit. Hydrophobic side seams as a leak barrier and a dampness indicator to show when it's time to change.

## Usage

Suitable for both bed-bound and mobile incontinence patients.

**Note on product combination:**To fix the pad in place, we recommend using our washable fixation briefs, PANTY or LIGHT.



REF	Size	Colour code	Absorbency (ISO 11948-1 standard)	Pack contents	Shipping unit	PZN	Assistive technology no.
20481	Plus	White/blue	1,250 ml (approx.)	30 pcs	5 packs	08801840	15.25.01.0098
20491	Extra	White/yellow	1,400 ml (approx.)	30 pcs	4 packs	08801857	15.25.01.1125
20501	Super	White/green	1,800 ml (approx.)	20 pcs	4 packs	08801863	15.25.01.2141
20511	Maxi	White/violet	2,200 ml (approx.)	20 pcs	4 packs	08834851	15.25.01.2142

## **BeeSana**<sup>®</sup> Absorbent cushion

## Pad for incontinence and nursing care

## Features

Latex-free absorbent cushion from cellulose flock with non-woven sleeve and outer film. This superabsorbent material protects underwear and clothing from soiling. Simple to apply in combination with Light or Panty fixation briefs.

## Usage

For incontinence/loose stools, for catheter wearers, bed-ridden and mobile patients. For use as a universal absorbent cushion in nursing care.





	REF	Size	Pack	Pack contents	PZN
With adhesive strip	2022	12 x 38 cm	50 pcs	8 packs	02214142
- Without adhesive strip	2024	18 x 40 cm	60 pcs	3 packs	02214165

## BeeSana® Absorbent pad

## Our large pad for incontinence care

## Features

Latex-free absorbent cushion from cellulose flock with non-woven sleeve. Simple to apply in combination with Light or Panty fixation briefs.

#### Usage

For incontinence for bed-ridden people.



REF	Size	Pack contents	Shipping unit	PZN
2119	16 x 55 cm	20 pcs	5 packs	02214219

## **BeeSana**<sup>®</sup> Rectangular pad

## Our absorbent pad for incontinence

#### **Features**

This latex-free absorbent cushion from cellulose flock with a non-woven sleeve is simple to apply in combination with Light or Panty fixation briefs.

#### Usage

For incontinence for bed-ridden people. For post-natal women in the period after giving birth.

Suitable for protection after vaginal procedures.

REF	Size	Pack contents	Shipping unit	PZN
2101	35 x 11 cm	25 pcs	8 packs	12147808



## **BeeSana**<sup>®</sup> Male incontinence pad

## Anatomically formed pad for men

## Features

Latex-, PVC- and DEHP-free. Ergonomic cut to fit the male anatomy. Not noticeable when worn. With adhesive strips for fixing to underwear.

## Usage

For men with urinary incontinence, such as after surgery in the genital area.

REF	Size	Pack contents	Shipping unit	PZN
2020	8.5 x 40 cm	21 pcs	7 packs	05880164

## **BeeSana**<sup>®</sup> Female incontinence pad

## Our hygienic female incontinence pad

#### **Features**

Latex-free absorbent cushion from absorbent cellulose flock with non-woven sleeve.

## Usage

As protection after vaginal procedures or as a pad for post-natal women after childbirth.

REF	Size	Pack contents	Shipping unit	PZN
2118	7 x 23 cm (approx.)	20 pcs	18 packs	01091275





REF 2118

20 St. / pcs.

Meditrade

# Light fixation briefs

For keeping incontinence pads in place

## Features

Latex-free fixation briefs with underwear-like fabric structure. Size indicated by coloured strips. Washable at 60  $^{\circ}$ C and can be tumble-dried on the delicates setting.

## Usage

For keeping incontinence pads and liners in place.





REF	Size	Colour	Hip measurement	Pack contents	Shipping unit	PZN
2221	Medium	Green	80–110 cm (approx.)	50 pcs	6 packs	11513151
2222	Large	White	90–130 cm (approx.)	50 pcs	6 packs	11513211
2223	X-Large	Blue	110–170 cm (approx.)	50 pcs	6 packs	11513228
2224	XX-Large	Orange	120–180 cm (approx.)	50 pcs	6 packs	11513234

## Panty fixation briefs

With legs for keeping incontinence pads in place

## Features

Latex-free fixation briefs, size indicated by coloured strips. Washable at 60  $^{\circ}$ C and can be tumble-dried on the delicates setting.

## Usage

For keeping incontinence pads and liners in place.



REF	Size	Colour	Hip measurement	Pack contents	Shipping unit	PZN
0711500301	Medium	Blue	70–100 cm (approx.)	100 pcs	6 packs	09081879
0711500501	Large	Brown	85—120 cm (approx.)	100 pcs	6 packs	09081885
0711500701	X-Large	Green	105–130 cm (approx.)	100 pcs	6 packs	09081891
0711500901	XX-Large	Green/blue	115–140 cm (approx.)	100 pcs	6 packs	09081916

## Layered undersheet

## An industry-standard undersheet for meeting the patient's needs

## Features

Undersheet with absorbent core made from environmentally friendly, 100% non-bleached recycled cellulose for fluid absorption. The dry non-woven material with outer film makes use easy in any situation.

#### Usage

For live-in nursing, undersheets help to make nursing easier by simplifying day-to-day care work and, last but not least, reducing the 'laundry mountain' by protecting the bed covers. We can supply the sheets as layered or flock undersheets in a range of different qualities.



REF	Size	Number of layers	Shipping unit	PZN
104606M	40 x 60 cm	6-ply	200 pcs	02426826
104608M	40 x 60 cm	8-ply	200 pcs	00559612
104610M	40 x 60 cm	10-ply	200 pcs	00559629
104612M	40 x 60 cm	12-ply	100 pcs	06178087
104620M	40 x 60 cm	20-ply	100 pcs	00559693
106606M	60 x 60 cm	6-ply	150 pcs	00560213
106610M	60 x 60 cm	10-ply	100 pcs	00560236
106906M	60 x 90 cm	6-ply	100 pcs	00560259
106908M	60 x 90 cm	8-ply	100 pcs	00560331
106910M	60 x 90 cm	10-ply	100 pcs	00560609
106912M50	60 x 90 cm	12-ply	50 pcs	00561046
106920M	60 x 90 cm	20-ply	50 pcs	00561052

## BeeSana<sup>®</sup> Comflock Super

## Undersheet made from cellulose flock

## Features

Undersheet with absorbent core made from cellulose flock for fluid absorption. The dry non-woven material with outer film makes use easy in any situation.



## Usage

For live-in nursing, undersheets help to make nursing easier by simplifying day-to-day care work and, last but not least, reducing the 'laundry mountain' by protecting the bed covers. We can supply the sheets as layered or flock undersheets in a range of different qualities.

REF	Size	Fluid absorbency	colour	Pack contents	Shipping unit	PZN
2064	40 x 60 cm	750 ml ±10%	white	25 pcs	6 packs	04136185
2065	60 x 60 cm	900 ml ±10%	white	25 pcs	6 packs	04136104
2066	60 x 90 cm	1,590 ml ±10%	blue	25 pcs	6 packs	04136156

# Nursing pad

## Our pad to protect underwear during breastfeeding

## Features

Anatomically formed nursing pad. The cellulose core absorbs leaks of milk while offering a discreet fit.





## Usage

For absorbing fluids and to protect underwear during post-natal breastfeeding.

## Non-sterile

REF	Size	Pack contents	Shipping unit	PZN
2068	10.5 cm	30 pcs	12 packs	02097937



# Disinfection and hygiene products with a broad spectrum of action

The core elements of our broad range of products reduce the risk of infection, standardise prevention methods and provide lasting skin care. Please note: Our product recommendations are based on scientific testing and are offered according to the best of our knowledge. More detailed recommendations – such as in terms of material compatibility – can only be provided separately.

Our recommendations do not exempt our customers from testing our products independently to ensure they are suitable for the intended purposes and procedures. We cannot, therefore, assume any liability.

# **Ethasept**<sup>®</sup>

## Our alcohol-based hand sanitiser

## Product type

Ethasept<sup>®</sup> is a ready-to-use, alcohol-based, hypoallergenic disinfectant without perfumes or colourings, and is suitable for disinfection applications in human hygiene.

## **Applications**

For hygienic and surgical hand disinfection in all areas of human hygiene.

## Usage

Clean heavily soiled hands before disinfection and let wet hands dry before starting disinfection. Ready-to-use. Keeps hands/forearms moist during application. Observe standard procedures for hygienic (EN 1500) and surgical (EN 12791) hand disinfection.

## **Product properties**

- fast-acting and broad-spectrum, works after just 15 s
- an effective bactericide, yeasticide, tuberculocide and virucide
- effective against MRSA
- dermatologically tested: "Excellent"
- listed by the German Association for Applied Hygiene (VAH)

## Shelf life

3 years (see label). Shelf life after breaking seal: 12 months. Keep containers tightly closed. Warehousing and storage under standard storage conditions (protect from moisture, direct sunlight, heat and frost).

## Composition

Active substance in 100 g: 73.5 g ethanol.

## **Regulatory requirements**

Ethasept<sup>®</sup> is not a medicinal product and is subject to biocide legislation. Use biocides safely. Always read the label and product information before use.

Type 1 biocide product ("biocide products for human hygiene"), biocide reg. no. N-44236.

Activities	15 sec	30 sec	1 min / 60 sec	1.5 min / 90 sec	2 min / 120 sec
Hygienic hand disinfection: (EN 1500)	+				
Surgical hand disinfection: (EN 12791)				+	
Virucide					+
Virucide (EN 14476)					
Adeno (high load)					+
Polio (high load)			+		

REF	Pack contents	Shipping unit	PZN
01001UK	100 ml	42 bottles	-
01004UK	500 ml	18 bottles	-
01005UK	11	10 bottles	-





# Ethasept<sup>®</sup> gel

## Hand disinfection

## Composition

73.5g Ethanol Ethasept Gel is ready-to-use alcohol-based hand sanitiser gel for hygienic hand disinfection.

## Product Type

Ethasept<sup>®</sup> gel is a ready-to-use, alcohol-based, hypoallergenic disinfectant without perfumes or colourings, and is suitable for disinfection applications in human hygiene.

## Application

Ready-to-use. Thoroughly rub hands with at least 2 x 3 ml onto the hands during 60 seconds, undiluted and without adding water. Observe standard methods for hygienic hand disinfection (EN 1500).

## **Product properties**

- Works quickly and comprehensively after just 60 seconds.
- Bactericidal, yeasticidal and limited virucidal activity
- Effective against MRSA
- Dermatologically tested and rated "Excellent"
- Listed by the German Association for Applied Hygiene (VAH)

## Shelf life

3 years (see label). Shelf life after breaking seal: 12 months. Keep containers tightly closed. Warehousing and storage under standard storage conditions (protect from moisture, direct sunlight, heat and frost).

## **Regulatory requirements**

Ethasept® gel is not a medicinal product and is subject to biocide legislation. Use biocides safely. Always read the label and product information before use.

REF	Box contents	Shipping unit	PZN
02001UK	100 ml	42 bottles	-
02002UK	500 ml	18 bottles	-
02003UK	1	10 bottles	-







## **Alcohol-based surface disinfection**

## Medizid<sup>®</sup> Rapid+

## High-performance, ready-to-use rapid disinfectant for alcohol-resistant surfaces and objects

## Product type

Medizid<sup>®</sup> Rapid+ is a ready-to-use, alcohol-based disinfectant for application using a cloth, with broad-spectrum effectiveness for the disinfection of alcohol-resistant surfaces in healthcare facilities. Ideal for rapid disinfection between treatments and procedures, especially when practical situations mean contact times need to be kept short.

#### **Applications**

For the rapid disinfection of alcohol-resistant surfaces in the patient's immediate environment in hospitals, general practice, outpatient units, nursing homes and manufacturing.

#### Usage

Remove visible soiling before performing disinfection. Wet the surfaces completely with the solution and use a hygienic, clean disposable cloth to distribute the disinfectant solution as

required. Remember to observe the contact times as necessary. For professional use only. Not suitable for the disinfection of invasive medical devices. Test material compatibility on an unobtrusive portion if necessary.

#### **Product properties**

- aldehyde-free
- effectiveness:
- effective bactericide (incl. MRSA), yeasticide
- effective fungicide
- effective mycobactericide
- effective virucide
- listed by the German Association for Applied Hygiene (VAH)

## Shelf life

3 years (see label). Shelf life after breaking seal: 28 days. Keep containers tightly closed. Use moist wipes only. Warehousing and storage under standard storage conditions (protect from moisture, direct sunlight, heat and frost).

REF

Pack contents

Shipping unit

PZN

## Composition

Active substances in 100 g: 0.05 g didecyldimethylammonium chloride **Regulatory requirements** 

Class IIa medical device

			ruen contento	sinpping and	. =	
Activities	fully wetted	01002+	250 ml	12 bottles	16015216	
As per VAH guidelines:		00975+	500 ml	18 bottles	16015222	
Bactericide, yeasticide incl. 4-field test:	1 min. (low load)	00976+	11	12 bottles	16015239	
Bactericide, yeasticide without	2 min. (high load)	00977+	5 l	2 canisters	16015245	
mechanical system:		Accessories (see page 149)				
Mycobactericide without	Aycobactericide without 5 min. (high load)		Canister tap			
mechanical system:		Canister key				
As per EN 16615 (4-field test):	2 min. (high load)	Spray noz	Spray nozzle for 1 l bottle			
		1				



# **Medizid**<sup>®</sup> Rapid+ disinfectant wipes



Disinfectant wipes for alcohol-resistant surfaces on medical devices

## Product type

Medizid<sup>®</sup> Rapid+ is a ready-to-use, alcohol-based disinfectant for application using a cloth, with broad-spectrum effectiveness for the disinfection of alcohol-resistant surfaces in healthcare facilities. Ideal for rapid disinfection between treatments and procedures, especially when practical situations mean contact times need to be kept short.

## **Applications**

For the rapid disinfection of alcohol-resistant surfaces in the patient's immediate environment in hospitals, general practice, outpatient units, nursing homes, manufacturing and public organisations.

## Usage

Remove visible soiling before performing disinfection. Then wipe down the surfaces carefully with the impregnated cloth. Ensure complete wetting of the surface and observe the contact times as necessary. For professional use only. Not suitable for the disinfection of invasive medical devices. Test material compatibility on an unobtrusive portion if necessary.

## **Product properties**

- aldehyde-free
- effective bactericide, yeasticide, mycobactericide, fungicide, virucide
- effective against MRSA
- listed by the German Association for Applied Hygiene (VAH)

## Shelf life

3 years (see label). Shelf life after breaking seal: 28 days. Keep containers tightly closed. Use moist wipes only. Warehousing and storage under standard storage conditions (protect from moisture, direct sunlight, heat and frost).



## Composition

Active substances in 100 g: 0.05 g didecyldimethylammonium chloride, 63.5 g ethanol

## **Regulatory requirements**

Class IIa medical device.

Activities	fully wetted	As per EN 16615 (4-field test) - practical test (for wipes):		
As per VAH guidelines:		Bio wipes	1 min. (high load)	
Bactericide, yeasticide incl. 4-field test:	1 min. (low load)	Standard wipe	2 min. (high load)	
Bactericide, yeasticide without mechanical system:	2 min. (high load)	Long-term bio wipe test	1 min. (high load)	
Mycobactericide without mechanical system:	5 min. (high load)	(28, 60 and 90 days)		

REF	Article	Wipe size	Pack contents	Shipping unit	PZN
00941+	Wipes in a PE dispenser tin	15 x 18 cm (approx.)	150 wipes	12 tins	16015268
00942+	Refill tin	15 x 18 cm (approx.)	150 wipes	20 bags	16015251
00946+	Wipes in a container	29 x 30 cm (approx.)	70 wipes	1 container	16015274
00947+	Refill container	29 x 30 cm (approx.)	70 wipes	6 rolls	16015280

## **Medizid**<sup>®</sup> Rapid QF Flowpack

Alcohol-based hygiene wipes for the disinfection of medical device surfaces and objects

## Product type

Medizid<sup>®</sup> Rapid QF Flowpack is impregnated with the very fast-acting fluid from the alcohol-based product Medizid<sup>®</sup> Rapid. The hygiene wipes are therefore a good choice for disinfecting medical device surfaces and objects used in medical applications, and in the patient's immediate environment between treatments and interventions.

#### Applications

For all alcohol-resistant medical devices, surfaces and objects, especially for device surfaces/objects in the patient's immediate environment in hospitals, general practice, outpatient units and nursing units.

#### Usage

Medizid<sup>®</sup> Rapid QF disinfectant wipes (flowpack) are ready-to-use. The surfaces and objects to be disinfected are wiped down with the Medizid<sup>®</sup> Rapid wipe, ensuring that the surface/object is wetted thoroughly. Exercise care with alcohol-susceptible materials. Do not use for the disinfection of invasive devices.

## **Product properties**

- aldehyde-free
- effective bactericide, yeasticide, mycobactericide and fungicide, "Limited virucide PLUS" rating
- effective against MRSA
- listed by the German Association for Applied Hygiene (VAH)

#### Shelf life

3 years (see label). Shelf life after breaking seal: 28 days. Keep containers tightly closed. Use moist wipes only. Warehousing and storage under standard storage conditions (protect from moisture, direct sunlight, heat and frost).

#### Composition

Active substances in 100 g: 50 g ethanol

### **Regulatory requirements**

Medizid<sup>®</sup> Rapid QF Flowpack is not a medicinal product and is subject to biocide legislation. Use biocides safely. Always read the label and product information before use. Type 1 biocide product ("biocide products for human hygiene"), biocide reg. no. N-75574.

REF	Wipe size	Pack contents	Shipping unit	PZN
00976FP	22 x 20 cm (approx.)	80 pcs	12 packs	_





Activities	30 sec	1 min	5 min
Activity as per EN 16615 practical test (for wipes), bactericide, yeasticide incl. 4-field test (high load): Bio wipes:		+	
Activities for immersion solution: Contact times for hygiene applications			
Yeasticide (EN 1650, low load)		+	
Bactericide, yeasticide (EN 13727, EN 13624, high load)		+	
Mycobactericide as per VAH method 9 (high load)			+
Tuberculocide (EN 14368, low load)		+	
Contact times for food applications			
Bactericide (EN 1276, low load)	+		
Yeasticide (EN 1650, low load)	+		
Bactericide/yeasticide (EN 13697, low load)			+
Virucidal activities (EN 14476)			
Limited virucide PLUS		+	
Vaccinia		+	
Adenovirus		+	
Norovirus		+	

## **Alcohol-free surface disinfection**

# Medizid<sup>®</sup> AF+

Concentrate for surface disinfection

## Product type

Medizid<sup>®</sup> AF+ is a liquid, aldehyde- and alcohol-free concentrate for the disinfection of surfaces not resistant to alcohol on medical devices, and especially device surfaces and objects in the patient's immediate environment.

## **Applications**

Ideal for areas with heavy soiling or odours, such as operating theatre tables, bed frames and surfaces in daily use in medical facilities, manufacturing and public organisations. Also suitable for surfaces not resistant to alcohol. Test material compatibility on an unobtrusive portion if necessary.

## Usage

Remove visible soiling before disinfection; wet surfaces completely and keep surfaces moist during the entire contact time. Use dosing aids and follow instructions; example of use: 10 l of a 1% solution = 9.9 l water + 100 ml Medizid<sup>®</sup> AF+.

## **Product properties**

- alcohol-free concentrate
- aldehyde-free
- effective bactericide, yeasticide, tuberculocide, "Limited virucide PLUS" rating
- high material compatibility
- listed by the German Association for Applied Hygiene (VAH)

## Contact time

Activity as per EN standards and VAH guidelines:	1 min	5 min	30 min	60 min
Bactericide, yeasticide (EN 13626, EN 16615; low load)		0.5%		
Bactericide, yeasticide (EN 16615; high load)	0.5%			
Tuberculocide (EN 14348; high load)				5%
Mycobactericide (EN 14348; incl. M. terrae and M. avium; high load)			2%	1%
Virucidal activities (EN 14476)				
Limited virucide PLUS				4%
Vaccinia Ankara(incl. HIV, HBV, HCV, influenza viruses, high load)		3%	1%	
Adeno-/noroviruses (high load)				4%
Rotaviruses (no load)			1%	
SV40 (high load)				2%



## Shelf life

3 years (see label). Shelf life after breaking seal: 12 months. Keep container tightly closed. Warehousing and storage under standard storage conditions (protect from moisture, direct sunlight, heat and frost).

## Composition

Active substances in 100 g: 2 g 2-phenoxyethanol, 7 g alkyl (C12-16) dimethyl benzyl ammonium chloride, 5 g N-(3-aminopropyl)-N-dodecylpropane-1,3-diamine

## **Regulatory requirements**

Class IIa medical device

REF	Pack contents	Shipping unit	PZN	
00972+	2	6 bottles	16015185	
00971+	5	2 canisters	16015191	
Accessorie	S (see page 149)			
Refillable dosing bottle				
Canister pump for dosing				
Canister ta	р			
Canister ke	ey .			

## Medizid<sup>®</sup> Alpha+

Alcohol-free, ready-to-use disinfection for wipe-resistant surfaces on medical devices

## Product type

Medizid<sup>®</sup> Alpha+ is a ready-to-use, aldehyde- and alcohol-free spray and wipe disinfectant solution with broad-spectrum effectiveness for the disinfection of surfaces in healthcare facilities.

## **Applications**

For healthcare facilities, sanitary/bathroom use and nursing homes.

#### Usage

Remove visible soiling before performing disinfection. Wet the surfaces completely with the solution and use a hygienic, clean disposable cloth to distribute the disinfectant solution as required. Keep the surface moist during the entire contact time. Wipe away residues of the disinfectant after the expiry of the contact time, either with a cloth soaked in water (at least drinking water quality) or a hygienically clean cloth. For professional use only. Always use undiluted.

#### **Product properties**

- Comprehensive bactericide and yeasticide
- "Limited virucide PLUS" rating
- Material compatibility for surfaces sensitive to alcohol
- Listed by the German Association for Applied Hygiene (VAH)

#### Shelf life

3 years (see label). Shelf life after breaking seal: 12 months. Keep container tightly closed. Warehousing and storage under standard storage conditions (protect from moisture, direct sunlight, heat and frost).

#### Contact time

Activities	
As per EN/methods from VAH/DGHM	
Bactericide, yeasticide (EN 16615; high load)	1 min
Bactericide, yeasticide (EN 13624, EN 13727, high load)	1 min
Virucidal activities (EN 14476)	
Limited virucide PLUS	60 min
Vaccinia Ankara (incl. HIV, HBV, HCV, high load)	1 min
Noroviruses (high load	30 min
Adenoviruses (high load)	60 min

REF	Pack contents	Shipping unit	PZN barcode
00964+	1	12 bottles	16015156
Accessorie	S (see page 149)		

Spray head





## Composition

Active substances in 100 g: 0.1 g 2-phenoxyethanol, 0.35 g alkyl (C12-16)dimethyl benzyl ammonium chloride, 0.25 g N-(3-aminopropyl)-N-dodecylpropane-1,3-diamine

**Regulatory requirements** 

Class IIa medical device

# **Medizid**<sup>®</sup> Alpha+ XXL disinfectant wipes

## Alcohol-free disinfectant wipes for wipe-resistant surfaces on medical devices

## Product type

Medizid<sup>®</sup> Alpha+ is a ready-to-use, aldehyde- and alcohol-free wipe disinfectant solution with broad-spectrum effectiveness for the disinfection of surfaces in healthcare facilities.

## Applications

For healthcare facilities, sanitary/bathroom use and nursing homes.

## Usage

Remove visible soiling before performing disinfection. Wet the surfaces completely with the solution and use a hygienic, clean disposable cloth to distribute the disinfectant solution as required. Keep the surface moist during the entire contact time. Wipe away residues of the disinfectant after the expiry of the contact time, either with a cloth soaked in water (at least drinking water quality) or a hygienically clean cloth. For professional use only. Always use undiluted.

## **Product properties**

- ready-to-use
- bactericide and yeasticide, "Limited virucide PLUS" rating
- material compatibility for surfaces sensitive to alcohol
- listed by the German Association for Applied Hygiene (VAH)

## Shelf life

3 years (see label). Shelf life after breaking seal: 12 months. Keep container tightly closed. Use moist wipes only. Warehousing and storage under standard storage conditions (protect from moisture, direct sunlight, heat and frost).

## Contact time

Activities	
As per EN/methods from VAH/DGHM	
Bactericide, yeasticide (EN 16615; high load)	1 min
Bactericide, yeasticide (EN 13624, EN 13727, high load)	1 min
Virucidal activities (EN 14476)	
Limited virucide PLUS	60 min
Vaccinia Ankara (incl. HIV, HBV, HCV, high load)	1 min
Noroviruses (high load	30 min
Adenoviruses (high load)	60 min

REF	Article	Wipe size	Pack contents	Shipping unit	PZN
00969+	Wipes in a container	26 x 29 cm	70 wipes	1 container	16015162
00970+	Refill pack	26 x 29 cm	70 wipes	6 rolls	16015179

## Composition

Active substances in 100 g: 0.1 g 2-phenoxyethanol, 0.35 g alkyl (C12-16)dimethyl benzyl ammonium chloride, 0.25 g N-(3-aminopropyl)-N-dodecylpropane-1,3-diamine

## **Regulatory requirements**

Class IIa medical device







# Medizid<sup>®</sup> Alpha+ Flowpack

Surface disinfection



## Product type

Medizid<sup>®</sup> Alpha+ Flowpack are ready-to-use, aldehyde- and alcohol-free hygiene wipes with a broad-based mechanism of action for the disinfection of surfaces not resistant to alcohol, medical devices and objects in healthcare facilities. Ideal when contact times need to be kept short for practical reasons.

## Applications

Suitable for the rapid disinfection of all surfaces and medical devices not resistant to alcohol in the patient's immediate environment in hospitals, general practice, outpatient units and nursing homes, as well as the food industry, manufacturing and public organisations. Test material compatibility on an unobtrusive portion if necessary. Not for use in the final disinfection of semi-critical and invasive medical devices

## Usage

Remove visible soiling before performing disinfection. Then wipe down the surfaces carefully with the impregnated cloth. Ensure wetting is complete and keep moist during the entire contact time.

## **Product properties**

- comprehensive bactericide and yeasticide
- "Limited virucide PLUS" rating
- material compatibility for surfaces sensitive to alcohol
- wipes made from 100% biodegradable cellulose fibres
- listed by the German Association for Applied Hygiene (VAH)

## Shelf life

3 years (see label). Shelf life after breaking seal: 28 days. Keep container tightly closed. Use moist wipes only. Warehousing and storage under standard storage conditions (protect from moisture, direct sunlight, heat and frost).

## Composition

Active substances in 100 g: 0.10 g 2-phenoxyethanol 0.10 g, 0.35 g alkyl (C12-16)dimethyl benzyl ammonium chloride, 0.25 g N-(3-aminopropyl)-N-dodecylpropane-1,3-diamine

## **Regulatory requirements**

Class IIa medical device

REF	Wipe size	Pack contents	Shipping unit	PZN
00964+FP	22 x 20 cm (approx.)	80 pcs	12 packs	17198644



As per VAH guidelines	
Bactericide, yeasticide (EN 16615; high load)	1 min
Bactericide, yeasticide (EN 13624, EN 13727) (high load)	1 min
Virucidal activities (EN 14476)	
Limited virucide PLUS	60 min
Vaccinia Ankara (incl. HIV, HBV, HCV, high load)	1 min
Noroviruses (high load)	30 min
Adenoviruses (high load)	60 min

## Instrusol<sup>®</sup> AF +

## Instrument disinfection for manual reprocessing

## Product type

Instrusol<sup>®</sup> AF + is a highly effective instrument disinfectant for the reprocessing of industry-standard medical and surgical instruments, as well as thermally sensitive anaesthetic accessories and rigid endoscopes, either by an immersion method or in an ultrasonic bath.

#### Applications

For all applications where medical instruments are reprocessed by the use of immersion procedures and in ultrasonic cleaning systems.

#### Usage

Always follow the usage instructions and reprocessing guidelines offered by instrument manufacturers; note in particular information about the preparatory steps to take, cleaning and disinfection, rinsing and drying, and storage.

## **Product properties**

- high-performance concentrate
- aldehyde- and phenol-free
- bactericide, yeasticide, tuberculocide
- "Limited virucide" rating
- listed by the German Association for Applied Hygiene (VAH)

#### Shelf life

3 years (see label). Shelf life after breaking seal: 12 days. Keep container tightly closed. Warehousing and storage under standard storage conditions (protect from moisture, direct sunlight, heat and frost).

#### Composition

Active substances in 100 g: 5 g alkyl dimethyl ethyl ammonium ethyl sulphate, 1 g polyhexamethylene biguanide hydrochloride, 1.8 g cocospropylene guanidinium diacetate

## **Regulatory requirements**

Class IIb medical device

## Contact time

Activities	5 min	15 min	30 min	60 min
Bactericide, yeasticide (as per VAH) (high load)		4.0 %	3.0 %	2.0 %
Tuberculocide (as per DGHM 2001) (high load)		4.0 %	3.0 %	2.0 %
Virucidal activity:				
Adenovirus (EN 14476) (Iow load)				3.0 %
Murine norovirus (EN 14476) (low load)				4.0 %
Rotavirus (DVV/RKI) (high load)	2.0 %			
SV 40 (DVV/RKI)			3.0 %	
BVDV/Vaccinia (DVV/RKI) (no load)	1.0 %			

REF	Pack contents	Shipping unit	PZN
00979+	500 ml	12 bottles	-
00973+	2	6 bottles	11169699
00974+	5	2 canisters	11169707
Accessorie	S (see page 149)		
Canister pu	Imp for dosing		
Canister ta	р		
Constanting			

Canister key



#### Reprocessing instructions Immersion method

Immerse instruments after use. All surfaces and cavities must be covered and/or filled.

#### Ultrasonic bath

The prescribed frequency range and the prescribed disinfection duration must be maintained. Water temperatures above 45 °C are to be avoided.

#### Product changeover

Instrusol<sup>®</sup> AF + is not compatible with disinfectant products containing aldehyde. On first use and when changing products, instruments and accessories, baths, reprocessing equipment and any application aids used must be thoroughly cleaned and rinsed off.

Instrusol<sup>®</sup> AF + ready-to-use solution is protein-resistant and can be used more than once if required. The ready-to-use solution is usable for up to 7 days if instruments immersed are properly pre-cleaned. If clouding or visible contamination occurs, we recommend renewing the ready-to-use solution earlier.

#### Notes

As a ready-to-use solution, Instrusol<sup>®</sup> AF + is not a hazardous formulation according to the German Ordinance on Hazardous Substances. The ready-to-use solution is not caustic. However, we do recommend taking the normal safety precautions for handling chemicals and wearing the usual kinds of protective clothing.



# Gentle Med<sup>®</sup> wash lotion

#### pH skin neutral and gentle washing lotion

#### Features

Dermatologically tested and rated a "Excellent" washing lotion for daily, gentle all-over body use, with a fresh scent.Prevents the skin drying out. Paraben-free. Alkali-free, free of silicone and paraffin oil.Can be used sparingly thanks to the optional metering pump.

#### Usage

Hand washing, body washing, showering and bathing.

REF	Pack contents	Shipping unit	PZN
00980D	500 ml	12 bottles	11170314
00981D	11	10 bottles	11170320
00982D	5	2 canisters	11170337



**Gentle Med** 

# **Gentle Med**<sup>®</sup> skin care lotion (O/W)

The skin care lotion to protect and care for stressed skin, with almond and jojoba oil

#### Features

Skin care lotion for sensitive and dry skin, pH skin neutral, enriched with urea. Creates a non-greasy protective film and enriches stressed skin with valuable natural oils (almond and jojoba oil) plus vitamin E. The ingredients urea and hyaluronic acid care for the skin and provide moisture.

Dermatologically tested as "Excellent". PEG-free. Silicone- and paraffin oil-free.

ALL N	100	
9		
~		
	at a l	

REF	Pack contents	Shipping unit	PZN
01011	100 ml	30 tubes	11170308
00990D	500 ml	12 bottles	11170231





## Gentle Med<sup>®</sup> cream lotion (W/O)

Skin protection balsam with panthenol and vitamin E

#### Features

Cream lotion with an advanced skincare formula to protect and promote the regeneration of extremely stressed skin. The Gentle Med<sup>®</sup> cream lotion in this product has been dermatologically tested and rated "Excellent": suitable for all-over body care, it is pH skin neutral and enriched with panthenol and vitamin E. Free from silicone oil, PEG and parabens, with marigold extract.

REF	Pack contents	Shipping unit	PZN
01013	100 ml	30 tubes	11169794
00988D	500 ml	12 bottles	11170053



# Gentle Med<sup>®</sup> active gel

## Vitalising and cooling active gel with essential oils

#### Features

Active rubbing gel to stimulate the circulation and enliven the skin, while being skin-friendly. Has a deodorising and cooling effect. Well-suited for massage to relieve tension. Dermatologically tested as "Excellent". Alkali-free. Paraben-free. PEG-free. Silicone- and paraffin oil-free. With essential oils from medicinal plants: menthol, camphor and rosemary.

REF	Pack contents	Shipping unit	PZN
00991D	500 ml	12 bottles	11169788





# Gentle Med<sup>®</sup> oil bath

#### Oil bath with natural oils

#### Features

Medical oil bath enriched with high-quality soya oil. Promotes the physiological regeneration of the skin. Offers relief for sensitive and dry skin. Suitable for use in showers and baths.



REF	Pack contents	Shipping unit	PZN
00992D	500 ml	12 bottles	11170260

Care

## **Gentle Med**<sup>®</sup> skin cleansing and care foam

Gentle cleaning and care with reliable skin protection

#### Features

Fast and gentle cleaning and care of the skin, also for incontinence, patients with a stoma, and for bed-ridden patients. This skin-cleansing care foam contains moisturising and odour absorbing substances. Paraben-free. PEG-free. Silicone- and paraffin oil-free. Contains no allergenic perfumes.

#### Usage

To lather, hold the can (valve) down at an angle, thus ensuring that the container is completely emptied. Leave to act briefly and then wipe off with a disposable cloth. A protective film continues to care for the skin and provide a deodorising effect. Shake the can vigorously before use.

REF	Pack contents	Shipping unit	PZN
00994D	500 ml	12 bottles	11170248

## **Gentle Med**<sup>®</sup> moist care wipes

#### Our care wipes for normal and sensitive skin

#### Features

Moist care wipes with a witch hazel extract texture, based on mild, coconut oil surfactants. Gentle Med<sup>®</sup> moist care wipes are PEG-, paraben-, paraffin- and silicone-free.

#### Usage

Ready-to-use, high-quality care wipes for the effective, rapid and especially mild care and cleansing of soiled areas of the skin. Disposable to prevent germs spreading any further. For use without water. Please re-seal the pack carefully after removing wipes.

P		
Santte Reinigung und Priege mit zuverba	Minimu Musecture Minimu Musecture 84 Stück	
Startle Winder Control Control Contro	oRIGI excelle	V P. ent

entle Med

mat

REF	Size	Pack contents	Shipping unit	PZN
2356	18 x 20 cm (approx.)	80 pcs	12 packs	12357658





# **Application aids**

### Our application aids for the metering and dosing of disinfectants and care products

#### Metering pumps/attachments

Enable products to be pumped out of the original containers comfortably and hygienically.

_	REF	Article		PZN
	00929	Metering pump	500 ml	_
	00930	Metering pump	1	_
	00931	Canister metering pump	2  , 5   and 10	_

#### (Foam) spray heads

Ensure even wetting, even in difficult-to-access areas, such as for extensive spraying with surface disinfectants. Also suitable for overhead use.

REF	Article	PZN
00934	Foam spray head 1 l	_

#### Canister tap

Ensure products can be removed from cans cleanly and without dripping.

REF	Article	PZN
00996	Canister tap 2 l, 5 l and 10 l	-

#### Canister key

For effortlessly opening cans.

REF	Article
00022	Contractor

00932 Canister key 5 | and 10 |



#### Dosing bottle/graduated flask

For removing the product in a user-friendly way when decanting or for making ready-to-use solutions or concentrates.

REF	Article		PZN
00943	Dosing bottle	11	-
00945	Graduated flask	1,000 ml	-
00927	Graduated flask	2,000 ml	_



# **Disinfectant dispenser**

For the dosed dispensing of alcohol-based hand disinfectants, liquid soaps and lotions

#### Aluminium disinfectant dispenser

Soap/disinfectant dispenser for disposable bottles and refill containers, aluminium housing, short operating lever, 500 ml or 1,000 ml. Suitable for wall mounting. Models with reusable or disposable pump, can be replaced from the front.

	REF Product name		Model
1417549 Aluminium dispenser 500 ml with rea		with reusable pump	
	4402102 Aluminium dispenser 500 ml		with DHP disposable pump
1417589 Aluminium dispenser 1,000 ml		Aluminium dispenser 1,000 ml	with reusable pump
-	4402107	Aluminium dispenser 1,000 ml	with DHP disposable pump



#### Touchless disinfectant dispenser

Touchless soap/disinfectant dispenser for an even greater level of hygiene. With drip tray, made from rugged plastic, 500 ml, with rear wall, lockable. Suitable for wall mounting. Models with reusable or disposable pump, can be replaced from the front.

REF		Product name	Model
	4402237	touchless plastic dispenser 500 ml	with reusable pump
	4402283	touchless plastic dispenser 500 ml	with DHP disposable pump



#### Plastic disinfectant dispenser

Manual dispenser for disinfectants, liquid soaps and low-viscosity lotions, suitable for wall mounting, white plastic housing with white drip tray, stainless steel operating lever. Model with DHP disposable pump or reusable pump, can be replaced from the front.

REF	Product name	Model
4401841	Plastic dispenser 500 ml	with reusable pump
4401974	Plastic dispenser 500 ml	with DHP disposable pump



#### DHP disposable pump

For even greater hygiene and safety! Recyclable plastic disposable pump for alcohol-based hand disinfectants, liquid soaps and lotions. Suitable for use with all our dispensers.

#### Good to know:

- packed individually for good hygiene, quick and easy use
- simple and hygienic installation: packaging is used as a glove
- no need to clean the pump or the dispenser
- cost optimisation: time-consuming dispenser reprocessing is not necessary
- no chance of accidental reuse: pump is destroyed when container is changed
- good for the environment: completely recyclable and produced using 100% renewable energy

REF	Product name
4401857	manual dispenser 500 ml
4401858	manual dispenser 1,000 ml
4401860	touchless dispenser 500 ml





# Our hygiene and cleaning products provide you with support in any application

The maintenance of a professional standard of cleaning and hygiene is key in healthcare sector as well as other industries. We aim to help you to achieve these hygiene standards at all times and in all areas with our range of cleaning and hygiene products.

Please note: Recommendations – such as in terms of material compatibility – can be provided separately. Our recommendations do not exempt our customers from testing our products independently to ensure they are suitable for the intended purposes and procedures. We cannot, therefore, assume any liability.

## **beopur** WC fresh

### Acidic, ready-to-use, viscous toilet cleaner

#### Features

Toilet fresh cleaner for cleaning toilets, urinals, bathtubs, shower cubicles and all acid-resistant surfaces in sanitary facilities.

#### Usage

\_

For limescale deposits, for removing limescale and intensive cleaning in the case of heavy soiling – can also be used undiluted. Then rinse thoroughly with water.

REF	Pack contents	Shipping unit	PZN
00926	750 ml	15 bottles	-



# **MEDIzid**<sup>®</sup>

### S-series sanitary cleaning

#### **Features**

- cleans WCs, urinals and sinks, tiles, bath and shower facilities, as well as all fittings with a chrome or stainless steel finish
- also usable as a WC cleaner, thanks to its viscosity
- removes soiling such as rust, urine stains, limescale and cement films
- complexing agent avoids the further formation of limescale
- inhibits further soiling
- usable on acid-resistant surfaces and floors. Not suitable for materials sensitive to acids, such as marble
- for metering, the use of our GASTROzid® KP10 canister pump, GASTROzid® AH10 outlet tap or GASTROzid® DK1000 dosing ball (for 1 | bottle) is recommended

#### **MEDIzid<sup>®</sup> S1 SANITARY SURFACE CLEANER** gel

1 | Spray bottle: S1-01 Pack: 12 btl. 10 I canister: S1-10




11 10 I







### Reliable hygiene in the sanitary sector



### **MEDIman**<sup>®</sup> **H-series hand cleaning**





#### MEDIman<sup>®</sup> H1 WASCHLOTION • perfume-free hand and all-over body wash lotion hypoallergenic perfume-free

500 ml Euro bottle: H1-500 1 | Euro bottle: H1-1000 Pack: 18 btl./pack: 10 btl.



• cleans stressed skin and hands

• hypoallergenic pH value maintains skin's natural protective acidic mantle



#### MEDIman<sup>®</sup> H2 HAND DISINFECTION "Limited virucide" rating

Short contact time 500 ml Euro bottle: H2-500 1 | Euro bottle: H2-1000 Pack: 18 btl./pack: 10 btl.



500 ml

٦Ì



#### MEDIman<sup>®</sup> H3 **SKIN PROTECTION LOTION** waterproof

1 | Euro bottle: H3-1000 Pack: 10 btl.





- waterproof skin protection and hand care cream
- for use <u>BEFORE</u> and <u>AFTER</u> work
- silicone- and perfume-free

#### • broad spectrum of action

- dose for hygienic hand disinfection is 3 ml (+15 s contact time)
- evaporates without leaving a residue

### Reliable hygiene in patient wards







### **MEDIzid® O-series surface cleaning**



#### MEDIzid<sup>®</sup> 01 **MULTI-PURPOSE CLEANER** broad spectrum

1 | Spray bottle: 01-01 Pack: 12 btl.

10 | canister: 01-10



10 I


ЭFГ

11

- for waterproof surfaces and floors, suitable for the one-stage wet and moist mopping/wiping of all types of coverings, does not result in layering
- removes soiling, also including soiling from rubber, ink and ballpoint pens
- preserves the absorbency of textile fibres in mop coverings and wiping cloths
- for metering, the use of our GASTROzid® PW4 bucket filler, GASTROzid® KP10 canister pump, GASTROzid® AH10 outlet tap or GASTROzid® DK1000 dosing ball (for 1 l bottle) is recommended

#### MEDIzid<sup>®</sup> O2 **ALL-PURPOSE CLEANER** AND DEGREASER

750 ml spray bottle: 02-750 Pack: 12 btl.



- usable as a universal cleaner for kitchen cleaning
- also suitable as a universal cleaner for industrial cleaning work
- non-aggressive also suitable for use on aluminium surfaces
- contains no perfumes or dyes

#### MEDIzid<sup>®</sup> 03 RAPID DISINFECTION

10

٦ř 750 ml

750 ml spray bottle: 03-750 Pack: 12 btl.



• ready-to-use, alcohol-based rapid disinfection • for rapid and effective surface disinfection

- for alcohol-resistant surfaces
- also for use in the food industry, in manufacturing and in public organisations (PT2/PT4)

#### MEDIzid<sup>®</sup> 04 **GLASS CLEANER**

1 | Spray bottle: 04-01 Pack: 12 btl. 10 | canister: 04-10



11

- ready-to-use rapid surface cleaner
- for rapid, daily routine cleaning
- for all waterproof smooth surfaces such as windows, glass, mirrors, tiles and painted surfaces, as well as stainless steel and chrome-finished surfaces
- for metering, the use of our GASTROzid® KP10 canister pump or GASTROzid® AH10 outlet tap is recommended

### **O-series surface cleaning**




#### MEDIzid<sup>®</sup> 06 STAINLESS STEEL CARE

500 ml spray bottle: O6-500 Pack: 12 btl.



- for the cleaning and care of stainless steel and aluminium surfaces
- removes minor dirt deposits, such as fingerprints
- creates a maintenance film that protects against moisture, corrosion and re-soiling (fingerprints)
- ideal for the upkeep of counter tops, work surfaces, equipment, fittings, lifts, etc.
- for heavily soiled surfaces, pre-cleaning with e.g. MEDIzid<sup>®</sup> O2 all-purpose cleaner and degreaser is necessary





MEDIzid<sup>®</sup> O6 STAINLESS STEEL CARE



MEDIzid<sup>®</sup> O2 ALL-PURPOSE CLEANER AND DEGREASER



MEDIzid<sup>®</sup> O4 GLASS CLEANER




MEDIzid<sup>®</sup> 03 RAPID DISINFECTION




MEDIzid<sup>®</sup> O1 MULTIPURPOSE CLEANER

### **Cleaning in commercial kitchens**



Gastrozid K1 kitchen cleaner – alkaline Gastrozid K2 kitchen cleaner – acidic

Medizid® Rapid+



Stovetop and surrounding area Gastrozid K1 kitchen cleaner – alkaline Gastrozid K2 kitchen cleaner – acidic Hot water devices Gastrozid K11 appliance descaler General surface disinfection Medizid® AF+

### **GASTROzid® GS-series dishwashing**



#### GASTROzid<sup>®</sup> GS1 **DISHWASHING DETERGENT** chlorinated

10 l /12 kg canister: GS1-10 20 l /25 kg canister: GS1-20

10 I

20 I

- alkaline detergent for commercial dishwashers with metering facilities
- removes starch and soiling such as coffee, tea or juice residues as well as lipstick
- · complexing agent reduces the further formation of limescale
- non-foaming
- helps to avoid deposits of grease in the dishwasher
- suitable for use with GASTROzid<sup>®</sup> GS2 rinse aid (acidic)
- for metering, use of the GASTROzid® PR4 metering pump is recommended



#### GASTROzid<sup>®</sup> GS2 **RINSE AID**

#### acidic

10 l /12 kg canister: GS2-10 20 l /25 kg canister: GS2-20







10 I



- alkaline dishwashing detergent and glassware cleaner for commercial dishwashers with metering facilities
- chlorine- and phosphonate-free.
- removes starch and soiling such as coffee, tea or juice residues as well as lipstick
- · helps to avoid deposits of grease in the dishwasher
- suitable for use with GASTROzid® GS2 rinse aid (acidic)
- for metering, use of the GASTROzid® PR4 metering pump is recommended



#### GASTROzid<sup>®</sup> GS3 **DISHWASHING DETERGENT** AND GLASSWARE CLEANER chlorine-free

10 l /12 kg canister: GS3-10 20 l /25 kg canister: GS3-20





#### GASTROzid<sup>®</sup> GS4 **RINSE AID** neutral

10 | canister: GS4-10



- neutral rinse aid for commercial and household machines as well as convection ovens or combined steamers if osmosis plants are used
- only suitable for use with fully demineralised water from an osmosis plant with max. 2-3 °dH; with hard water, always ensure you use GASTROzid® GS2 acidic rinse aid
- ensures dishes and glassware dry without spotting to a sparkling finish
- suppresses foaming
- suitable for use with GASTROzid® GS1 dishwashing detergent
- for metering, use of the GASTROzid® PR1 metering pump is recommended

- acidic universal rinse aid for commercial and household dishwashers • ensures dishes and glassware dry without spotting to a sparkling finish
- prevents limescale from forming
- suitable for use with GASTROzid® GS1 dishwashing detergent
- for metering, use of the GASTROzid® PR1 metering pump is recommended

### **GS-series dishwashing**



#### GASTROzid<sup>®</sup> GS5 **DISHWASHING DETERGENT** chlorine- and phosphate-free

10 l /12 kg canister: GS5-10



- alkaline detergent for commercial dishwashers with metering facilities
- free from substances that pollute wastewater, such as phosphate, NTA, NTDA and chlorine compounds
- removes starch and soiling such as coffee, tea or juice residues as well as lipstick
- complexing agent reduces the further formation of limescale
- non-foaming
- helps to avoid deposits of grease in the dishwasher
- suitable for use with GASTROzid® GS2 rinse aid (acidic)
- for metering, use of the GASTROzid® PR4 metering pump is recommended



#### GASTROzid<sup>®</sup> GS10 **HAND WASHING-UP LIQUID** powerful against grease

1 | Easy-grip bottle: GS10-01 Pack: 12 btl.

10 | canister: GS10-10



• rinses dishes, glassware, cutlery and cookware

- removes grease and soiling such as coffee, tea or juice residues as well as lipstick
- also suitable for use as a neutral surface cleaner
- no streaking and self-drying
- for manual dosing, the use of our GASTROzid<sup>®</sup> KP10 canister pump or GASTROzid® AH10 outlet tap is recommended

### **Metering system**



11



#### 1 GASTROzid<sup>®</sup> PR1 metering pump 230 V for rinse aid

- for under-sink and hood-type dishwashers
- usable with GASTROzid® GS2 or GS4 rinse aid

#### 1 pc. #PR1

#### 2 GASTROzid<sup>®</sup> PR4 metering pump 230 V for dishwashing detergent

- for under-sink and hood-type dishwashers
- usable with GASTROzid<sup>®</sup> GS1, GS3 or
  - GS5 dishwasher detergent
- 1 pc. #PR4

#### **3** GASTROzid<sup>®</sup> T50 combined metering unit for rinse aid and dishwasher detergent

- for hood-type and conveyor dishwashers
- usable with the GASTROzid<sup>®</sup> GS series
- 1 pc. #T50

#### – GASTROzid<sup>®</sup> SPG3 internal hose pump (not shown)

- for dishwashers with an internal dosing unit (as a replacement)
- optimum results obtained when used with the GASTROzid<sup>®</sup> GS series
- 1 pc. #SPG3



### **GASTROzid®** K-series kitchen cleaning



#### GASTROzid<sup>®</sup> K1 **KITCHEN CLEANER** ALKALINE

10 | canister: K1-10



10 I

- alkaline high-performance cleaner and degreaser
- for routine cleaning
- for water- and alkali-resistant floors and surfaces (e.g. work surfaces, stainless steel surfaces and fittings, kitchen appliances) in commercial kitchens or for universal use in the food processing industry
- removes oils, grease, blood, protein and other soiling
- also suitable as a universal cleaner for industrial cleaning work
- also suitable for use on aluminium surfaces
- · contains no perfumes or dyes

• acidic foam cleaner

• for basic and routine cleaning

• contains no perfumes or dyes

• additional note:

• removes limescale, rust and urine stains

• limescale inhibitor blocks further soiling

• also suitable for sanitary facilities or swimming pools

- additional note:
- a) For routine and basic cleaning: GASTROzid® K2 kitchen cleaner acidic
- b) For surface disinfection: GASTROzid® K3 disinfectant cleaner
- for metering, the use of the GASTROzid<sup>®</sup> KP10 canister pump, GASTROzid<sup>®</sup> PW4 bucket filler or GASTROzid<sup>®</sup> PW2 spray/foam unit is recommended

• suitable for acid-resistant floors (e.g. tiles) and surfaces (e.g. work surfaces, stainless steel surfaces and fittings, kitchen appliances) in commer-cial kitchens or for universal use in the food processing industry



#### GASTROzid<sup>®</sup> K2 **KITCHEN CLEANER** ACIDIC

10 | canister: K2-10



Di 10 I



#### GASTROzid<sup>®</sup> K3 DISINFECTANT **CLEANER**

slightly alkaline

10 | canister: K3-10





 grill and conventional/convection oven cleaning • suitable for all kitchen appliances, such as grills, combi-steamers, conventional and convection ovens, trollies and oven racks, baking trays, conveyor belts,

- GASTROzid® GS2 rinse aid acidic
- for metering, the use of the GASTROzid<sup>®</sup> KP10 canister pump, GASTROzid® PW4 bucket filler or GASTROzid® PW2 spray/foam unit is recommended



#### GASTROzid<sup>®</sup> K10 OVEN CLEANER

breaks down encrustations

750 ml spray bottle: K10-750 Pack: 12 btl. 10 | Canister: K10-10



10 I

750 ml

#### • slightly alkaline disinfectant cleaner • for disinfection and cleaning in one work step • suitable for water-resistant floors (e.g. tiles) and surfaces

b) For surface disinfection: GASTROzid® K3 disinfectant cleaner

• also suitable for use on aluminium surfaces, with corrosion protection

a) For alkaline routine cleaning: GASTROzid® K1 alkaline kitchen cleaner

for metering, the use of the GASTROzid<sup>®</sup> KP10 canister pump, GASTROzid<sup>®</sup> PW4 bucket filler or GASTROzid<sup>®</sup> PW2 spray/foam unit is recommended

- (e.g. work surfaces, stainless steel surfaces and fittings, kitchen appliances) in commercial kitchens or for use in the food processing industry
- bactericidal, limited-spectrum virucide and yeasticide
- also removes slight soiling involving oil and grease
- non-aggressive
- additional note:

a) For alkaline routine cleaning: GASTROzid® K1 alkaline kitchen cleaner b) For acidic routine and basic cleaning: GASTROzid® K2 kitchen cleaner acidic

- for metering, the use of the GASTROzid® KP10 canister pump, GASTROzid® PW4 bucket filler or GASTROzid® PW2 spray/foam unit is recommended
- appliances, machines, containers, worktops and floors • removes burnt-on encrustations
  - do not use on soft metals (aluminium)!
  - additional note for metering units built into the oven:

### **K-series kitchen cleaning**



#### GASTROzid<sup>®</sup> K11 APPLIANCE DESCALER

1 | Spray bottle: K11-01 Pack: 12 btl.



- acidic appliance descaler
- for descaling appliances and surfaces
- suitable for acid-resistant kitchen appliances (e.g. hot water appliances, coffee machines or pans) and surfaces (e.g. stainless steel surfaces and fittings) and floors
- removes limescale deposits from stainless steel, aluminium, brass, copper, tin and (to an extent) zinc
- caution! Do not use on marble

### Metering system and aids



#### 1 GASTROzid<sup>®</sup> PW2 spray/foam device

- for 2 products
- with 15 m hose and spray gun, canister brackets for standard water connection/<u>NO</u> electrical connection necessary.
- $\bullet$  for metering the entire GASTROzid  $^{\otimes}$  K series
- 1 pc. #PW2

#### 2 GASTROzid<sup>®</sup> PW4 bucket filler

- for 4 products
- with plunger, flow rate 16 l /min
- for standard water connection/<u>NO</u> electrical connection necessary.
- for metering the entire GASTROzid® K series
- 1 pc. #PW4

#### **3** GASTROzid<sup>®</sup> AH10 outlet tap

• for 10 l canister (51 mm DM) 1 pc. #AH10

#### 4 GASTROzid<sup>®</sup> KP10 canister pump

• for bucket metering (1 stroke = approx. 20 ml) 1 pc. #KP10

#### 5 GASTROzid<sup>®</sup> KS10 canister key

• for 10 | canister 1 pc. #KS10

#### 6 GASTROzid<sup>®</sup> DK1000 dosing ball

- for 1 I round bottles
- for bucket metering (reservoir = 20 ml)
- 1 pc. #DK1000


#### 7 GASTROzid<sup>®</sup> WH10/WH20 stainless steel wall bracket 101


• for wall storage of all GASTROzid<sup>®</sup> 10 | products 1 pc. #WH10

#### 20 I

- for wall storage of all GASTROzid<sup>®</sup> 20 | products
- 1 pc. #WH20

### Cleaning in other areas





### Cleaning and hygiene plan

what	with what		how	hazardous substances
Disinfect hands		MEDIman H2 hand disinfection	hygienic: Rub in 3 ml / 30 sec undiluted, lather with a little water, clean hands, rinse with water.	
Dishes, cutlery		Gastrozid GS5 dishwashing detergent	machine dishwashing	
Dishes, cutlery		Gastrozid GS2 rinse aid acidic	machine dishwashing	
Dishes, pots		Gastrozid GS10 hand washing-up liquid	machine dishwashing	
Stove, extractors, surroundings		Gastrozid K1 kitchen cleaner alkaline	Apply cleansing solution evenly, allow to act, work in, rinse off and, if necessary, wipe dry	
Water heaters		Gastrozid K11 appliance descaler	Add approx. 100 ml descaler per litre of water. Fill the water tank of the appliance with the solution 2/3. Heat up the appliance. Important! After each descaling, rinse the appliance sufficiently.	A A A A A A A A A A A A A A A A A A A
Surfaces, floors		Medizid® AF+ Concentrate for surface disinfection, VAH listed	Pre-clean surfaces after each use, apply the disinfectant solution, clean, reaction time min.15 min	
Floor cleaning		Gastrozid K1 kitchen cleaner alkaline Gastrozid K2 kitchen cleaner acidic Wechselreinigung	Apply product solution, work on it with a scrubber, rinse and wipe dry	
Stainless steel surfaces		Medizid O6 stainless steel care only vertical surfaces, <u>not for</u> food processing surfaces	Apply with clean, dry cloth; distribute the product evenly on the surface, do not polish again; <b>do not use on surfaces for food processing</b>	
Convection oven, grill		Gastrozid K10 OVEN CLEANER Grill cleaner with metal protection	If necessary, pre-heat the surface, pre-rinse, pre-steam, spray on the product, allow to act	
Surfaces, appliances		<b>Medizid®</b> Rapid+ rapid disinfection VAH listed	After each use, pre-clean and disinfect equipment and surfaces, spray on clean and apply evenly, reaction time 60 sec	

Surfaces, appliances



# Everyday products — and therefore essential!

Work processes in the healthcare sector should primarily be patient-centred. For this reason, we have a comprehensive range of ward and patient supplies in our portfolio to properly meet this everyday need. From daily physical hygiene for bed-ridden patients with our care products and washing gloves, to essential items such as cotton buds, wooden tongue depressors or protective towels.

### **BagBath**<sup>®</sup> haircare Washing hair without water

#### Features

A water-free hair washing system with no rinsing. Can be used either heated up or cooled down. Contains no perfumes. Can also be used in combination with the BagBath body wash system.

#### Usage

Hair washing system for bed-ridden patients and people requiring care.

REF	Pack contents	Shipping unit	PZN
APO021012	100 ml	30 bottles	_

## **BagBath**<sup>®</sup>

#### Disposable body wash system without water or soap

#### **Features**

The wipes in the BagBath<sup>®</sup> have a pleasant and subtle scent. Moisturising ingredients and clinically tested formulas for these ingredients work to protect the skin from irritation.



#### Usage

The packaging can be opened easily and then resealed, even with gloves. Wash bowls, soap, water, flannels, towels and care lotion are no longer required.

REF	Size	Pack contents	Shipping unit	PZN
2359	25 x 20 cm	8 wipes	30 packs	06149252






BagBath<sup>®</sup> enables hospitals and nursing care centres to improve the quality of skin care for their patients, and also while saving time and money.

BagBath<sup>®</sup> makes it possible to wash the body without soap, washcloths and towels in situations where there is either no bath available or it is not possible to get the patient to the bath.

#### All body parts are covered:

BagBath<sup>®</sup> ensures full and complete body hygiene. A pack contains 8 wipes – one for each area of the body.



#### Benefits for the patient:

- Avoidance of contact infections caused by traditional washing utensils such as flannels and towels.
- Prevention of cross-contamination by dispensing with reusable washing materials.
- Washing with BagBath<sup>®</sup> has fewer risks for patients who have undergone surgery, been intubated, suffered trauma or who are immobile, since there is no need to move the patient as is the case when the bed clothing becomes wet.

#### Special pack that can be warmed up:

#### Microwave + BagBath<sup>®</sup>

- Most patients prefer the feeling of a warm bath.
- Warmth prevents the patient losing body heat.



#### Special pack that can be chilled:

- Refrigerator + BagBath<sup>®</sup>
- A refreshing and cooling effect on patients



#### Simple to dispose of:

#### 🛥 rubbish bin

As adisposable product, BagBath<sup>®</sup> can be disposed of with normal household waste. When washing patients who have been quarantined, internal disposal regulations must be respected.



# BeeSana® Wash glove

Hygienic, disposable wash gloves

#### Features

• disposable wash glove in various designs





#### Usage

For skin cleansing and personal care. Can be used in all areas of nursing, with a particular focus on the care of the elderly and sick.

REF	Product name	Model	Pack contents	Shipping unit	PZN
3709	Film	Ultrasonically welded	50 pcs	20 packs	05461042
3723	Molton/soft	Ultrasonically welded	50 pcs	20 packs	11703466
3724	Molton/soft	Sewn	25 pcs	40 packs	05529001
3725	Molton 80 g	Sewn	50 pcs	20 packs	09336931

## **Protective towel**

#### Features

- protective towels made from a layer of cellulose and PE film
- PE inner coating lessens wetness leaking through onto the clothing
- with punched straps for fixing around the neck, plus a movable catch pouch lower down
- thanks to the soft cellulose outer layer, can also be used as a serviette after the meal
- latex-, PVC- and DEHP-free
- colour: Apricot

#### Usage

Disposable.

REF	Size	Pack contents	Shipping unit	PZN
2997	38 x 68 cm	100 pcs	6 packs	11727314



# Protective bib

#### Features

- protective bib made from medical gauze
- 6-ply
- white

REF	Size	Pack contents	Shipping unit	PZN
4239	20 x 30 cm	10 pcs	10 packs	12724401



## **BeeSana**<sup>®</sup> Disposable briefs

Briefs to give a feeling of secure privacy

#### Features

- disposable briefs made from non-woven material
- with rubber drawstring
- white with coloured waistband for easy sizing

#### Usage

To protect the privacy of patients during examinations and after surgery.

REF	Size	Waist colour	Waist size	Pack contents	Shipping unit	PZN
2251	М	Blue	51–55 cm	20 pcs	30 packs	11857318
2252	L	Brown	54–58 cm	20 pcs	30 packs	11857324
2253	XL	Green	58–62 cm	20 pcs	30 packs	11857330
2254	XXL	Red	62–64 cm	20 pcs	30 packs	11857301





# Anti-slip sock universal

### Our sock to prevent slips and falls

#### Features

- disposable ABS sock (laid flat)
- with grips as an anti-slip coating
- made from cotton
- black
- universal size
- washable at 30 °C



#### Usage

Disposable, usable for many different applications in hospitals and nursing homes.

REF	Size	Pack contents	Shipping unit	PZN
9605	40 x 9 cm	20 pcs	8 packs	12671811

# Surgical sock

#### Features

- disposable sock (laid flat), without heel
- top with rubber drawstring
- white
- washable at 30 °C



Suitable for many applications. Disposable, for use in the operating theatre, lab or clean rooms.

REF	Size	Pack contents	Shipping unit	PZN
9837	39 x 8 cm	40 pcs	20 packs	12671805
9836	44 x 8 cm	40 pcs	20 packs	12671797



## PE shoe covers

Shoe covers made from CPE with rubber drawstring

#### Features

- with rubber drawstring
- disposable
- colour: blue

Dimensions (approx.):	Length:	41 cm
	Height:	15 cm
	Thickness:	0.03 mm

REF	Pack contents	Shipping unit	PZN
70-151	1000 pcs	1 box	11363007



## Cellulose sponge sheets

#### Features

• highly bleached (German pharmacopoeia (DAB) standard)

absorbent

#### Usage

Cleaning and hygiene cloth for medical applications.

REF	Size	Pack contents	Shipping unit	PZN
12431	20 x 20 cm	300 g	26 packs	13655576

## Layered cellulose

Crepe cuts in loose layers for multiple uses

#### Features

- layered cellulose, unbleached, made from 100% recycled cellulose
- many different uses

#### Usage

Can be used as an underlay for non-sterile activities, as an absorbent material or for cleaning and wiping.

REF	Size	Pack contents	Shipping unit	PZN
102200D	20 x 20 cm	5 kg	3 bags	06803349
102300D	20 x 30 cm	5 kg	3 bags	06803355
104600D	40 x 60 cm	5 kg	3 bags	06803361



# Rollicel<sup>®</sup>/medical crepe

#### Exam table covers on a roll

#### Features

- exam table covering made from cellulose
- absorbent and micro-perforated
- several available variants
- colour: white

#### Usage

For treatment tables, and partially-clothed examinations and treatments.

REF	Size	Layers	Coating	Shipping unit	PZN
3739-50	39 cm x 50 m	2-ply	No	9 rolls	11350476
3739-50B	39 cm x 50 m	2-ply	Yes	6 rolls	11350571
3739-100	39 cm x 100 m	2-ply	No	6 rolls	11350536
3745-100B	45 cm x 100 m	2-ply	Yes	6 rolls	11350884
3750-50	50 cm x 50 m	2-ply	No	9 rolls	11350482
3750-50B	50 cm x 50 m	2-ply	Yes	6 rolls	11350855
3750-100	50 cm x 100 m	2-ply	No	6 rolls	11350542
3755-50	55 cm x 50 m	2-ply	No	9 rolls	11350499
3755-50B	55 cm x 50 m	2-ply	Yes	6 rolls	11350861
3755-100	55 cm x 100 m	2-ply	No	6 rolls	11350559
3759-50	59 cm x 50 m	2-ply	No	9 rolls	11350513
3759-50B	59 cm x 50 m	2-ply	Yes	6 rolls	11350878
3759-100	59 cm x 100 m	2-ply	No	6 rolls	11350565
3759-100B	59 cm x 100 m	2-ply	Yes	6 rolls	11350890

## **Protective sheets**

#### Disposable product as a surface for sitting or lying

#### Features

- 2-ply cellulose with transparent PE film that prevents fluids from seeping through
- latex-, DHEP- and PVC-free.
- reinforced lengthwise with rugged polyester stitching and side barriers.

#### Usage

Disposable protective sheets to cover stretchers, examination tables and gurneys, work tables and mattresses.

REF	Size (cm)	Shipping unit	PZN
1988	80 x 175	100 pcs	04353787
1992	80 x 200	100 pcs	04353793
1995	80 x 210	100 pcs	04353801



## **Protective covering**

Disposable product as a surface for sitting or lying, the cost-effective alternative to protective sheets

#### Features

- 1-ply cellulose with transparent PE film that prevents fluids from seeping through
- latex-, DHEP- and PVC-free.

#### Usage

Disposable protective sheets to cover stretchers, examination tables and gurneys, work tables and mattresses.

REF	Size (cm)	Side barriers	Shipping unit	PZN
1993	38 x 40	No	1,000 pcs	04353681
1990	80 x 90	Yes	200 pcs	04353735
2982	80 x 140	Yes	100 pcs	01082431
2984	80 x 210	Yes	100 pcs	01082454

## Multi-purpose cloths

#### Features

- multi-purpose cloth made from thick felted non-woven
- tear-resistant, absorbent and lint-free

#### Usage

For all cleaning work in the kitchen, laboratory, hospitals and hotel rooms.

REF	Cloth format	Pack contents	Shipping unit	PZN
3539	32 x 38 cm	50 pcs	20 packs	_



## Protective bed coverings

### Transparent HDPE

#### Features

protective HDPE (high-density polyethylene) cover

#### Usage

For hospital and nursing beds, to protect the fresh, newly-made bed.

REF	Cloth format	Pack contents	Shipping unit	PZN
1985	320 x 95 cm	200 pcs	1 roll	_



## **Universal cloth**

#### Features

- disposable product for cleaning and hygienic care of intact skin
- high-viscose version also optionally available

#### Usage

Disposable product for use in all areas of care.

**Universal cloth** 

REF	Size	Pack contents	Shipping unit	PZN
6051	30 x 30 cm (approx., unfolded) 15 x 15 cm (approx., folded)	50 pcs	10 packs	04355415

### Universal soft wipe (thicker quality)

REF	Size	Pack contents	Shipping unit	PZN
6052	30 x 30 cm (approx., unfolded) 15 x 15 cm (approx., folded)	50 pcs	10 packs	11722759

## Wipe cloths


### Recyclable wipes made from cellulose

#### Features

- multi-layered, semi-bleached wipes made from cellulose
- latex-, PVC- and DEHP-free

#### Usage

Disposable for use as a cleaning material.







## **Cosmetic wipes**

#### Features

- made from 2-ply cellulose
- colour: brilliant white

REF	Size	Pack contents	Shipping unit	PZN
6701	20.5 x 21 cm	100 pcs	35 packs	-



# V-fold paper hand towels

#### Features

• dimensions: 25 x 23 cm

1-ply, nat	ural		
REF	Pack contents	Shipping unit	PZN
5722	20 x 250 pcs	1 pack	_
2-ply, wh	ite		
REF	Pack contents	Shipping unit	PZN
5721	20 x 160 pcs	1 pack	_



## **BaSick Bag**

### Sick bag

#### Features

- for use with the onset of nausea
- partially transparent, tear-proof sick bag made from polyethylene
- with graduation and seal ring

#### Usage

Disposable, for use with nausea.







REF	Size, height	Capacity	Pack contents	Shipping unit	PZN	
9401	16 x 40 cm	1,500 ml	25 pcs	20 bags	13818825	
SP9401	47 cm	3 packs	1 pc.	1 рс.	_	_

## Kidney dish

### Non-sterile, paperboard

REF	Size	Colour	Shipping unit	PZN
9400	255 x 150 x 48 mm	Grey	300 pcs	01082164

### Double sterile, plastic, pack of 1 pc.

REF	Size	Pack contents	Shipping unit	PZN
4249	700 ml	12 x 1 pc.	2 packs	_





## **Disposable comb**

Plastic comb, individually wrapped in cellophane

0670 Prown 100 x 1 nc 20 nocks 119576	Shipping unit PZN	Pack contents	Colour	REF
9079 DIOWII 100 X I JUC. 20 JACKS 110570	20 packs 11857672	100 x 1 pc.	Brown	9679



## **Disposable razor**

#### Features

• stainless steel, non-sterile disposable razor

REF	Razor blades	Pack contents	Shipping unit	PZN
9610	Single-edged	100 pcs	10 packs	02501189

# Urine bottle

### Urine bottle for men 1,000 ml

#### Features

• opaque, graduated urine bottle with clip fastener

REF	Shipping unit	PZN
3006	25 pcs	_





## Urine bottle holder

### Urine bottle holder with chain and lid

#### Features

• urine bottle holder for securely picking up the urine bottle

#### Usage

The urine bottle holder should attached (hung on) within the reach of the patient to the bed frame of the patient's bed. Clean as necessary. <u>IMPORTANT:</u> When using disinfectants, always respect the manufacturer's instructions for use.

REF	Unit	Shipping unit	PZN	
3001	1 pc.	20 pcs	_	_



## Nurse's stethoscope

### Nurse's stethoscope for auscultation

#### Features

- acoustic flat-head stethoscope
- mechanical
- ear hooks with earpieces, tube and diaphragm
- tube length: approx. 60 cm
- latex-free

#### Usage

The stethoscope must only be used by authorised medical personnel. The stethoscope is placed on the part of the body relevant for the exam and used by exerting only a little pressure.

REF	Pack contents	Shipping unit	PZN
9805	1 pc.	100 pcs	11857459

## Tourniquet

## High-quality tourniquet with innovative method of use

- usable with one hand, no "reaching around"
- with safety fastener
- contains no latex or metal



REF	Article	Colour	PZN
5020	Meditrade tourniquet	Colourful designer band	06177892
N20200	Jetpull	Blue designer band	-



## **Dosicair**<sup>®</sup>

### Our infusion regulator with injection adapter

#### Features

- luer taper connection system
- can be connected to infusion systems

#### Usage

For infusion treatments.

REF	Pack contents	Shipping unit	PZN
DF100	50 pcs	8 packs	_

## Urine beaker

REF	Size	Pack contents	Shipping unit	PZN
3009	125 ml	100 pcs	5 packs	04354427
3010	100 ml	500 pcs	1 pack	-





# Ultrasound scan gel

#### Features

• ultrasound gel for scans using ultrasound probes

#### Usage

For all types of ultrasound scans.

REF	Pack contents	Shipping unit	PZN
6831	250 g	30 bottles	01904778



## **Cold-Warm-Swab**

### Pad for cold and warm applications

#### Features

- permanently gel-filled pad
- can be easily shaped to fit
- suitable for the freezer (stays flexible to -30 °C), water bath and microwave
- follow information on label

#### Usage

For strains, pulled muscles and swelling. Also for fever, headaches and bruising.



	REF	Size	Pack contents	Shipping unit	PZN
Gel colour: transparent	9668	13 x 14 cm	50 pcs	1 pack	04352500
	9870	12 x 29 cm	25 pcs	1 pack	04352463
-	9666	16 x 26 cm	25 pcs	1 pack	04352517
	9661	30 x 40 cm	6 pcs	1 pack	04352546
	9663	7.5 x 35 cm	40 pcs	1 pack	04352552
Gel colour: blue	9465	16 x 25 cm	20 pcs	1 pack	-

## Non-woven cover for pads

### Non-woven cover for cooling/heating pads

- non-woven covers for cooling/heating pads
- used to protect the skin, since direct contact between the pad and skin must be avoided
- for disposable use only

REF	Size	Pack contents	Shipping unit	PZN
9682	13 x 14 cm	500 pcs	1 pack	_
9683	12 x 29 cm	500 pcs	1 pack	_
9684	16 x 26 cm	500 pcs	1 pack	_
9687	21 x 38 cm	500 pcs	1 pack	_
9693	30 x 40 cm	500 pcs	1 pack	_

## **Rescue blanket**

#### Features

- colour: gold/silver
- reflects bodily heat or solar radiation
- reduces the risk of hypothermia or works to shield the body from the heat of the sun's rays
- packed individually with instructions for use.

#### Usage

Cold and heat protection for injured parties. Usable for accidents in windy, rainy or hot conditions.



REF	Size	Pack contents	Shipping unit	PZN
4070	160 x 210 cm	50 pcs	5 packs	13510895

## **Medical beaker**

### Medical beaker, 30 ml, graduated

- made from transparent polypropylene
- contents visible
- with graduation at side

	REF	Colour	Pack contents	Shipping unit	PZN
Medical beaker					
	9729	Transparent	75 pcs	50 packs	_
	9727	Blue	75 pcs	50 packs	-
	9725	Red	75 pcs	50 packs	-
	9726	Yellow	75 pcs	50 packs	_
	9728	Green	75 pcs	50 packs	-
Medical beaker lid					
	9700	White	750 pcs	5 packs	_

## **Medication dispenser**

#### Features

- tablet dispenser made from PP/PVC with transparent sliding lid
- helps patients to take their medication regularly and on time in the morning, at lunchtime, in the evening and at night

REF	Colour	Pack contents	PZN
9707	White	200 pcs	_
9708	Blue	200 pcs	_



## Medi-Box

### Container for securely disposing of needles and syringes

#### Features

- round disposal container
- with lid and round opening
- flat closure
- pull-off and screw-off aids

#### Usage

For disposing of needles and syringes.

REF	Size	Pack contents	PZN
9694	0.5	1 pc.	-
9702	2.5 l	1 pc.	_



## Wall holder

### Wall holder for glove boxes of all sizes

- acrylic or polystyrene
- easy to fit









REF	Article	Product name	Dimensions (in mm)	Shipping unit	PZN
4040	Dracula	PS wall holder, without overprint (white)	Universal	1 pc.	-
4050	Ghost	Acrylic wall holder with overprint	250 x 145 x 100	1 pc.	_
4059	Ghost	Acrylic wall holder without overprint	250 x 140 x 70	1 pc.	_
4060	Ghost	Acrylic wall holder without overprint	250 x 145 x 100	1 pc.	_

## Rubbish bags HD

### PE rubbish bags

#### Features

- rubbish bags made from HDPE (high-density polyethylene)
- for collecting general waste
- disposable

REF	Size	Colour	Pack contents	Shipping unit	PZN
9747	49 x 60 cm (approx. 24 l)	Transparent	50 pcs	40 rolls	-
9748	63 x 74 cm (approx. 60 l)	Transparent	50 pcs	40 rolls	-

# Rubbish bags LD

### PE rubbish bags

- rubbish bags made from LDPE (low-density polyethylene)
- for collecting general waste
- disposable

REF	Size	Colour	Pack contents	Shipping unit	PZN
9651	70 x 110 cm (approx. 120 l)	Blue	25 pcs	10 rolls	-
9746	50 x 60 cm (approx. 30 l)	Grey	50 pcs	20 rolls	-
9749	63 x 72 cm (approx. 60 l)	Transparent	50 pcs	20 rolls	-
9750	63 x 85 cm (approx. 65 l)	White	40 pcs	9 rolls	-
9760	57.5 x 100 cm (approx. 70 l)	Blue	25 pcs	10 rolls	-
9761	70 x 110 cm (approx. 120 l)	Blue	25 pcs	10 rolls	-
9762	70 x 110 cm (approx. 120 l)	Transparent	25 pcs	10 rolls	-
9764	60 x 80 cm (approx. 60 l)	Grey	25 pcs	10 rolls	-
9765	70 x 110 cm (approx. 120 l)	Red	25 pcs	10 rolls	-
9766	57.5 x 100 cm (approx. 70 l)	Yellow	25 pcs	10 rolls	_
9768	70 x 110 cm (approx. 120 l)	Blue	15 pcs	10 rolls	-



## Kitchen towel

#### Features

- brilliant white
- made from cellulose
- 50 sheets per roll
- roll length: 11.2 m
- sheet width: 25.7 cm
- sheet length: 22.4 cm

### 3-ply

REF	Pack contents	Shipping unit	PZN
57301	4 rolls	8 packs	-

## **Toilet** paper

#### Features

- brilliant white
- made from cellulose
- 250 sheets per roll
- roll length: 28 m
- sheet width: 9.5 cm
- sheet length: 11.2 cm

### 2-ply

REF	Pack contents	Shipping unit	PZN
9200	8 rolls	8 packs	_
3-ply			
REF	Pack contents	Shipping unit	PZN







Notes:	

Notes:	
	<u> </u>

### Index

#### Α

ABE® Absorbent adhesive dressing 81 ABE<sup>®</sup> body bandage 99 ABE<sup>®</sup> cannula retention dressing 83 ABE<sup>®</sup> derm (plus) 82 ABE<sup>®</sup> drainal 71 ABE<sup>®</sup> extra 94 ABE<sup>®</sup> extra adhesive 95 ABE<sup>®</sup> fix color 89 ABE<sup>®</sup> gauze bandage 86 ABE<sup>®</sup> Gilchrist's bandage 98 ABE<sup>®</sup> heart catheter pressure bandage 100 ABE® Ideal 100 91 ABE<sup>®</sup> injection plaster 80 ABE<sup>®</sup> Last 87 ABE<sup>®</sup> Last adhesive 90 ABE<sup>®</sup> Last crepe 88 ABE<sup>®</sup> lastic 93 ABE<sup>®</sup> light 80 ABE<sup>®</sup> limb holder 100 ABE<sup>®</sup> non-woven fabric 83 ABE<sup>®</sup> nose bandage 74 ABE<sup>®</sup> padded bandage 94 ABE® pressure bandage roll 101 ABE<sup>®</sup> ribbon gauze 76 ABE<sup>®</sup> silk 84 ABE® ribbon gauze X-ray contrast 109 ABE® tracheal 71 ABE<sup>®</sup> tubular bandage 96 ABE® tubular net bandage 97 ABE<sup>®</sup> universal 92 ABF<sup>®</sup> Ven 95 Absorbent cushions 121 Absorbent dressing pads 72 Absorbent pads 122 Active gel 147 Alcohol swabs 77 All-purpose cleaners 160-161 Anti-slip sock universal 176 Application aids 149-151 Arm bandages 98 Astronaut-shaped surgical cap 55

#### В

BagBath<sup>®</sup> (Haircare) 172 Bandage change sets 114 Bandages 86-99 BaSick Bag 182 BeeSana<sup>®</sup> wash gloves 174 BeeSana® absorbent cushion 121 BeeSana<sup>®</sup> absorbent dressing pad 72 BeeSana<sup>®</sup> absorbent pad 122 BeeSana<sup>®</sup> bandage change set 114 BeeSana<sup>®</sup> cellulose swab 77 BeeSana<sup>®</sup> Comflock Super 127 BeeSana<sup>®</sup> disposable brief 175 BeeSana<sup>®</sup> eye pad 73 BeeSana<sup>®</sup> female incontinence pad 123 BeeSana® gauze balls 75 BeeSana<sup>®</sup> gauze ball with X-ray contrast 106 / 109 BeeSana<sup>®</sup> gauze swab with X-ray contrast 104 - 105 BeeSana<sup>®</sup> Lap sponge with X-ray contrast 102-103 BeeSana® male incontinence pad 123 BeeSana<sup>®</sup> medical gauze 64 BeeSana<sup>®</sup> Medi-Pants 119 BeeSana® non-woven swab 69 BeeSana<sup>®</sup> PEG bandage set 113 BeeSana<sup>®</sup> pointed swab 76 BeeSana® pointed swab with X-ray contrast 108 BeeSana<sup>®</sup> PP gown 56 BeeSana® PP scrubs 56 BeeSana<sup>®</sup> PP/PE gown 57 BeeSana<sup>®</sup> pre-cut gauze 64 BeeSana<sup>®</sup> preparation balls with X-ray contrast 107 BeeSana<sup>®</sup> rectangular pad 122 BeeSana® SafeCut 85 BeeSana® Slit gauze swab 68 BeeSana® Slit non-woven swab 70 BeeSana® SMS gown 57-58 BeeSana<sup>®</sup> suture removal set 114 BeeSana<sup>®</sup> ward catheter set 112 BeeSana® ward protection set 115 BeeSana® Wash glove 174

beopur WC fresh 154 Body bandages 99 Body wash systems 172 Box holder 195

### С

Canister keys 149 Canister metering pump 149 Canister/outlet taps 149 Cannula fixations/plasters 83 Care foam 148 Cellulose sponge sheets 177 Cellulose swabs 77 Chemical protective suit 58 Cleaning products 154-169 Cleaning products and disinfectants 156-169 Cold-Warm-Swab 186 Combs 182 Comflock Super 127 Compression bandages 91-95 Copolymed<sup>®</sup> 30 Copolymer gloves 30 Cosmetic wipes 181 Cotton buds 78 Cotton gloves 26 Cream lotion 146

### D

Dental masks 48-54 Disinfectant cleaner 162-175 Disinfectant dispenser 150-151 Disinfectant wipes 135-136, 139-140 Disinfection (products) 130-143 Disposable aprons 59 Disposable body wash system 172 Disposable briefs 175 Disposable comb 182 Disposable gloves 9-47 Disposable razors 183 Disposal container 194 Dosicair<sup>®</sup> 185 Dosing aids 149 Dosing bottle 149 Dressing material 64-109 Dressing pads 71

#### E

Ethasept<sup>®</sup> 132 Ethasept<sup>®</sup> Gel 133 Exam table covers 178 Examination gloves, (non-)sterile 9-25, 30-32 Eye pads 73

#### F

Face masks 49-54 Female incontinence pad 123 Female pads 123 Film wound bandages 82 Finessis<sup>®</sup> Aegis<sup>®</sup> 45 Finessis<sup>®</sup> Corium<sup>®</sup> 47 Finessis® Zero® 46 Finger cots 28 Fixation bandages 85-90 Fixation briefs 124-125 Fixation dressings 97 Fixation plasters 83-84 Flexylon gloves 44-47 Flock (patient) undersheets 128-129 Flock incontinence briefs 124 Foam spray heads 149 Froggy<sup>™</sup> 27

#### G

Gauze 64 Gauze balls 75 / 109 Gauze bandage 86 Gauze swabs 65-73, 104-105 Gentle Med<sup>®</sup> active gel 147 Gentle Med<sup>®</sup> care foam 154 Gentle Med<sup>®</sup> cream lotion (W/O) 146 Gentle Med<sup>®</sup> moist care wipes 148 Gentle Med® oil bath 147 Gentle Med<sup>®</sup> skin care lotion (O/W) 146 Gentle Med<sup>®</sup> skin cleansing and care foam 148 Gentle Med® wash lotion 145 Gentle Skin<sup>®</sup> Aloecare<sup>™</sup> 13 Gentle Skin<sup>®</sup> black 14 Gentle Skin<sup>®</sup> classic (x-long) 10 Gentle Skin<sup>®</sup> compact + 11 Gentle Skin<sup>®</sup> grip 11 Gentle Skin® Isopretex® OP 42 Gentle Skin® Micro OP® 39 Gentle Skin<sup>®</sup> Premium OP<sup>™</sup> 36 Gentle Skin® Securitex OP 41 Gentle Skin<sup>®</sup> sensitive 9 Gentle Skin® Smooth OP 38 Gentle Skin® sterile 31 Gentle Skin<sup>®</sup> Superior OP 37 Gentle Tec<sup>®</sup> 12 Gilchrist's bandage 98 Glove box holder 195 Glove liners/undergloves 26, 30, 39, 47 Gloves 8-47 Gowns 56-58 Graduated flasks 149

#### Η

Hair wash system 172 Hand disinfection 132 / 133 Hand spray head 149 Heart catheter pressure bandage 100 Household gloves 29 Hygiene gloves 27 Hygiene products 154-169

### I

Ideal bandages 91 Incontinence briefs 118-120 Incontinence briefs 118-119 Incontinence briefs (with leg) 126-127 Incontinence lines 120 Incontinence pads 120-123 Incontinence products 118-127 Incontinence undersheets 126-127 Infusion regulators 187 Injection plasters 80 Instrument disinfection 142-143

#### K

Kidney dishes 182 Kitchen towel 191

#### L

Lap sponge with X-ray contrast 102-103 Latex gloves 9-14, 28-29, 31, 35-41 Layered cellulose 177 Layered undersheets 126 Leg cuffs 100 Liberty Light 120 Light fixation briefs 124 Limb holder 100 Long-stretch bandages 97

## Index

#### Μ

Male incontinence pads 123 Masks 50-54 Maxitex<sup>®</sup> UG 40 Medi-Box 188 Medi-Pants 119 Medical beakers 187 Medical crepe 180 Medical gauze 64 Medication dispenser 188 Medium-stretch bandages 92 Medizid<sup>®</sup> AF+ 137 Medizid<sup>®</sup> Alpha+ 138 Medizid<sup>®</sup> Alpha+ Flowpack 140 Medizid<sup>®</sup> Alpha+ XXL (disinfectant wipes) 139 Medizid<sup>®</sup> Rapid QF Flowpack 136 Medizid<sup>®</sup> Rapid+ 134 Medizid<sup>®</sup> Rapid+/XXL (disinfectant wipes) 135 MEDIzid<sup>®</sup> sanitary cleaner 155 Metering pumps 149 Moisturising wipes 148 Mouth tongue depressor 79 Multi-purpose cloths 179 MyClean® easy touch 29

#### Ν

Nasal slings 74 Needle collection boxes 188 Neoprene gloves 42 Neopretex<sup>®</sup> OP 43 Nitril<sup>®</sup> 3000 15 Nitril<sup>®</sup> 3000 X-Long Blue 18 Nitril<sup>®</sup> BestGen<sup>®</sup> 17 Nitril<sup>®</sup> black 19 Nitril<sup>®</sup> LIOX 22 / 23 Nitril<sup>®</sup> Magenta 21 Nitril<sup>®</sup> NextGen<sup>®</sup> 15 Nitril<sup>®</sup> Sensory<sup>®</sup> 16 Nitril<sup>®</sup> sterile 32 Nitril<sup>®</sup> Viola 20 Nitrile gloves 15-23, 32 Non-slip socks 176

Non-woven (slit) swab 69-70 Non-woven cover for dressing pads 186 Non-woven fabric 83 Nose bandages 74 Nurse's caps 55 Nurse's stethoscopes 184 Nursing pads 128

### 0

Oil bath 147

#### Ρ

Padded bandages 94 Pads for post-natal women 122-123 Panty fixation briefs 125 Paper hand towels 181 PE aprons 59 PE disposable aprons 59 PE gloves 27, 30 PE shoe covers 177 PEG bandage set 113 Plasters 80-84 Pointed swabs 76, 108 Polychloroprene gloves 43 Polyisoprene gloves 42 PP/(PE) gowns/scrubs/aprons 56-58 Preparation balls 107 Pressure bandages 101 Protective bed coverings 179 Protective bibs 175 Protective coverings 179 Protective scrubs/gowns/aprons 56-61 Protective sheets 178 Protective towels 174

### R

Razors 183 Ready-to-use finger bandages 85 Rectangular pads 122 Reference<sup>™</sup> 9 Reference<sup>™</sup> OP 35 Rescue blankets 187 Ribbon gauze with X-ray 107 Roll plasters 84 Rollicel<sup>®</sup>/medical crepe 178 Rubbish bags 190

### S

SafeCut 85 Scrotum roll 101 Set solutions 112-115 Shoe covers 177 Short-stretch bandages 91, 93-95 Shoulder bandage 98 Sick bags (BaSick Bag) 182 Skin care lotion 146 Skin cleansing foam 148 Skin disinfection 136 Slit gauze swab 68 Slit (hole) dressing pads 71 SMMS surgical scrubs 60-61 SMS gowns 57-58 Soap dispensers 150-151 Socks 176 Stethoscopes 184 Suavel<sup>®</sup> Antifluid 53 Suavel<sup>®</sup> Antifog 53 Suavel<sup>®</sup> Apollo 55 Suavel<sup>®</sup> Astrid 55 Suavel<sup>®</sup> Comfort 52 Suavel<sup>®</sup> Comfort Plus 52 Suavel<sup>®</sup> OP PRO Expert 61 Suavel<sup>®</sup> OP PRO Standard 60 Suavel<sup>®</sup> Protec 51 Suavel<sup>®</sup> Protec Black 51 Suavel<sup>®</sup> Protec Plus 50 Suavel<sup>®</sup> Sensima 54 Suavel<sup>®</sup> Yeti Coverall 58 Superform<sup>®</sup> 120 Superform® Air 121

Superslip<sup>®</sup> 118 Superslip<sup>®</sup> Air Comfort 119 Support bandages 96-101 Surface disinfection 137-143 Surgical caps 55 Surgical caps 55 Surgical dental masks 49-54 Surgical face masks 49-54 Surgical gloves 35-47 Surgical gowns/scrubs/aprons 56-61 Surgical masks 49-54 Surgical sets 110-115 Surgical socks 176 Suture removal sets 114 Swabs 75-77, 106-109

#### T

Tablet dispensers 188 Toilet cleaner 154 Toilet paper 191 Tourniquet 184 Tube bandages 98-99 Tubular bandage 96 Tubular net bandage 97

### U

Ultrasound scan gel 185 Undersheets 126-127 Universal bandages 92-93 Universal wipes 180 Urine beaker 185 Urine bottle (holder) 183

#### V

V-fold paper hand towels 181 Vinyl 2000 24 Vinyl 2000 Stretch PF 25 Vinyl Blue Gen PF 25 Vinyl gloves 24-25 Visitor gowns 56-58

#### W


Wall holders 189 Ward catheter sets 112 Ward gowns/aprons 56-58 Ward protection set 115 Wash gloves 174 Wash lotions 145 Wet (moisturising) wipes 150 Wipe cloths 180 Wipes 135-136, 139-140 Wooden tongue depressor 79 Wound bandages 80-82 Wound plasters 80-82

### Х

X-ray wound care 102-109

### Y

Yarn glove 26



#### Senator Hans-Hermann Rösner-Mautby

(\* 27 June 1917 in Beuthen, † 24 June 1999 in Oberaudorf)

For his services towards the establishment of the hospital in Oberaudorf, he was made honorary senator of Ludwig-Maximilian University, Munich. He was also awarded the German Federal Cross of Merit and the Bavarian Order of Merit for his entrepreneurial achievements in the healthcare sector.



### Touch. Help. Create opportunities.

With its wide-ranging activities and carefully targeted public relations work, the Senator Rösner Foundation collects funds for people who have fallen into financial difficulty because of cancer. Founders Peter Rösner and his wife Eveline Rösner use their foundation to support cancer research and the training of medical personnel.

As a successor to the "Hilfswerks der Klinik Bad Trissl" (Friends of Bad Trissl Hospital), the Senator Rösner Foundation can now look back on over 25 years of experience as a charitable organisation.

#### The following projects are currently receiving long-term support.

#### **HEAT** study

The Senator Rösner Foundation is supporting the clinical research study led by Prof. Dr. med. Dipl.-Biochem. Rolf D. Issels at Munich University Hospital with its innovative approaches to the field of pancreatic cancer; funding volume here is several hundred thousand euros.

#### Rapid help where it's needed most

Cancer can leave some people facing unexpected financial and social difficulties. Here, the Senator Rösner Foundation "Leben mit Krebs" (Living with Cancer) provides quick and straightforward help.

### Promoting skills and professional development

The Senator Rösner Foundation takes advantage of its proximity to all the major oncological centres in Munich and Bavaria. The Foundation maintains regular contact with scientists and researchers while also providing funds for research and development.

#### We run to help

The charity run organised by the Senator Rösner Foundation is held under this banner, plus the patronage of FISWorld Champion Felix Neureuther, and is now well-known throughout Germany. In 2013 alone, almost 800 runners together covered the amazing distance of 6,756 km, and managed to amass a record-breaking sponsorship figure of  $\in$  105,000.

#### WSG-Swarovski-Wattens

#### **SOS Kinderdörfer**

Apioneering aid project that goes beyond national boundaries. Thanks to the bonds of friendship between the Rösner family and the charity president, Diana Langes-Swarovski, children and young people are helped via SOS Kinderdörfer to get through difficult life crises, including by other means than sporting development.

www.wsg-wattens-fussball.at

Contact: senator-roesner-stiftung.de

## **Company development**

- 1982 Formation of Hermes Trans Atlantik Trading as a subsidiary of our Canadian companies
- 1988 Formation of Meditrade GmbH, with its present-day business focus and slogan: "360° protection" – hand protection, face masks, caps, operating scrubs, disinfection, incontinence products, bandages, dressing pads, swabs, abdominal dressings, etc.
- 1989 50% participation in Meditrade Austria (divested 1992), and shareholding in Meditrade Switzerland (divested 1992)
- 1990 Turnover transfer from Cosmotec and Sauter Trading
- **1990–1992** Acquisition of first glove factory in Thailand, Hycare Ltd. Hatyai
- 1992 Formation of Meditrade Dresden (now Praximed Meditrade DD)
- 1994 Acquisition of AiA Leinberger GmbH Großräschen
- 1992 Relocation to new company site at Thierseestr. 196–198, Kiefersfelden
- 1993 Acquisition of Hese GmbH in Feucht (Nuremberg), cellulose and undersheet maker, plus VP-Schickedanz/Temca factory
- **1992–1994** Formation of new companies Hycare Ltd. UK/Barrier Healthcare Ltd. UK in Cambridge
- 1995 Asset deal with MCM Landhut



- 1994–1995 Breakeven point reached
- 1996 Shareholding acquired in Eurocare GmbH and Mediclean Bayern GmbH
- 1996 Formation of Meditrade International Thailand Ltd. with subsequent purchase of shareholding in joint venture Siamount Plastic Ltd. Bangkok with Mitsubishi Group (vinyl factory)
- Up to 2002 Range continuously extended with exclusivity agreements and minority shareholdings in Italy
- 2003 Construction of another new factory in Thailand Hatyai, Rösner-Mautby Thailand Ltd. (DVD)
- Up to 2006 Sale of all factories in Thailand
- 2008 Share investments in Advent Health Konzern AG, Malaysia (global leader in surgical glove production)
- 2009–2011 Acquisition of Beese GmbH & Co. KG in Hamburg with production subsidiary in China
- 2012 Move to new company premises at Medipark 1, Kiefersfelden
- 2016 Majority shareholding in MaiMed Group
- 2018 Expansion of logistics centre in Neuenkirchen, Meditrade GmbH celebrates 30 years in business

2020

Formation of Meditrade UK, a Meditrade subsidiary based in London: Meditrade UK Ltd., 176 – 178 Brompton Road, SW3 1HQ, London, England www.meditradeuk.com



#### Meditrade UK Ltd.

176 – 178 Brompton Road SW3 1HQ, London England Phone: +44 (0) 20 4513 2390 Email: info@meditradeuk.com

Please contact our Customer Service if you have any questions about our products or require further information: Phone: +44 (0) 20 4513 2390 Email: info@meditradeuk.com

You can also visit our website at www.meditradeuk.com The website provides full details of our product portfolio as well as other information.

#### **Meditrade GmbH**

Medipark 1 83088 Kiefersfelden Germany Phone: +49(0)8033/9760-0 Fax: +49(0)8033/9760-70





